



Pacific Medical Catalogue

TABLE OF CONTENTS

Blood Chairs & Stools 2

 Blood Collection Chair 2

 Round Stools 3

 Saddle Stools 6

 Step Up Stools 8

 Portable Massage Chair 10

BLOOD CHAIRS & STOOLS

BLOOD COLLECTION CHAIR



BLOOD COLLECTION CHAIR

Comfortable and supportive blood collection chair including arm rests and a foldable foot rest.

Specifications:

Comfortable soft, supportive high density foam

Extra wide Seat: 55cm

Width including arm rests: 91.5cm

Arm rests width: 18.5cm

Length: 65cm

Length including foot rests: 84cm

Height: 104cm

Foldable foot rest

Weight capacity: 175kg

ORDER CODE:

BCCGR - GREY

BCCBL - BLACK

[Click For More Information](#)

ROUND STOOLS



ROUND STOOL WITH BACKREST

High quality round stool with backrest. It has adjustable backrest angle along with height. Suitable for many applications.

Features:

- Adjustable backrest tilt
- Height adjustable backrest
- Comfortable foam seat

Specifications:

- Adjustable height: 47cm - 63cm
- Seat diameter: 40cm
- Backrest size: 36.5 x 26.5cm
- Backrest height range: 26.5 - 39cm
- Backrest angle adjustment: 25 degrees
- Foam thickness: 5cm
- Weight Capacity: 160kg

ORDER CODES:

RSNBRR - Navy Blue

RSBLBR - Black

RSGRBR - Grey

[Click For More Information](#)

ROUND STOOL WITH ARMREST

High quality round stool with armrest. It has a swing around armrest with height adjustment. Suitable for many applications.

Features:

- Swing around armrest
- Height adjustable armrest
- Comfortable foam seat

Specifications:

- Adjustable height: 48cm - 64cm



Seat diameter: 40cm
Armrest size: 53 x 8 x 5cm
Armrest height range: 23 - 31cm
Foam thickness: 7cm
Weight Capacity: 160kg

ORDER CODES:**RSNBAR - Navy Blue****RSBLAR - Black****RSGRAR - Grey**[Click For More Information](#)

ROUND STOOL STANDARD

A great value round stool suitable for many applications and practices.

Specifications:

Adjustable height: 51cm - 67cm
Cushion diameter: 35cm
Foam thickness: 6cm
Weight Capacity: 160kg

ORDER CODES:**RSNB - Navy Blue****RSBL - Black****RSGR - Grey**[Click For More Information](#)

ROUND STOOL PREMIUM WITH FOOT RING

A premium round stool with a height adjustable foot ring and strong base.

Specifications:

Adjustable height: 52cm - 72cm
Cushion diameter: 35cm
Foam thickness: 6cm
Weight Capacity: 160kg
Base: 56cm
Adjustable height foot ring

ORDER CODES:**PRSNB - Navy Blue**

PRSBL - Black
PRSGR - Grey

[Click For More Information](#)

SADDLE STOOLS



SADDLE STOOL WITH BACKREST

High quality saddle stool with backrest. It has adjustable seat and backrest angle along with backrest height and angle. Suitable for many applications.

Features:

- Adjustable seat tilt
- Height and angle adjustable backrest
- Moulded foam seat

Specifications:

- Adjustable height: 47cm - 63cm
- Seat diameter: 42cm
- Backrest size: 36.5 x 26.5cm
- Backrest height range: 23 - 38cm
- Seat angle range: -3 to +6 degrees
- Backrest angle adjustment: 25 degrees
- Foam thickness: 5cm
- Weight Capacity: 160kg

ORDER CODES:

SSNBBR - Navy Blue

SSBLBR - Black

SSGRBR - Grey

[Click For More Information](#)

SADDLE STOOL STANDARD

A great value saddle stool suitable for many applications and practices.

Specifications:

Adjustable height: 51cm - 67cm
Cushion diameter: 35cm
Foam thickness: 6cm
Weight Capacity: 160kg

ORDER CODES:

SSNB - Navy Blue

SSBL - Black

SSGR - Grey

[Click For More Information](#)



SADDLE STOOL PREMIUM WITH FOOT RING

A premium saddle stool with a height adjustable foot ring and strong base.

Specifications:

Adjustable height: 51cm - 71cm
Cushion diameter: 35cm
Foam thickness: 6cm
Weight Capacity: 160kg
Base: 56cm
Adjustable height foot ring

ORDER CODES:

PSSNB - Navy Blue

PSSBL - Black

PSSGR - Grey

[Click For More Information](#)

STEP UP STOOLS

SINGLE STEP STOOL



Perfect for any examination couch this sturdy single step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Dimensions: 455x245x200 (LxWxH)
Weight Capacity: 150KG

ORDER CODE: SST

[Click For More Information](#)

DOUBLE STEP STOOL

Perfect for any examination couch this sturdy double step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

Specifications:

Stainless steel frame
Non-slip ribbed rubber platform
Reinforced rubber tips
Manufactured with 1" steel tubing
Overall Dimensions: 400x550x470mm (LxWxH)
First Step: 250mm
Second Step: 220mm
Weight Capacity: 150KG

ORDER CODE: DST

[Click For More Information](#)



PORTABLE MASSAGE CHAIR



PORTABLE MASSAGE CHAIR

Perfect for on the go practitioners our strong & luxurious portable massage chair requires no tools for assembly and folds up easily.

Features:

- Strongly made steel frame
- Removable chest pad, seat pad, knee support, face rest
- Comfortable design with High density foam
- Luxurious washable & durable PU vinyl
- No tools needed for assembly
- Folds up easily

Specifications:

- Open dimensions: 54x76x120-127cm
- Folded dimensions: 54x29x118cm
- Chair weight: 15kg
- Weight capacity: 180kg
- Warranty: 2 Years on frame

ORDER CODE:
PMCBL - BLACK

[Click For More Information](#)

INDEX

B

Blood Collection Chair 2

D

Double Step Stool 8

P

Portable Massage Chair 10

R

Round Stool Premium With Foot Ring 4

Round Stool Standard 4

Round Stool With Armrest 3

Round Stool With Backrest 3

S

Saddle Stool Premium With Foot Ring 7

Saddle Stool Standard 7

Saddle Stool With Backrest 6

Single Step Stool 8



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au