



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Round Stools** ..... 2

# ROUND STOOLS



## ROUND STOOL WITH BACKREST

High quality round stool with backrest. It has adjustable backrest angle along with height. Suitable for many applications.

### Features:

Adjustable backrest tilt  
Height adjustable backrest  
Comfortable foam seat

### Specifications:

Adjustable height: 47cm - 63cm  
Seat diameter: 40cm  
Backrest size: 36.5 x 26.5cm  
Backrest height range: 26.5 - 39cm  
Backrest angle adjustment: 25 degrees  
Foam thickness: 5cm  
Weight Capacity: 160kg

### ORDER CODES:

**RSNBBR - Navy Blue**

**RSBLBR - Black**

**RSGRBR - Grey**

[Click For More Information](#)

## ROUND STOOL WITH ARMREST

High quality round stool with armrest. It has a swing around armrest with height adjustment. Suitable for many applications.

### Features:

Swing around armrest  
Height adjustable armrest  
Comfortable foam seat

### Specifications:

Adjustable height: 48cm - 64cm



Seat diameter: 40cm  
Armrest size: 53 x 8 x 5cm  
Armrest height range: 23 - 31cm  
Foam thickness: 7cm  
Weight Capacity: 160kg

**ORDER CODES:****RSNBAR - Navy Blue****RSBLAR - Black****RSGRAR - Grey**[Click For More Information](#)**ROUND STOOL STANDARD**

A great value round stool suitable for many applications and practices.

**Specifications:**

Adjustable height: 51cm - 67cm  
Cushion diameter: 35cm  
Foam thickness: 6cm  
Weight Capacity: 160kg

**ORDER CODES:****RSNB - Navy Blue****RSBL - Black****RSGR - Grey**[Click For More Information](#)**ROUND STOOL PREMIUM WITH FOOT RING**

A premium round stool with a height adjustable foot ring and strong base.

**Specifications:**

Adjustable height: 52cm - 72cm  
Cushion diameter: 35cm  
Foam thickness: 6cm  
Weight Capacity: 160kg  
Base: 56cm  
Adjustable height foot ring

**ORDER CODES:****PRSNB - Navy Blue**

**PRSB** - Black  
**PRSG** - Grey

[Click For More Information](#)

# INDEX

## R

Round Stool Premium With Foot Ring 3

Round Stool Standard 3

Round Stool With Armrest 2

Round Stool With Backrest 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)