



Pacific Medical Catalogue

TABLE OF CONTENTS

Saddle Stools 2

SADDLE STOOLS



SADDLE STOOL WITH BACKREST

High quality saddle stool with backrest. It has adjustable seat and backrest angle along with backrest height and angle. Suitable for many applications.

Features:

- Adjustable seat tilt
- Height and angle adjustable backrest
- Moulded foam seat

Specifications:

- Adjustable height: 47cm - 63cm
- Seat diameter: 42cm
- Backrest size: 36.5 x 26.5cm
- Backrest height range: 23 - 38cm
- Seat angle range: -3 to +6 degrees
- Backrest angle adjustment: 25 degrees
- Foam thickness: 4cm
- Weight Capacity: 160kg

ORDER CODES:

SSNBBR - Navy Blue

SSBLBR - Black

SSGRBR - Grey

[Click For More Information](#)

SADDLE STOOL STANDARD

A great value saddle stool suitable for many applications and practices.

Specifications:

Adjustable height: 51cm - 67cm

Cushion diameter: 35cm

Foam thickness: 4cm

Weight Capacity: 160kg

ORDER CODES:

SSNB - Navy Blue

SSBL - Black

SSGR - Grey

[Click For More Information](#)



SADDLE STOOL PREMIUM WITH FOOT RING

A premium saddle stool with a height adjustable foot ring and strong base.

Specifications:

Adjustable height: 51cm - 71cm

Cushion diameter: 35cm

Foam thickness: 4cm

Weight Capacity: 160kg

Base: 56cm

Adjustable height foot ring

ORDER CODES:

PSSNB - Navy Blue

PSSBL - Black

PSSGR - Grey

[Click For More Information](#)

INDEX

S

Saddle Stool Premium With Foot Ring 3

Saddle Stool Standard 3

Saddle Stool With Backrest 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au