



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Step Up Stools** ..... 2

# STEP UP STOOLS

## SINGLE STEP STOOL



Perfect for any examination couch this sturdy single step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

### Specifications:

Stainless steel frame  
Non-slip ribbed rubber platform  
Reinforced rubber tips  
Manufactured with 1" steel tubing  
Dimensions: 455x245x200 (LxWxH)  
Weight Capacity: 150KG

### ORDER CODE: SST

[Click For More Information](#)

## DOUBLE STEP STOOL

Perfect for any examination couch this sturdy double step stool is made from 1" stainless steel tubing with reinforced rubber tips & a non slip ribbed rubber platform.

### Specifications:

Stainless steel frame  
Non-slip ribbed rubber platform  
Reinforced rubber tips  
Manufactured with 1" steel tubing  
Overall Dimensions: 400x550x470mm (LxWxH)  
First Step: 250mm  
Second Step: 220mm  
Weight Capacity: 150KG

### ORDER CODE: DST

[Click For More Information](#)





# INDEX

## **D**

Double Step Stool 2

## **S**

Single Step Stool 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)