



Pacific Medical Catalogue

TABLE OF CONTENTS

Consumables 2

CONSUMABLES

PILLOW CASE WITH FLAP BLUE BX200



Easy to use disposable pillow case with flap that provide reliable protection and help maintain a clean, sanitary patient environment.

Specifications:

75x50cm in size
SMS 40GSM
With flap
Blue in colour
Box of 200

ORDER CODES: PCFLAP

[Click For More Information](#)

PVC PILLOW CASE WITH ZIP LIGHT BLUE BX100

Water-resistant, easy-to-use disposable PVC pillow cases with a zip that provide reliable protection and help maintain a clean, sanitary patient environment.

Specifications:

75x50cm in size
Lightweight PVC
Zip for easy closure
Light blue in colour
Box of 100

ORDER CODES: PCPVCZIP

[Click For More Information](#)





PATIENT GOWN WITH SHORT SLEEVES BX100

High-quality disposable patient gowns offering lightweight comfort and dependable coverage. Perfect for maintaining cleanliness and convenient patient care in healthcare environments.

Specifications:

M (110x125cm) & XL (120x140cm) sizes available
Non-woven SMS
Neck & back closed with ties
Dark blue in colour
45 GSM
Box of 100

ORDER CODES:

PGSSLM - Medium

PGSSLXL - Extra Large

[Click For More Information](#)

MODESTY SHEET 100X70CM BX200

A dark blue, durable & lightweight modesty sheet designed to provide comfort and privacy during treatments and procedures.

Specifications:

100x70cm in size
SMS 40GSM
Dark blue in colour
Box of 200

ORDER CODES:

MSH100

[Click For More Information](#)





MODESTY SHEET 200X120CM BX100

A dark blue, durable & lightweight modesty sheet designed to provide comfort and privacy during treatments and procedures.

Specifications:

200x120cm in size
SMS 25GSM
Dark blue in colour
Box of 100

ORDER CODES: MSH200

[Click For More Information](#)

PACIFIC TOWEL ROLL CTN8

Pacific towel roll an easy and economical absorbent perforated paper roll. Many applications as a hand towel, liner, underlay, surface & device protection.

Features:

Easy & economical roll format
Perforated
Absorbent
Multi Use

Specification:

2 Ply
100 sheets per roll
Individually shrink wrapped rolls
Carton of 8 rolls
Length: 41.5m
Width: 49cm

ORDER CODES: PMPTL

[Click For More Information](#)





PACIFIC BED ROLL CTN6

Pacific towel roll an easy and economical absorbent perforated paper roll. Many applications as a hand towel, liner, underlay, surface & device protection.

Features:

- Soft yet thick & durable
- Continuous roll
- Strong & hygenic cover

Specification:

- Individually shrink wrapped rolls
- Carton of 6 rolls
- Length: 80m
- Width: 54cm

ORDER CODES: PMPBRL

[Click For More Information](#)

UNDERPAD BX300

A high quality waterproof underpad suitable for use in the most demanding healthcare setting. 5 Ply for greater absorbency and waterproof backing allow better containment of fluids.

Features:

- Waterproof
- Multi layer construction
- Light blue in colour
- Super absorbent
- Comfortable
- Environmentally friendly

Specification:

- Fluff pulp mix
- Absorbent polymer
- Waterproof PE Film
- Non skid backing



Box of 300 pieces
Packaged in 6 bags of 50 pieces
Sealed on all four sides
Length: 60cm
Width: 40cm

ORDER CODES:
PMUPL

[Click For More Information](#)



HALF SIZE UNDERPAD BX600

A high quality waterproof half size underpad suitable for use in the most demanding healthcare setting. 5 Ply for greater absorbency and waterproof backing allow better containment of fluids.

Features:

- Waterproof
- Multi layer construction
- Light blue in colour
- Super absorbent
- Comfortable
- Environmentally friendly

Specification:

- Fluff pulp mix
- Absorbent polymer
- Waterproof PE Film
- Non skid backing
- Box of 600 pieces
- Packaged in 12 bags of 50 pieces
- Sealed on all four sides
- Length: 40cm
- Width: 30cm

ORDER CODES:
PMUPHF

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 57X135MM BX200

A simple high quality self seal autoclave pouch 57x135mm in size. It features internal and external sterilisation indicators and a strong triple seal.



Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP57x135

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 70X260MM BX200

A simple high quality self seal autoclave pouch 70x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP70x260

[Click For More Information](#)



SELF SEALING AUTOCLAVE POUCH 90X160MM BX200

A simple high quality self seal autoclave pouch 90x160mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP90x160

[Click For More Information](#)



SELF SEALING AUTOCLAVE POUCH 90X260MM BX200

A simple high quality self seal autoclave pouch 90x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP90x260

[Click For More Information](#)





SELF SEALING AUTOCLAVE POUCH 135X285MM BX200

A simple high quality self seal autoclave pouch 135x285mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP135x285

[Click For More Information](#)



SELF SEALING AUTOCLAVE POUCH 190X360MM BX200

A simple high quality self seal autoclave pouch 190x360mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP190x360

[Click For More Information](#)



SELF SEALING AUTOCLAVE POUCH 300X380MM BX200

A simple high quality self seal autoclave pouch 300x380mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:

SSP300x380

[Click For More Information](#)



FITTED DISPOSABLE BEDSHEET BX100

The best form of protection for your examination couch is a disposable fitted bed sheet. Suitable for various treatment couches and stretchers, made to size and packaged in four packs of 25 for easy distribution it is the preferred solution for many healthcare facilities.

Specifications:

- 2000 x 750mm
- Non-woven polypropylene
- Available in navy blue & white
- 40 GSM
- Box of 100 (4 x 25)

ORDER CODES:

BSHFITWH - White

BSHFITNB - Navy Blue

[Click For More Information](#)



NON-FITTED DISPOSABLE BEDSHEET BX100

The best form of protection for your examination couch is a disposable flat bed sheet. Suitable for various treatment couches, made to size and packaged in ten packs of 10 for easy distribution it is the preferred solution for many healthcare facilities.



Specifications:

2400 x 700mm
Non-woven polypropylene
Available in navy blue & white
40 GSM
Box of 100 (4 x 25)

ORDER CODE:

BSHWH - White

BSHNB - Navy Blue

[Click For More Information](#)

INDEX

F

Fitted Disposable Bedsheet BX100 10

H

Half Size Underpad BX600 6

M

Modesty Sheet 100x70cm BX200 3

Modesty Sheet 200x120cm BX100 4

N

Non-Fitted Disposable Bedsheet BX100 11

P

Pacific Bed Roll CTN6 5

Pacific Towel Roll CTN8 4

Patient Gown With Short Sleeves BX100 3

Pillow Case With Flap Blue BX200 2

PVC Pillow Case With Zip Light Blue BX100 2

S

Self Sealing Autoclave Pouch 135x285mm BX200 9

Self Sealing Autoclave Pouch 190x360mm BX200 9

Self Sealing Autoclave Pouch 300x380mm BX200 10

Self Sealing Autoclave Pouch 57x135mm BX200 7

Self Sealing Autoclave Pouch 70x260mm BX200 7

Self Sealing Autoclave Pouch 90x160mm BX200 8

Self Sealing Autoclave Pouch 90x260mm BX200 8

U

Underpad BX300 5



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au