



Pacific Medical Catalogue

TABLE OF CONTENTS

Bipolar Foreceps & Accessories 2

BIPOLAR FORCEPS & ACCESSORIES



BOVIE REUSABLE BIPOLAR FORCEPS CORD FOR NEW UNITS

This cable allows connection of current generation Bovie electrosurgical units to bipolar forceps.

This cable allows connection of current Bovie electrosurgical units to bipolar forceps. The cable is 3M long and suited to devices such as the Derm 101, Derm 102, Derm 942 (A942), Bantam Pro (A952) & Specialist Pro (A1250S)

ORDER CODES:
BOVA827V

[Click For More Information](#)

BOVIE 184MM FINE TIP BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Bovie 184mm fine tip bipolar forcep are guaranteed for twenty autoclave cycles. Bipolar forceps require an electrosurgery footswitch for activation, see all Bovie accessories for the available range. For care and handling, refer to the Instructions for Use in the package insert.

Features:

- High quality materials
- Guaranteed for 20 autoclave cycles
- Require an electrosurgery footswitch for activation

ORDER CODES:
BOVA820

[Click For More Information](#)



BOVIE 177MM STRAIGHT FINE TIP BIPOLAR FORCEP



Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Bovie 177mm straight fine tip bipolar forcep are guaranteed for twenty autoclave cycles. Bipolar forceps require an electrosurgery footswitch for activation, see all Bovie accessories for the available range. For care and handling, refer to the Instructions for Use in the package insert.

Features:

- High quality materials
- Guaranteed for 20 autoclave cycles
- Require an electrosurgery footswitch for activation

ORDER CODES:
BOVA822

[Click For More Information](#)

BOVIE 114MM STRAIGHT MICRO TIP BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Bovie 114mm straight micro tip bipolar forcep are guaranteed for twenty autoclave cycles. Bipolar forceps require an electrosurgery footswitch for activation, see all Bovie accessories for the available range. For care and handling, refer to the Instructions for Use in the package insert.

Features:

- High quality materials
- Guaranteed for 20 autoclave cycles
- Require an electrosurgery footswitch for activation

ORDER CODES:
BOVA826

[Click For More Information](#)



OLSEN 178MM CUSHING SMOOTH STRAIGHT 1.5MM TIPS INSULATED BIPOLAR FORCEP



Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Olsen 178mm cushing smooth straight 1.5mm tips are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials
Insulated
Suitable for use with all electrosurgical devices

ORDER CODES:
BOV10-11621

[Click For More Information](#)

OLSEN 127MM ADSON SMOOTH 1.0MM TIPS INSULATED BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Olsen 127mm adson smooth 1.0mm tips are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials
Insulated
Suitable for use with all electrosurgical devices

ORDER CODES:
BOV10-10821

[Click For More Information](#)





BOVIE REUSABLE BIPOLAR FORCEPS CORD FOR OLD UNITS

This cable allows connection of older Bovie electrosurgical units to bipolar forceps.

This cable allows connection of older Bovie electrosurgical units to bipolar forceps. It is 3M long and suited to devices such as the A900 series.

ORDER CODES:
BOVA827BP

[Click For More Information](#)

OLSEN 184MM GERALD SMOOTH 1.0MM TIPS INSULATED BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Olsen 184mm Gerald Smooth 1.0mm are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials
Insulated
Suitable for use with all electrosurgical devices

ORDER CODES:
BOV10-11801

[Click For More Information](#)



OLSEN 140MM SEMKEN SMOOTH 0.75MM TIPS INSULATED BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Olsen 140mm semken smooth 0.75mm tips are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials
Insulated
Suitable for use with all electrosurgical devices

ORDER CODES:
BOV10-10911

[Click For More Information](#)

SYMMETRY 101MM JEWELERS 0.4MM TIPS INSULATED BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Symmetry 101mm jewelers 0.4mm tips are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials
Insulated
Suitable for use with all electrosurgical devices

ORDER CODES:
BOV08-00011

[Click For More Information](#)





OLSEN 127MM ADSON SERRATED 1.0MM TIPS INSULATED BIPOLAR FORCEP

Reusable Bipolar Forceps are available in micro .5mm, 1.0mm, and 2.0mm tip size and in lengths ranging from 3.5" to 8".

Olsen 127mm adson serrated 1.0mm tips are insulated bipolar forceps manufactured to industry gold standard utilising the highest quality materials.

Features:

High quality materials

Insulated

Suitable for use with all electrosurgical devices

ORDER CODES:

BOV10-10811

[Click For More Information](#)

INDEX

B

Bovie 114mm Straight Micro Tip Bipolar Forcep 3
Bovie 177mm Straight Fine Tip Bipolar Forcep 3
Bovie 184mm Fine Tip Bipolar Forcep 2
Bovie Reusable Bipolar Foreceps Cord For New Units 2
Bovie Reusable Bipolar Foreceps Cord For Old Units 5

O

Olsen 127mm Adson Serrated 1.0mm Tips Insulated Bipolar Forcep 7
Olsen 127mm Adson Smooth 1.0mm Tips Insulated Bipolar Forcep 4
Olsen 140mm Semken Smooth 0.75mm Tips Insulated Bipolar Forcep 6
Olsen 178mm Cushing Smooth Straight 1.5mm Tips Insulated Bipolar Forcep 4
Olsen 184mm Gerald Smooth 1.0mm Tips Insulated Bipolar Forecep 5

S

Symmetry 101mm Jewelers 0.4mm Tips Insulated Bipolar Forcep 6



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au