



Pacific Medical Catalogue

TABLE OF CONTENTS

Cautery 2

 Power Handles & Replacement Tips 2

 High Temperature Cauteries 7

 Low Temperature Cauteries 11

 Sleeves 13

CAUTERY

POWER HANDLES & REPLACEMENT TIPS

BOVIE POWER HANDLE LOW-TEMP



Bovie manufactures more cauteries than anyone in the world. Therefore, Bovie has the ability to put the right cautery in the physician's hand regardless of the procedure. High temperature cauteries are popular for many types of procedures including, pinpoint hemostasis, controlling diffuse bleeders, and many others.

Bovie power handle low-temp replaceable battery and replaceable tip cauteries are the next advancement in battery operated cauteries where reusability is a must. These products are designed for the office, ER, and clinical procedures. The Bovie power handle low-temp can be purchased individually in high or low temperature models with a variety of tips designed to fit your procedural needs.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Reusable handle with switch
- Replaceable batteries
- Replaceable tips

Specification:

Dimensions

Length: 205mm
Diameter: 19mm

Battery

1 x AA

Kit Inclusions

HIT0 power handle low-temp x 1
H100 high-temp fine tip x 1
AA battery x 1

ORDER CODES:

BOVHIT1

MAINTENANCE: Please make sure your device is checked and serviced yearly by an authorised technician.

[Click For More Information](#)

BOVIE POWER HANDLE HIGH-TEMP

Bovie manufactures more cauteries than anyone in the world. Therefore, Bovie has the ability to put the right cautery in the physician's hand regardless of the procedure. High temperature cauteries are popular for many types of procedures including, pinpoint hemostasis, controlling diffuse bleeders, and many others.



Bovie power handle high-temp replaceable battery and replaceable tip cauteries are the next advancement in battery operated cauteries where reusability is a must. These products are designed for the office, ER, and clinical procedures. The Bovie power handle high-temp can be purchased individually in high or low temperature models with a variety of tips designed to fit your procedural needs.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Reusable handle with switch
- Replaceable batteries
- Replaceable tips

Specification:

Dimensions

Length: 205mm
Diameter: 19mm

Battery

2 x AA

Kit Inclusions

HIT1 power handle high-temp x 1
H101 high-temp fine tip x 1
AA battery x 2

ORDER CODES:

BOVHIT1

MAINTENANCE: Please make sure your device is checked and serviced yearly by an authorised technician.

[Click For More Information](#)

BOVIE REPLACEMENT TIP LOW-TEMP FINE BX10



Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Features:

- Easy to use
- Sterile drape included
- Suitable for power handles

ORDER CODES: BOVH100

[Click For More Information](#)

BOVIE REPLACEMENT TIP HIGH-TEMP FINE BX10

Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Replacement tip high-temp fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Features:

- Easy to use
- Sterile drape included
- Suitable for power handles

ORDER CODES: BOVH101

[Click For More Information](#)



BOVIE REPLACEMENT TIP HIGH-TEMP LOOP BX10



Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Replacement tip high-temp loop are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Features:

- Easy to use
- Sterile drape included
- Suitable for power handles

ORDER CODES: BOVH103

[Click For More Information](#)

BOVIE REPLACEMENT TIP HIGH-TEMP FINE ELONGATED BX10

Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.


Replacement tip high-temp elongated fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Features:


- Easy to use
- Sterile drape included
- Suitable for power handles

ORDER CODES: BOVH104

[Click For More Information](#)



BOVIE REPLACEMENT TIP HIGH-TEMP VASECTOMY BX10



Replacement Tip Low-Temp Fine are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Replacement tip high-temp vasectomy are sterile cautery tips including a sterile drape to cover the non-sterile handle. Tips are individually packaged, 10 per box. The replacement tips are used in conjunction with the power handle cauteries.

Features:

- Easy to use
- Sterile drape included
- Suitable for power handles

ORDER CODES: BOVH121

[Click For More Information](#)

HIGH TEMPERATURE CAUTERIES

BOVIE CAUTERY HIGH-TEMP FINE TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

ORDER CODES: BOVAA01

[Click For More Information](#)



BOVIE CAUTERY HIGH-TEMP LOOP TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp loop tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.



Features:

Portable & easy to use
Safe design with recessed power button
Easy snap design

**ORDER CODES:
BOVAA03**

[Click For More Information](#)

**BOVIE CAUTERY HIGH-TEMP LOOP TIP EXTENDED
SHAFT BX10**

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp loop tip extended shaft has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

Features:

Portable & easy to use
Safe design with recessed power button
Easy snap design

**ORDER CODES:
BOVAA05**

[Click For More Information](#)

BOVIE CAUTERY HIGH-TEMP FINE TIP EXTENDED SHAFT BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.



Bovie cautery high-temp fine tip extended shaft has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to "snap off" the end of the cautery for safe and effective removal of the battery power source.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

ORDER CODES: BOVAA17

[Click For More Information](#)

BOVIE CAUTERY HIGH-TEMP ELONGATED TIP BX10



Bovie has been the number one manufacturer of battery-operated cauterics for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp elongated tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to "snap off" the end of the cautery for safe and effective removal of the battery power source.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

ORDER CODES:
BOVAA21

[Click For More Information](#)

LOW TEMPERATURE CAUTERIES

BOVIE CAUTERY LOW-TEMP FINE TIP BX10



Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery low-temp fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

ORDER CODES: BOVAA00

[Click For More Information](#)

BOVIE CAUTERY LOW-TEMP MICRO FINE TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery low-temp micro fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.



Features:

Portable & easy to use
Safe design with recessed power button
Easy snap design

ORDER CODES:**BOVAA90**[Click For More Information](#)

SLEEVES

BOVIE CAUTERY STERILE SHEATH BX10

Bovie cautery sterile sheath suitable for use with power handles.



Bovie cautery sterile sheath suitable for use with power handles.

ORDER CODES:
BOVHISL

[Click For More Information](#)

INDEX

B

Bovie Cautery High-Temp Elongated Tip BX10 9
Bovie Cautery High-Temp Fine Tip BX10 7
Bovie Cautery High-Temp Fine Tip Extended Shaft BX10 9
Bovie Cautery High-Temp Loop Tip BX10 7
Bovie Cautery High-Temp Loop Tip Extended Shaft BX10 8
Bovie Cautery Low-Temp Fine Tip BX10 11
Bovie Cautery Low-Temp Micro Fine Tip BX10 11
Bovie Cautery Sterile Sheath BX10 13
Bovie Power Handle High-Temp 3
Bovie Power Handle Low-Temp 2
Bovie Replacement Tip High-Temp Fine BX10 4
Bovie Replacement Tip High-Temp Fine Elongated BX10 5
Bovie Replacement Tip High-Temp Loop BX10 5
Bovie Replacement Tip High-Temp Vasectomy BX10 6
Bovie Replacement Tip Low-Temp Fine BX10 4



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au