



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**High Temperature Cauteries ..... 2**

# HIGH TEMPERATURE CAUTERIES

## BOVIE CAUTERY HIGH-TEMP FINE TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

### Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

**ORDER CODES:**  
**BOVAA01**

[Click For More Information](#)



## BOVIE CAUTERY HIGH-TEMP LOOP TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp loop tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.



**Features:**

Portable & easy to use  
Safe design with recessed power button  
Easy snap design

**ORDER CODES:  
BOVAA03**

[Click For More Information](#)

**BOVIE CAUTERY HIGH-TEMP LOOP TIP EXTENDED  
SHAFT BX10**

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp loop tip extended shaft has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

**Features:**

Portable & easy to use  
Safe design with recessed power button  
Easy snap design

**ORDER CODES:  
BOVAA05**

[Click For More Information](#)

## BOVIE CAUTERY HIGH-TEMP FINE TIP EXTENDED SHAFT BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.



Bovie cautery high-temp fine tip extended shaft has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

### Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

### ORDER CODES:

**BOVAA17**

[Click For More Information](#)

## BOVIE CAUTERY HIGH-TEMP ELONGATED TIP BX10



Bovie has been the number one manufacturer of battery-operated cauterics for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery high-temp elongated tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

### Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

**ORDER CODES:**  
**BOVAA21**

[Click For More Information](#)

# INDEX

## B

Bovie Cautery High-Temp Elongated Tip BX10 4

Bovie Cautery High-Temp Fine Tip BX10 2

Bovie Cautery High-Temp Fine Tip Extended Shaft BX10 4

Bovie Cautery High-Temp Loop Tip BX10 2

Bovie Cautery High-Temp Loop Tip Extended Shaft BX10 3



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)