



Pacific Medical Catalogue

TABLE OF CONTENTS

Low Temperature Cauteries 2

LOW TEMPERATURE CAUTERIES

BOVIE CAUTERY LOW-TEMP FINE TIP BX10



Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery low-temp fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.

Features:

- Portable & easy to use
- Safe design with recessed power button
- Easy snap design

ORDER CODES:
BOVAA00

[Click For More Information](#)

BOVIE CAUTERY LOW-TEMP MICRO FINE TIP BX10

Bovie has been the number one manufacturer of battery-operated cauteries for decades. The safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal.

Bovie cautery low-temp micro fine tip has a safety-enhanced, snap-design & is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. The cautery units are also modified to enable practitioners to “snap off” the end of the cautery for safe and effective removal of the battery power source.



Features:

Portable & easy to use
Safe design with recessed power button
Easy snap design

ORDER CODES:**BOVAA90**[Click For More Information](#)

INDEX

B

Bovie Cautery Low-Temp Fine Tip BX10 2

Bovie Cautery Low-Temp Micro Fine Tip BX10 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au