



Pacific Medical Catalogue

TABLE OF CONTENTS

Disposable Electrodes 2

DISPOSABLE ELECTRODES



BOVIE DISPOSABLE DERMAL ELECTRODE TIPS SHARP NON-STERILE BX100

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVA804

[Click For More Information](#)

BOVIE DISPOSABLE DERMAL ELECTRODE TIPS SHARP STERILE BX50

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVA805

[Click For More Information](#)





BOVIE DISPOSABLE DERMAL ELECTRODE TIPS BLUNT ANGLED NON-STERILE BX100

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVA806**

[Click For More Information](#)

BOVIE DISPOSABLE DERMAL ELECTRODE TIPS BLUNT ANGLED STERILE BX50

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVA807**

[Click For More Information](#)





BOVIE DISPOSABLE DERM-ELITE ROUND ELECTRODE TIPS BLUNT ANGLED NON-STERILE BX100

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These round Derm-Elite electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVA806DE**

[Click For More Information](#)

BOVIE DISPOSABLE DERM-ELITE ROUND ELECTRODE TIPS BLUNT ANGLED STERILE BX50

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These round Derm-Elite electrodes are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVA807DE**

[Click For More Information](#)



BOVIE DISPOSABLE STANDARD BLADE ELECTRODE BX50



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable blade electrode are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES01**

[Click For More Information](#)

BOVIE DISPOSABLE RESTICK II COATED BLADE ELECTRODES BX12

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable Restick II coated blade electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES01T**

[Click For More Information](#)



BOVIE DISPOSABLE EXTENDED BLADE ELECTRODES BX25



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 15.24cm disposable extended blade electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES04**

[Click For More Information](#)

BOVIE DISPOSABLE STANDARD NEEDLE ELECTRODES BX25

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES02**

[Click For More Information](#)



BOVIE DISPOSABLE RESTICK II COATED NEEDLE ELECTRODES BX12



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable Restick II coated needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVES02T

[Click For More Information](#)

BOVIE DISPOSABLE MICRO DISSECTION NEEDLE ELECTRODE BX5

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 4.5cm disposable micro dissection needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVES61

[Click For More Information](#)



BOVIE DISPOSABLE RESTICK II COATED MODIFIED NEEDLE ELECTRODES BX12



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable Restick II coated modified needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES38T**

[Click For More Information](#)

BOVIE DISPOSABLE SUPER CUT MODIFIED FINE NEEDLE ELECTRODES BX5

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 4.5cm disposable super cur modified fine needle electrode are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES63**

[Click For More Information](#)



BOVIE DISPOSABLE IMA / ENT STANDARD BLADE ELECTRODES BX25



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable IMA / ENT standard blade electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVES37

[Click For More Information](#)

BOVIE DISPOSABLE IMA / ENT STANDARD NEEDLE ELECTRODES BX25

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable IMA / ENT standard needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

ORDER CODES:
BOVES38

[Click For More Information](#)



BOVIE DISPOSABLE IMA / ENT EXTENDED NEEDLE ELECTRODES BX25



All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These 6.99cm disposable IMA / ENT extended needle electrodes are suitable for use with all Bovie units including the Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVES40**

[Click For More Information](#)

BOVIE DISPOSABLE H-TYPE ELECTRODE TIPS SHARP ANGLED NON-STERILE BX100

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These H-Type electrodes sharp non-sterile are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVH10012**

[Click For More Information](#)





BOVIE DISPOSABLE H-TYPE ELECTRODE TIPS SHARP STERILE BX50

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These H-Type electrodes sharp sterile are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVH10008**

[Click For More Information](#)

BOVIE DISPOSABLE H-TYPE ELECTRODE TIPS BLUNT ANGLED NON-STERILE BX100

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These H-Type electrodes blunt non-sterile are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVH101112**

[Click For More Information](#)





BOVIE DISPOSABLE H-TYPE ELECTRODE TIPS BLUNT ANGLED STERILE BX50

All electrodes utilize the standard 2.3mm stainless steel shaft electrodes are manufactured to the highest standards and come individually packaged.

These H-Type electrodes blunt sterile are suitable for use with the Bovie Derm 101, Derm 102, Derm 942 (A942) & Bantam PRO (A952). They are also suitable for use with any pencil that accepts the standard 2.3mm stainless steel shaft.

**ORDER CODES:
BOVH10108**

[Click For More Information](#)

INDEX

B

Bovie Disposable Derm-Elite Round Electrode Tips Blunt Angled Non-Sterile BX100 4
Bovie Disposable Derm-Elite Round Electrode Tips Blunt Angled Sterile BX50 4
Bovie Disposable Dermal Electrode Tips Blunt Angled Non-Sterile BX100 3
Bovie Disposable Dermal Electrode Tips Blunt Angled Sterile BX50 3
Bovie Disposable Dermal Electrode Tips Sharp Non-Sterile BX100 2
Bovie Disposable Dermal Electrode Tips Sharp Sterile BX50 2
Bovie Disposable Extended Blade Electrodes BX25 6
Bovie Disposable H-Type Electrode Tips Blunt Angled Non-Sterile BX100 11
Bovie Disposable H-Type Electrode Tips Blunt Angled Sterile BX50 12
Bovie Disposable H-Type Electrode Tips Sharp Angled Non-Sterile BX100 10
Bovie Disposable H-Type Electrode Tips Sharp Sterile BX50 11
Bovie Disposable IMA / ENT Extended Needle Electrodes BX25 10
Bovie Disposable IMA / ENT Standard Blade Electrodes BX25 9
Bovie Disposable IMA / ENT Standard Needle Electrodes BX25 9
Bovie Disposable Micro Dissection Needle Electrode BX5 7
Bovie Disposable Restick II Coated Blade Electrodes BX12 5
Bovie Disposable Restick II Coated Modified Needle Electrodes BX12 8
Bovie Disposable Restick II Coated Needle Electrodes BX12 7
Bovie Disposable Standard Blade Electrode BX50 5
Bovie Disposable Standard Needle Electrodes BX25 6
Bovie Disposable Super Cut Modified Fine Needle Electrodes BX5 8



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au