



Pacific Medical Catalogue

TABLE OF CONTENTS

Smoke Evacuators & Accessories 2

SMOKE EVACUATORS & ACCESSORIES

BOVIE SMOKE SHARK II SMOKE EVACUATOR



The Bovie Smoke Shark II is designed for ease of use and effective smoke and particle removal at the operative site. The Smoke Shark II is designed for use with all electrosurgical generators and effectively controls smoke plume from lasers and other forms of advanced energy.

The Bovie Smoke Shark II is designed for ease of use and effective smoke and particle removal at the operative site. The Bovie Smoke Shark II is designed for use with all electrosurgical generators and effectively controls smoke plume from lasers and other forms of advanced energy.

Features:

- Easy to operate
- Quiet and effective smoke evacuation
- Long life 35 hour filter
- Low cost per procedure
- Excellent particle capture

Specification:

Recommended Accessories

- 786T - Vacuum hose 7/8" x 6"
- 786TS - Vacuum hose 7/8" x 6" hose to hose connect
- SS95 - Smoke evacuation pencil adaptor 3/8" x 20" hose
- SERS2 - Ez link smart remote activator switch

Dimensions

- Height: 15.2cm
- Width: 27.9cm
- Depth: 39.4cm
- Weight: 4.3KG

Line Voltage

220VAC

Line Frequency

50-60Hz

Warranty

2 year warranty

Kit Inclusions

- SE02-220 smoke shark 2 x 1
- FSSEP pneumatic footswitch x 1
- SF35 long life 35 hour filter x 1
- AC power cord x 1
- User guide x 1

ORDER CODES:
BOVSE02-220


MAINTENANCE: Please make sure your device is checked and serviced yearly by an authorised technician.

[Click For More Information](#)

BOVIE SMOKE SHARK III SMOKE EVACUATOR

The Bovie Smoke Shark III is designed for ease of use and effective smoke and particle removal at the operative site. The Smoke Shark III is designed for use with all electrosurgical generators and effectively controls smoke plume.

Bovie smoke shark III reduces biological, physical, and chemical hazards present in electrosurgical smoke plume. Uniquely designed for the physician's office with compact size for wall mounting. Specifically intended for shorter procedures (15 minute or less) with a 4-Stage ULPA filtration 99.999% efficient in trapping 0.1-0.2 micron particulates and micro-organisms. A quiet design offers excellent performance with less noise. Includes a shark skin adapter (SS95) for out of the box usability, a smart shark remove activation device (SERS2) for maximizing filter life & a 6 hour filter (SF06).



Features:

- Compact for the physician office
- Quiet and effective smoke evacuation
- Multi procedure 6 hour filters
- Excellent particle capture

Specification:

Dimensions

- Width: 23.9cm
- Depth: 15.5cm
- Height: 18cm
- Weight: 2.8kg

Line Voltage

220VAC

Line Frequency

50-60Hz

Warranty

2 year warranty

Kit Inclusions

- SE03-220 smoke shark 3 x 1
- SERS2 ezi link remote activator x 1
- SS95 smoke evacuation pencil adaptor x 1
- SF06 6 hour filter x 1
- Wall mounting kit x 1
- AC power cord x 1

User guide x 1

ORDER CODES:
BOVSE03-220

MAINTENANCE: Please make sure your device is checked and serviced yearly by an authorised technician.

[Click For More Information](#)



BOVIE SMOKE SHARK II 35 HOUR FILTER NON-STERILE

Bovie smoke shark II 35 hour filter suitable for use with the Smoke Shark II. Designed for ease of use and effective smoke and particle removal at the operative site.

Bovie smoke shark II 35 hour filter suitable for use with the Smoke Shark II. Designed for ease of use and effective smoke and particle removal at the operative site.

Features:

Easy replacement
35 hour filter life

ORDER CODES:
BOVSF35

[Click For More Information](#)

BOVIE SMOKE SHARK III FILTER NON-STERILE

Bovie smoke shark III 6 hour filter suitable for use with the Smoke Shark III. Designed for ease of use and effective smoke and particle removal at the operative site.

Bovie smoke shark III 6 hour filter suitable for use with the Smoke Shark III. Designed for ease of use and effective smoke and particle removal at the operative site.

Features:

Easy replacement
6 hour filter life

ORDER CODES:
BOVSF06

[Click For More Information](#)



BOVIE SMOKE EVACUATOR PNEUMATIC FOOT SWITCH

Smoke Evacuator Pneumatic Footswitch to suit use with smoke shark devices.

Smoke Evacuator Pneumatic Footswitch to suit use with smoke shark devices.

ORDER CODES:
BOVFSSEP

[Click For More Information](#)

BOVIE VACUUM HOSE 7/8" X 6" BX24

Vacuum hose to hose connect 7/8" x 6" to suit use with smoke shark devices.

Bovie vacuum hose 7/8" x 6" to suit use with smoke shark devices.

ORDER CODES:
BOV786T



[Click For More Information](#)

BOVIE REDUCER FITTER 7/8" TO 1/4" X 24" BX10

Bovie Reducer Fitter 7/8" to 1/4" x 24" to suit use with smoke shark devices.

Bovie Reducer Fitter 7/8" to 1/4" x 24" to suit use with smoke shark devices.

ORDER CODES:
BOVSERF


[Click For More Information](#)

BOVIE SMOKE EVACUATION PENCIL ADAPTOR 3/8" X 20" HOSE BX20

Smoke evacuation pencil adaptor 3/8" X 20" hose to suit use with smoke shark devices.

Smoke evacuation pencil adaptor 3/8" X 20" hose to suit use with smoke shark devices.

**ORDER CODES:
BOVSS95**



[Click For More Information](#)

BOVIE EZ LINK

Bovie ez link remote activator switch to suit use with smoke shark devices.

Bovie ez link remote activator switch to suit use with smoke shark devices. Simply clip onto the handpiece cable to allow automatic activation of the Smoke Shark smoke evacuator.

**ORDER CODES:
BOVSERS2**

[Click For More Information](#)



BOVIE SMOKE EVACUATOR ARTICULATED ARMSTAND

Smoke Evacuator Articulated Armstand to suit use with smoke shark devices.

Smoke evacuator articulated armstand to suit use with smoke shark devices.

ORDER CODES:
BOVFSEAS



[Click For More Information](#)

BOVIE HANDHELD VACUUM HOSE 7/8" X 6" WITH WAND & SPONGE GUARD BX24

Handheld vacuum hose 7/8" x 6" with wand & sponge guard to suit use with smoke shark devices.

Handheld vacuum hose 7/8" x 6" with wand & sponge guard to suit use with smoke shark devices.

ORDER CODES:
BOVSETWS

[Click For More Information](#)



**BOVIE VACUUM HOSE 7/8" X 6" HOSE TO HOSE
CONNECT BX24**

Vacuum hose to hose connect 7/8" x 6" to suit use with smoke shark devices.

Vacuum hose to hose connect 7/8" x 6" to suit use with smoke shark devices.

ORDER CODES:
BOV786TS

[Click For More Information](#)



INDEX

B

Bovie Ez Link	7
Bovie Handheld Vacuum Hose 7/8" x 6" With Wand & Sponge Guard BX24	8
Bovie Reducer Fitter 7/8" to 1/4" x 24" BX10	6
Bovie Smoke Evacuation Pencil Adaptor 3/8" X 20" Hose BX20	7
Bovie Smoke Evacuator Articulated Armstand	8
Bovie Smoke Evacuator Pneumatic Foot Switch	5
Bovie Smoke Shark II 35 Hour Filter Non-Sterile	4
Bovie Smoke Shark II Smoke Evacuator	2
Bovie Smoke Shark III Filter Non-Sterile	5
Bovie Smoke Shark III Smoke Evacuator	3
Bovie Vacuum Hose 7/8" x 6" BX24	6
Bovie Vacuum Hose 7/8" x 6" Hose To Hose Connect BX24	9



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au