



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Hospital Beds** ..... 2

# HOSPITAL BEDS

## ES SERIES HOSPITAL BED WITH X-RAY BACKREST (COMING SOON)



A fully adjustable electric hospital bed specifically designed for standard hospital wards and long-term care units. The bed has a simple structure, a wide range of functions beyond the standard, and complex functions, providing patients and caregivers with cost effectiveness and ease of use. Uniquely equipped with an X-ray translucent backrest.

### ES Series Hospital Bed Features:

- Five function electric operation of backrest, knee break, height and trendelenburg
- Integrated side rail controls for patient & caregiver
- Integral function lockout
- X-Ray translucent backrest
- Electric & manual CPR release on both sides
- Removeable head & foot board
- Integrated mattress platform extension
- Undercarriage light
- Supports fowlers chair position
- Visual indicator for minimum height position
- Battery backup
- Emergency stop button
- Corner bumpers with self help & iv pole locking mechanism
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm central castors with directional lock
- Foam mattress

### ES Series Hospital Bed Specification:

- Overall length: 228cm
- Overall width: 105cm
- Height range: 40.5-75cm
- Hoist clearance: 15cm
- Mattress platform length: 200cm
- Mattress platform width: 90cm
- Mattress platform extension: 10cm
- Mattress size: 200x90x10cm
- Backrest angle: 0-70 degrees
- One touch backrest positions: 30 degrees then continuous to 70 degrees
- Knee break angle: 0-27 degrees
- Leg angle: 0-12 degrees
- Trendelenburg & reverse trendelenburg: 0-13 degrees
- SWL: 300KG

Warranty: 3 Years  
Certification: CE, FDA, ISO13485-2016

**ORDER CODES:**  
**HBES407EX**

**MAINTENANCE:** Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.

[Click For More Information](#)

## **ES SERIES HOSPITAL BED WITH X-RAY BACKREST & WEIGHT SCALE (COMING SOON)**

A fully adjustable electric hospital bed specifically designed for standard hospital wards and long-term care units. The bed has a simple structure, a wide range of functions beyond the standard, and complex functions, providing patients and caregivers with cost effectiveness and ease of use. Uniquely equipped with an X-ray translucent backrest & integrated weight scale.



### **ES Series Hospital Bed Features:**

- Five function electric operation of backrest, knee break, height and trendelenburg
- Integrated side rail controls for patient & caregiver
- Integral function lockout
- X-Ray translucent backrest
- Integrated weight scale
- Bed exit alarm
- Electric & manual CPR release on both sides
- Removeable head & foot board
- Integrated mattress platform extension
- Undercarriage light
- Supports fowlers chair position
- Visual indicator for minimum height position
- Battery backup
- Emergency stop button
- Corner bumpers with self help & iv pole locking mechanism
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm central castors with directional lock
- Foam mattress

### **ES Series Hospital Bed Specification:**

- Overall length: 228cm
- Overall width: 105cm
- Height range: 40.5-75cm
- Hoist clearance: 15cm
- Mattress platform length: 200cm

Mattress platform width: 90cm  
 Mattress platform extension: 10cm  
 Mattress size: 200x90x10cm  
 Backrest angle: 0-70 degrees  
 One touch backrest positions: 30 degrees then continuous to 70 degrees  
 Knee break angle: 0-27 degrees  
 Leg angle: 0-12 degrees  
 Trendelenburg & reverse trendelenburg: 0-13 degrees  
 Weight scale capacity: 350KG  
 Weight scale graduation: 0.1KG  
 SWL: 300KG  
 Warranty: 3 Years  
 Certification: CE, FDA, ISO13485-2016

**ORDER CODES:  
HBES408EXW**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)

## THREE FUNCTION HOSPITAL BED



A strong, reliable and high quality three function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for aged care and hospital facilities.

### Three Function Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Removeable head & foot board
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### Three Function Hospital Bed Specification:

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm
- Hoist clearance: 16cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Mattress size: 195x90x10cm
- Backrest angle: 0-73 degrees
- Knee break angle: 0-30 degrees

Leg angle: 0-16 degrees  
SWL: 225KG  
Warranty: 3 Years  
Certification: CE, FDA, ISO13485-2016

**ORDER CODES:  
HB3**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)

## FIVE FUNCTION HOSPITAL BED

A strong, reliable and high quality five function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for higher needs patients and longer term patient care.



### Five Function Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### Five Function Hospital Bed Specification:

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm
- Hoist clearance: 16cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Mattress size: 195x90x10cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-16 degrees
- Trendelenburg & reverse trendelenburg: 0-13 degrees
- SWL: 300KG
- Warranty: 3 Years
- Certification: CE, FDA, ISO13485-2016

### ORDER CODES:

**HB5**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)



## **TILTING ICU HOSPITAL BED**

A strong, reliable and high quality five function tilting ICU hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for high needs patients requiring long term treatment with the tilting function allowing treatment of the nervous system with weakness, pain, muscle spasm, and for auxiliary rehabilitation purposes such as hemiplegia, paraplegia and traumatic brain injury and other severe patients to resume standing training

### **Tilting ICU Hospital Bed Features:**

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Tilt inclinometer
- Emergency stop button
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Durable patient belts
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### **Tilting ICU Hospital Bed Specification:**

- Overall length: 215cm
- Overall width: 110cm
- Height range: 52-76cm
- Hoist clearance: 14cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-20 degrees
- Trendelenburg & reverse trendelenburg: 0-16 degrees & 0-13 degrees
- Tilting angle: 0-85 degrees
- SWL: 225KG
- Warranty: 3 Years
- Certification: CE, FDA, ISO13485-2016

**ORDER CODES:  
HBT**

**Please note these are an indent item and are manufactured to order.**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)

# INDEX

## **E**

ES Series Hospital Bed With X-Ray Backrest & Weight Scale (COMING SOON) 3

ES Series Hospital Bed With X-Ray Backrest (COMING SOON) 2

## **F**

Five Function Hospital Bed 5

## **T**

Three Function Hospital Bed 4

Tilting ICU Hospital Bed 6



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)