



Pacific Medical Catalogue

TABLE OF CONTENTS

Infection Control	2
Autoclaves	2
Self Seal Pouches	4
Ultrasonic Cleaners	8

INFECTION CONTROL

AUTOCLAVES



AURA TOUCHSCREEN AUTOCLAVE 17 / 22L - CLASS B & S

The new gold standard in benchtop sterilisation. The AURA Autoclave is loaded with features and European componentry, it will reliably serve any clinic time and time again.

Features:

- Reliable European Componentry
- Sleek Finish
- Easy To Operate User Interface
- Touch Screen & Buttons Allowing Fast Operation
- Inbuilt Printer & USB Functions
- Microprocessor Control System For Precise Sensor Management
- Automated Cleaning Program & Maintenance Reminders
- Automatic RO & Manual Fill Options
- S & B Class Cycles
- Delayed Start Function
- Automatic Component Diagnosis Tool
- Water Quality Sensor
- Five Selectable Languages

Specification:

- External Dimensions: 625 x 440 x 475mm
- Chamber Dimensions: 249 x 355mm (17L) / 249 x 450mm (22L)
- Weight: 52KG (17L) / 55KG (22L)
- Selectable 134°C Cycles: Universal Class B, Fast Class S, Prion, Custom
- Selectable 121°C Cycles: Universal Class B, Fast Class S, Custom
- Test Cycles: B&D Helix, Leak & Vacuum
- Jacket Heater With Six Different Settings
- Drying Settings Available: Normal, Extra Dry 5 min, Extra Dry 10 min
- Included Accessories: Tray Rack, 3 X Trays, 2 X Drain Hoses & Connectors, Tray Handle, Door Seal, Chamber & Water Tank Filter Mesh, Paper Roll, USB Drive, AC Power Cord
- Required Voltage: 120 - 230V AC
- Rated Frequency: 50 / 60 Hz
- Rated Power: 1800 VA
- Sterilising Temperature / Pressure: 121°C @ 100-130 Kpa, 134°C @ 200-240 Kpa
- Water Tank Capacity: 3.5L

Operating Temperature: +5°C - +40°C

Relative Humidity Range: 30% - 75%

Warranty: 2 Years

Certifications: CE, ISO13485-2016, ASME / Worksafe, EN13060

ORDER CODES:

AURA17 - 17L

AURA22 - 22L

NOTE: Professional installation & validation required. Please make sure your autoclave is serviced in accordance with the users manual and machine prompts to ensure it functions to its optimum. All servicing must be completed by our network of authorised agents to maintain machine warranty.

[Click For More Information](#)

SELF SEAL POUCHES

SELF SEALING AUTOCLAVE POUCH 57X135MM BX200



A simple high quality self seal autoclave pouch 57x135mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE: SSP57x135

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 70X260MM BX200

A simple high quality self seal autoclave pouch 70x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.



Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE: SSP70x260

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 90X160MM BX200



A simple high quality self seal autoclave pouch 90x160mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:
SSP90x160**

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 90X260MM BX200

A simple high quality self seal autoclave pouch 90x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.



Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:
SSP90x260**

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 135X285MM BX200



A simple high quality self seal autoclave pouch 135x285mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:
SSP135x285

[Click For More Information](#)

SELF SEALING AUTOCLAVE POUCH 190X360MM BX200

A simple high quality self seal autoclave pouch 190x360mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

ORDER CODE:
SSP190x360

[Click For More Information](#)





SELF SEALING AUTOCLAVE POUCH 300X380MM BX200

A simple high quality self seal autoclave pouch 300x380mm in size. It features internal and external sterilisation indicators and a strong triple seal.

Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:
SSP300x380**

[Click For More Information](#)

ULTRASONIC CLEANERS

ULTRASONIC CLEANER 2 LITRE



Our 2L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 2L
- Tank Dimensions: 150 x 140 x 100mm (LxWxH)
- External Dimensions: 175 x 165 x 230mm (LxWxH)
- Basket Dimensions: 138 x 128 x 73mm (LxWxH)
- Heating power: 100W
- Ultrasonic power: 50W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD2

[Click For More Information](#)

ULTRASONIC CLEANER 3 LITRE

Our 3L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 3L
- Tank Dimensions: 240 x 140 x 100mm (LxWxH)
- External Dimensions: 330 x 165 x 230mm (LxWxH)
- Basket Internal Dimensions: 229 x 125 x 70mm(LxWxH)
- Heating power: 100W
- Ultrasonic power: 100W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD3

[Click For More Information](#)



ULTRASONIC CLEANER 6 LITRE

Our 6L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:



Capacity: 6L
Tank Dimensions: 300 x 150 x 150mm (LxWxH)
External Dimensions: 400 x 175 x 305mm (LxWxH)
Basket Internal Dimensions: 265 x 116 x 110mm (LxWxH)
Heating power: 300W
Ultrasonic power: 150W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:

USD6

[Click For More Information](#)

ULTRASONIC CLEANER 9 LITRE

Our 9L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

Capacity: 9L
Tank Dimensions: 300 x 235 x 150mm (LxWxH)
External Dimensions: 400 x 265 x 305mm (LxWxH)
Basket Internal Dimensions: 264 x 204 x 96mm (LxWxH)
Heating power: 300W
Ultrasonic power: 200W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:

USD9

[Click For More Information](#)





ULTRASONIC CLEANER 13 LITRE

Our 13L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 13L
- Tank Dimensions: 330 x 300 x 150mm (LxWxH)
- External Dimensions: 360 x 330 x 310mm (LxWxH)
- Basket Dimensions: 294 x 264 x 96mm (LxWxH)
- Heating power: 300W
- Ultrasonic power: 400W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: Yes
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD13

[Click For More Information](#)

ULTRASONIC CLEANER 20 LITRE

Our 20L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 20L
- Tank Dimensions: 500 x 300 x 150mm (LxWxH)
- External Dimensions: 550 x 330 x 310mm (LxWxH)
- Basket Internal Dimensions: 458 x 263 x 96mm (LxWxH)
- Heating power: 500W
- Ultrasonic power: 400W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: Yes
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD20

[Click For More Information](#)



ULTRASONIC CLEANER 27 LITRE

Our 27L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation



Comes with basket and lid as standard

Specifications:

Capacity: 27L

Tank Dimensions: 500 x 300 x 200mm (LxWxH)

External Dimensions: 550 x 330 x 360mm (LxWxH)

Basket Internal Dimensions: 458 x 263 x 130mm (LxWxH)

Heating power: 500W

Ultrasonic power: 500W

Ultrasonic frequency: 40KHz

Temp Adjustable: Ambient to 80°C

Drainage Valve: Yes

Display: Digital

Timer: 1 - 99min

ORDER CODE:

USD27

[Click For More Information](#)

INDEX

A

AURA Touchscreen Autoclave 17 / 22L - Class B & S 2

S

Self Sealing Autoclave Pouch 135x285mm BX200 6

Self Sealing Autoclave Pouch 190x360mm BX200 6

Self Sealing Autoclave Pouch 300x380mm BX200 7

Self Sealing Autoclave Pouch 57x135mm BX200 4

Self Sealing Autoclave Pouch 70x260mm BX200 4

Self Sealing Autoclave Pouch 90x160mm BX200 5

Self Sealing Autoclave Pouch 90x260mm BX200 5

U

Ultrasonic Cleaner 13 Litre 11

Ultrasonic Cleaner 2 Litre 8

Ultrasonic Cleaner 20 Litre 12

Ultrasonic Cleaner 27 Litre 12

Ultrasonic Cleaner 3 Litre 9

Ultrasonic Cleaner 6 Litre 9

Ultrasonic Cleaner 9 Litre 10



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au