



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Autoclaves** ..... 2

# AUTOCLAVES



## AURA TOUCHSCREEN AUTOCLAVE 17 / 22L - CLASS B & S

The new gold standard in benchtop sterilisation. The AURA Autoclave is loaded with features and European componentry, it will reliably serve any clinic time and time again.

### Features:

- Reliable European Componentry
- Sleek Finish
- Easy To Operate User Interface
- Touch Screen & Buttons Allowing Fast Operation
- Inbuilt Printer & USB Functions
- Microprocessor Control System For Precise Sensor Management
- Automated Cleaning Program & Maintenance Reminders
- Automatic RO & Manual Fill Options
- S & B Class Cycles
- Delayed Start Function
- Automatic Component Diagnosis Tool
- Water Quality Sensor
- Five Selectable Languages

### Specification:

- External Dimensions: 625 x 440 x 475mm
- Chamber Dimensions: 249 x 355mm (17L) / 249 x 450mm (22L)
- Weight: 52KG (17L) / 55KG (22L)
- Selectable 134°C Cycles: Universal Class B, Fast Class S, Prion, Custom
- Selectable 121°C Cycles: Universal Class B, Fast Class S, Custom
- Test Cycles: B&D Helix, Leak & Vacuum
- Jacket Heater With Six Different Settings
- Drying Settings Available: Normal, Extra Dry 5 min, Extra Dry 10 min
- Included Accessories: Tray Rack, 3 X Trays, 2 X Drain Hoses & Connectors, Tray Handle, Door Seal, Chamber & Water Tank Filter Mesh, Paper Roll, USB Drive, AC Power Cord
- Required Voltage: 120 - 230V AC
- Rated Frequency: 50 / 60 Hz
- Rated Power: 1800 VA
- Sterilising Temperature / Pressure: 121°C @ 100-130 Kpa, 134°C @ 200-240 Kpa
- Water Tank Capacity: 3.5L
- Operating Temperature: +5°C - +40°C
- Relative Humidity Range: 30% - 75%
- Warranty: 2 Years
- Certifications: CE, ISO13485-2016, ASME / Worksafe, EN13060

### ORDER CODES:

**AURA17 - 17L**

**AURA22 - 22L**

**NOTE: Professional installation & validation required. Please make sure your autoclave is serviced in accordance with the users manual and machine prompts to ensure it functions to its optimum. All servicing must be completed by our network of authorised agents to maintain machine warranty.**

[Click For More Information](#)

# INDEX

## A

AURA Touchscreen Autoclave 17 / 22L - Class B & S 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)