



Pacific Medical Catalogue

TABLE OF CONTENTS

Ultrasonic Cleaners 2

ULTRASONIC CLEANERS

ULTRASONIC CLEANER 2 LITRE



Our 2L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 2L
- Tank Dimensions: 150 x 140 x 100mm (LxWxH)
- External Dimensions: 175 x 165 x 230mm (LxWxH)
- Basket Dimensions: 138 x 128 x 73mm (LxWxH)
- Heating power: 100W
- Ultrasonic power: 50W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD2

[Click For More Information](#)

ULTRASONIC CLEANER 3 LITRE

Our 3L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 3L
- Tank Dimensions: 240 x 140 x 100mm (LxWxH)
- External Dimensions: 330 x 165 x 230mm (LxWxH)
- Basket Internal Dimensions: 229 x 125 x 70mm(LxWxH)
- Heating power: 100W
- Ultrasonic power: 100W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: No
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD3

[Click For More Information](#)



ULTRASONIC CLEANER 6 LITRE

Our 6L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:



Capacity: 6L
Tank Dimensions: 300 x 150 x 150mm (LxWxH)
External Dimensions: 400 x 175 x 305mm (LxWxH)
Basket Internal Dimensions: 265 x 116 x 110mm (LxWxH)
Heating power: 300W
Ultrasonic power: 150W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD6**[Click For More Information](#)

ULTRASONIC CLEANER 9 LITRE

Our 9L ultrasonic cleaner with a sleek black titanium finish and plexiglass operation panel are sleek, modern & functional.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

Capacity: 9L
Tank Dimensions: 300 x 235 x 150mm (LxWxH)
External Dimensions: 400 x 265 x 305mm (LxWxH)
Basket Internal Dimensions: 264 x 204 x 96mm (LxWxH)
Heating power: 300W
Ultrasonic power: 200W
Ultrasonic frequency: 40KHz
Temp Adjustable: Ambient to 80°C
Drainage Valve: Yes
Display: Digital
Timer: 1 - 99min

ORDER CODE:**USD9**[Click For More Information](#)



ULTRASONIC CLEANER 13 LITRE

Our 13L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 13L
- Tank Dimensions: 330 x 300 x 150mm (LxWxH)
- External Dimensions: 360 x 330 x 310mm (LxWxH)
- Basket Dimensions: 294 x 264 x 96mm (LxWxH)
- Heating power: 300W
- Ultrasonic power: 400W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: Yes
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD13

[Click For More Information](#)

ULTRASONIC CLEANER 20 LITRE

Our 20L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation
- Comes with basket and lid as standard

Specifications:

- Capacity: 20L
- Tank Dimensions: 500 x 300 x 150mm (LxWxH)
- External Dimensions: 550 x 330 x 310mm (LxWxH)
- Basket Internal Dimensions: 458 x 263 x 96mm (LxWxH)
- Heating power: 500W
- Ultrasonic power: 400W
- Ultrasonic frequency: 40KHz
- Temp Adjustable: Ambient to 80°C
- Drainage Valve: Yes
- Display: Digital
- Timer: 1 - 99min

ORDER CODE:

USD20

[Click For More Information](#)



ULTRASONIC CLEANER 27 LITRE

Our 27L ultrasonic cleaners are quiet with variable rinse times, temperatures, stainless steel finish and a clear digital display they are fitting for any modern clinic.

Features:

- Turbo power - normal and soft functions
- Auto-degas function
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel for easy operation



Comes with basket and lid as standard

Specifications:

Capacity: 27L

Tank Dimensions: 500 x 300 x 200mm (LxWxH)

External Dimensions: 550 x 330 x 360mm (LxWxH)

Basket Internal Dimensions: 458 x 263 x 130mm (LxWxH)

Heating power: 500W

Ultrasonic power: 500W

Ultrasonic frequency: 40KHz

Temp Adjustable: Ambient to 80°C

Drainage Valve: Yes

Display: Digital

Timer: 1 - 99min

ORDER CODE:**USD27**

[Click For More Information](#)

INDEX

U

Ultrasonic Cleaner 13 Litre	5
Ultrasonic Cleaner 2 Litre	2
Ultrasonic Cleaner 20 Litre	6
Ultrasonic Cleaner 27 Litre	6
Ultrasonic Cleaner 3 Litre	3
Ultrasonic Cleaner 6 Litre	3
Ultrasonic Cleaner 9 Litre	4



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au