



Pacific Medical Catalogue

TABLE OF CONTENTS

Lift Chairs & Recliners 2

LIFT CHAIRS & RECLINERS



PACIFIC LIFT CHAIR / RECLINER

The Pacific Lift Chair / Recliner is a well featured dual motor chair with a generous 170KG SWL, removable armrests and 180 degree recline to accommodate most requirements.

Features:

- Four Castors With Two Lockable
- Two motor with lift function
- Reclining functions up to 180 degrees
- Light up remote control with USB charging port
- Quiet operation
- Removable armrests
- Removable wing headrests
- High quality charcoal grey vinyl
- Reliable TiMotion motors

Specification:

- Seat Width: 56cm
- Seat Depth: 54cm
- Seat Height: 53cm
- Overall Width (At Armrest): 74cm
- Overall Height: 109cm
- Total Length After Reclining: 170cm
- SWL: 170kg

ORDER CODE:

REC - Charcoal Grey

MAINTENANCE: Please make sure your chair is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

INDEX

P

Pacific Lift Chair / Recliner 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au