



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Liquid Nitrogen Dewars ..... 2**

# LIQUID NITROGEN DEWAR



## 10L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 10L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

- Made of a high-strength aluminium alloy
- Advanced powder coating technology
- Designed for low liquid nitrogen consumption
- 5 year vacuum warranty
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

- Capacity: 10L
- Weight (empty): 6.0kg
- Caliber (neck opening of dewar): 50mm
- Outer diameter of dewar: 300mm
- Height of dewar: 552mm
- Static evaporation rate: 0.12L/d
- Static holding time: 86 days
- Ladle outer diameter: 38mm
- Ladle height: 120mm
- No. ladles per dewar: 3~6
- Warranty: 5 Years Vacuum, 2 Years Container

- Straws holding capacity:
  - Single 120mm (0.5mL): 792pc
  - Single 120mm (0.25mL): 1788pc

### ORDER CODE: LNC10

[Click For More Information](#)

## 20L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 20L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

- Made of a high-strength aluminium alloy
- Advanced powder coating technology
- Designed for low liquid nitrogen consumption
- 5 year vacuum warranty
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

- Capacity: 20L
- Weight (empty): 9.5kg
- Caliber (neck opening of dewar): 50mm
- Outer diameter of dewar: 394mm
- Height of dewar: 672mm
- Static evaporation rate: 0.12L/d
- Static holding time: 168 days
- Ladle outer diameter: 38mm
- Ladle height: 120/276mm
- No. ladles per dewar: 3~6

- Straws holding capacity:
  - Single 120mm (0.5mL): 792pc
  - Single 120mm (0.25mL): 1788pc
  - Double 276mm (0.5mL): 1284
  - Double 276mm (0.25mL): 2832

### ORDER CODE:

**LNC20**

[Click For More Information](#)





## 30L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 30L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

- Made of a high-strength aluminium alloy
- Advanced powder coating technology
- Designed for low liquid nitrogen consumption
- 5 year vacuum warranty
- Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

- Capacity: 30L
- Weight (empty): 12.9kg
- Caliber (neck opening of dewar): 50mm
- Outer diameter of dewar: 462mm
- Height of dewar: 706mm
- Static evaporation rate: 0.12L/d
- Static holding time: 254 days
- Ladle outer diameter: 38mm
- Ladle height: 120/276mm
- No. ladles per dewar: 3~6

- Straws holding capacity:
  - Single 120mm (0.5mL): 792pc
  - Single 120mm (0.25mL): 1788pc
  - Double 276mm (0.5mL): 1284
  - Double 276mm (0.25mL): 2832

### ORDER CODE: LNC30

[Click For More Information](#)

## CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

### Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

### Order Code:

**LNPACK1**

[Click For More Information](#)





## CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

### Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

### Order Code:

**LNPACK2**

[Click For More Information](#)

# INDEX

## **1**

10L Liquid Nitrogen Dewar 2

## **2**

20L Liquid Nitrogen Dewar 3

## **3**

30L Liquid Nitrogen Dewar 4

## **C**

Cryotherapy Package One 5

Cryotherapy Package Two 6





**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)