



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Liquid Nitrogen Dewars** ..... 2

# LIQUID NITROGEN DEWAR



## 10L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 10L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

Made of a high-strength aluminium alloy  
Advanced powder coating technology  
Designed for low liquid nitrogen consumption  
5 year vacuum warranty  
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

Capacity: 10L  
Weight (empty): 6.0kg  
Caliber (neck opening of dewar): 50mm  
Outer diameter of dewar: 300mm  
Height of dewar: 552mm  
Static evaporation rate: 0.12L/d  
Static holding time: 86 days  
Ladle outer diameter: 38mm  
Ladle height: 120mm  
No. ladles per dewar: 3~6  
Warranty: 5 Years Vacuum, 2 Years Container

Straws holding capacity:  
Single 120mm (0.5mL): 792pc  
Single 120mm (0.25mL): 1788pc

### ORDER CODE: LNC10

[Click For More Information](#)

## 20L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 20L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

Made of a high-strength aluminium alloy  
Advanced powder coating technology  
Designed for low liquid nitrogen consumption  
5 year vacuum warranty  
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

Capacity: 20L  
Weight (empty): 9.5kg  
Caliber (neck opening of dewar): 50mm  
Outer diameter of dewar: 394mm  
Height of dewar: 672mm  
Static evaporation rate: 0.12L/d  
Static holding time: 168 days  
Ladle outer diameter: 38mm  
Ladle height: 120/276mm  
No. ladles per dewar: 3~6

Straws holding capacity:  
Single 120mm (0.5mL): 792pc  
Single 120mm (0.25mL): 1788pc  
Double 276mm (0.5mL): 1284  
Double 276mm (0.25mL): 2832

### ORDER CODE:

**LNC20**

[Click For More Information](#)





## 30L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 30L dewar with a safety lock lid and low liquid nitrogen consumption.

### Features:

Made of a high-strength aluminium alloy  
Advanced powder coating technology  
Designed for low liquid nitrogen consumption  
5 year vacuum warranty  
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

### Specifications:

Capacity: 30L  
Weight (empty): 12.9kg  
Caliber (neck opening of dewar): 50mm  
Outer diameter of dewar: 462mm  
Height of dewar: 706mm  
Static evaporation rate: 0.12L/d  
Static holding time: 254 days  
Ladle outer diameter: 38mm  
Ladle height: 120/276mm  
No. ladles per dewar: 3~6

Straws holding capacity:  
Single 120mm (0.5mL): 792pc  
Single 120mm (0.25mL): 1788pc  
Double 276mm (0.5mL): 1284  
Double 276mm (0.25mL): 2832

**ORDER CODE:**  
**LNC30**

[Click For More Information](#)

## CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

### Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

### Order Code:

**LNPACK1**

[Click For More Information](#)





## CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

### Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

### Order Code:

**LNPACK2**

[Click For More Information](#)

# INDEX

## **1**

10L Liquid Nitrogen Dewar 2

## **2**

20L Liquid Nitrogen Dewar 3

## **3**

30L Liquid Nitrogen Dewar 4

## **C**

Cryotherapy Package One 5

Cryotherapy Package Two 6





**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)