



Pacific Medical Catalogue

TABLE OF CONTENTS

Liquid Nitrogen Flask 2

LIQUID NITROGEN FLASK



NITROSPRAY CRYOTHERAPY FLASK

Premier Nitro Spray cryotherapy flask constructed of a double-walled vacuum stainless steel container to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage & disengage spray tips and other attachments as required.

Features:

- Double walled vacuum cryotherapy flask
- Easy connect luer lock connector
- Relief valve precisely regulates pressure
- Curved trigger provides unobstructed view of treatment site
- Easy to hold grooved head
- Durable stainless steel construction
- Spray tip protector prevents accidental damage and contact with patient

Specification:

- Holding Time: 6-8 Hours
- Intermittent Use: 3-6 Hours Treatment Time
- Package Inclusions: Flask, Set of five spray tips & Spray tip protector
- Capacity: 300 (Nitro Spray Lite) & 475ml (Nitro Spray Plus)
- Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

ORDER CODES:

LNSLITEPLUS - 300ml

LNSPLUS - 475ml

[Click For More Information](#)

CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK1

[Click For More Information](#)





CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK2

[Click For More Information](#)

NITROSPRAY CRYOTHERAPY SPRAY TIPS

Spray tips suitable for the Nitrospray cryotherapy flasks.

Features:

Suitable for use with the Nitrospray 300ml & 475ml flasks

Easy connect luer lock connector

Comes in 16, 17, 18, 19 & 20 gauge sizes

Available individually or as a set of five

ORDER CODES:

LNST16G - 16 Gauge

LNST17G - 17 Gauge

LNST18G - 18 Gauge

LNST19G - 19 Gauge



LNST20G - 20 Gauge**LNSTSET - 16, 17, 18, 19 & 20 Gauge**[Click For More Information](#)

NITROSPRAY SPRAY TIP PROTECTOR

A protector for the spray tip of the Nitrospray cryotherapy flask.

A protector for the spray tip of the Nitrospray cryotherapy flask.

Order Code:

LNSTP

[Click For More Information](#)

NITROSPRAY ACCESSORY ADAPTOR FOR SCREW TIPS

An adaptor that allows the connection of screw tips from brands such as Brymill.

An adaptor that has a luer lock connection on one end and a screw tip on the other. This allows the connection of spray tips from brands such as Brymill.

Order Code:

LNSABA

[Click For More Information](#)

NITROSPRAY SPRAY TIP EXTENSION BENT



A long curved needle spray tip which connects to the luer lock connection.

A long curved needle spray tip which connects to the luer lock connection. Allows better patient access when treating difficult areas.

Order Codes:

LNSTEB

[Click For More Information](#)

NITROSPRAY NEEDLE EXTENDER RIGID

A rigid extender for the Nitrospray spray tips.

A rigid steel extender for the existing luer lock connection. Allows the tip better reach and patient access when treating difficult areas.



Order Code:

LNSNER

[Click For More Information](#)

INDEX

C

Cryotherapy Package One 3

Cryotherapy Package Two 4

N

Nitrospray Accessory Adaptor For Screw Tips 5

Nitrospray Cryotherapy Flask 2

Nitrospray Cryotherapy Spray Tips 4

Nitrospray Needle Extender Rigid 6

Nitrospray Spray Tip Extension Bent 6

Nitrospray Spray Tip Protector 5



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au