



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Liquid Nitrogen Flask** ..... 2

# LIQUID NITROGEN FLASK



## NITROSPRAY CRYOTHERAPY FLASK

Premier Nitrospray cryotherapy flask constructed of a double-walled vacuum stainless steel container to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage & disengage spray tips and other attachments as required.

### Features:

- Double walled vacuum cryotherapy flask
- Easy connect luer lock connector
- Relief valve precisely regulates pressure
- Curved trigger provides unobstructed view of treatment site
- Easy to hold grooved head
- Durable stainless steel construction
- Spray tip protector prevents accidental damage and contact with patient

### Specification:

- Holding Time: 6-8 Hours
- Intermittent Use: 3-6 Hours Treatment Time
- Package Inclusions: Flask, Set of five spray tips & Spray tip protector
- Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)
- Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

### ORDER CODES:

**LNSLITEPLUS - 300ml**

**LNSPLUS - 475ml**

[Click For More Information](#)

## CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



### Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

### Order Code:

**LNPACK1**

[Click For More Information](#)

## CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

### Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



**Dewar Specification:**

Capacity: 20L  
Weight (empty): 9.5kg  
Caliber (neck opening of dewar): 50mm  
Outer diameter of dewar: 394mm  
Height of dewar: 672mm  
Static evaporation rate: 0.12L/d  
Static holding time: 168 days  
Ladle outer diameter: 38mm  
Ladle height: 120/276mm  
No. ladles per dewar: 3~6

**Order Code:****LNPACK2**[Click For More Information](#)**NITROSPRAY CRYOTHERAPY SPRAY TIPS**

Spray tips suitable for the Nitrospray cryotherapy flasks.

**Features:**

Suitable for use with the Nitrospray 300ml & 475ml flasks  
Easy connect luer lock connector  
Comes in 16, 17, 18, 19 & 20 gauge sizes  
Available individually or as a set of five

**ORDER CODES:****LNST16G - 16 Gauge****LNST17G - 17 Gauge****LNST18G - 18 Gauge****LNST19G - 19 Gauge****LNST20G - 20 Gauge****LNSTSET - 16, 17, 18, 19 & 20 Gauge**[Click For More Information](#)

## NITROSPRAY SPRAY TIP PROTECTOR

A protector for the spray tip of the Nitrospray cryotherapy flask.

A protector for the spray tip of the Nitrospray cryotherapy flask.

**Order Code:**

**LNSTP**

[Click For More Information](#)



## NITROSPRAY ACCESSORY ADAPTOR FOR SCREW TIPS

An adaptor that allows the connection of screw tips from brands such as Brymill.

An adaptor that has a luer lock connection on one end and a screw tip on the other. This allows the connection of spray tips from brands such as Brymill.

**Order Code:**

**LNSABA**

[Click For More Information](#)



## NITROSPRAY SPRAY TIP EXTENSION BENT



A long curved needle spray tip which connects to the luer lock connection.

A long curved needle spray tip which connects to the luer lock connection. Allows better patient access when treating difficult areas.

**Order Codes:**

**LNSTEB**

[Click For More Information](#)

## NITROSPRAY NEEDLE EXTENDER RIGID

A rigid extender for the Nitrospray spray tips.

A rigid steel extender for the existing luer lock connection. Allows the tip better reach and patient access when treating difficult areas.



**Order Code:**

**LNSNER**

[Click For More Information](#)

# INDEX

## C

Cryotherapy Package One 3

Cryotherapy Package Two 3

## N

Nitrospray Accessory Adaptor For Screw Tips 5

Nitrospray Cryotherapy Flask 2

Nitrospray Cryotherapy Spray Tips 4

Nitrospray Needle Extender Rigid 6

Nitrospray Spray Tip Extension Bent 6

Nitrospray Spray Tip Protector 5





**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)