



# Pacific Medical Catalogue

## TABLE OF CONTENTS

<b>Medical Carts</b> .....	2
Medical Carts .....	2
Computer Trolleys .....	11

# MEDICAL CARTS

## MEDICAL CARTS



### EMERGENCY TROLLEY

A rugged portable emergency trolley constructed with steel & ABS. Along with a heavy load rating and full set of accessories included it's a best seller.

**Note: Product pictured with central & seals lock. It can only come with one or the other.**

#### Features:

- Five drawer
- Extension table
- CPR board
- Defibrillator Shelf
- Stainless steel IV pole
- Waste bin
- 304 stainless steel wire basket
- Oxygen tank holder
- External power holder
- ABS drawer dividers
- Central or seals lock
- Noise resistant castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Smooth welding

#### Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-3 100mm, 4 - 125mm, 5 - 190mm
- Drawer load capacity: 30kg
- Overall load capacity: 200kg
- Weight: 55kg

All accessories included

#### ORDER CODE:

**EMT - Emergency Trolley With Central Lock**  
**EMTSL - Emergency Trolley With Seals Lock**

[Click For More Information](#)

## MEDICINE TROLLEY

For the safe transport of medications and supplies our four drawer medicine trolley is the perfect fit. Including all accessories and boasting a number of features this product represents great value.

### Features:

- Five drawers
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

### Specifications:

- Dimensions: 650 x 480 x 1070mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-4 110mm, 5 - 250mm
- Top loading: 75kg
- Weight: 42kg

All accessories included

**ORDER CODE:**  
**MEDT**

[Click For More Information](#)





## TREATMENT TROLLEY ONE DRAWER

This single drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

### Features:

- One drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

### Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 100mm
- Top loading: 75kg
- Weight: 35kg

All accessories included

### ORDER CODE: TTR01

[Click For More Information](#)

## TREATMENT TROLLEY TWO DRAWER

This two drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

### Features:

- Two drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer



**Specifications:**

Dimensions: 650 x 480 x 1000mm  
Drawer Depth: 100mm  
Top loading: 75kg  
Weight: 35kg

All accessories included

**ORDER CODE:  
TTR02**

[Click For More Information](#)

**TREATMENT TROLLEY THREE DRAWER**

This three drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

**Features:**

Three drawer  
Bottom shelf  
Extension table  
Waste bin  
304 stainless steel wire basket  
ABS drawer dividers  
Central lock  
Castors with brakes  
Non slip mat  
Luxurious push handle  
Soft closing drawer

**Specifications:**

Dimensions: 650 x 480 x 1000mm  
Drawer Depth: 100mm  
Top loading: 75kg  
Weight: 35kg

All accessories included

**ORDER CODE:  
TTR03**

[Click For More Information](#)

## ANAESTHESIA TROLLEY

Anesthesia trolley with an extension table, waste bin, stainless steel wire basket and soft closing drawers.

### Features:

- Five drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Storage rack

### Specifications:

Dimensions: 650 x 480 x 1000mm  
 Drawer Depth: 3 x 100mm, 1 x 130mm, 1 x 185mm  
 Weight: 55kg

All accessories included

### ORDER CODE: ANT

[Click For More Information](#)



## DRESSING CHANGE TROLLEY

A popular two drawer ABS dressing trolley loaded with features such as extension table, bins, two buckets & a drawer central lock.

### Features:

- Two drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Two mounted buckets

**Specifications:**

Dimensions: 650 x 480 x 1000mm

Drawer Depth: 100mm

Weight: 35kg

All accessories included

**ORDER CODE:****DCTR**[Click For More Information](#)**INFUSION TROLLEY**

A fully featured infusion trolley including bins, basket & adjustable height IV bag pole.

**Features:**

- One drawer
- Adjustable height rack
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

**Specifications:**

Dimensions: 650 x 480 x 1000mm

Drawer Depth: 100mm

Weight: 22kg

All accessories included

**ORDER CODE:****INFTR**[Click For More Information](#)



## INSTRUMENT TROLLEY NO SHELF

A strong ABS instrument trolley for transport of goods and equipment throughout healthcare facilities.

### Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

### Specifications:

Dimensions: 650 x 480 x 1000mm  
Weight: 22kg

**ORDER CODE:**  
**INST00**

[Click For More Information](#)

## INSTRUMENT TROLLEY ONE SHELF

A strong ABS instrument trolley with one shelf for transport of goods throughout healthcare facilities.

### Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

### Specifications:

Dimensions: 650 x 480 x 1000mm  
Weight: 22kg

**ORDER CODE:**  
**INST01**

[Click For More Information](#)





## DIFFICULT AIRWAY TROLLEY

Six drawer difficult airway trolley with extension table, luxurious push handle and all accessories included.

### Features:

- Six drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Lockable end cabinet

### Specifications:

Dimensions: 650 x 480 x 1070mm  
 Drawer Depth: 3 x 100mm, 2 x 125mm, 1 x 185mm  
 Weight: 60kg

All accessories included

### ORDER CODE: DATR

[Click For More Information](#)

## MEDICAL RECORDS TROLLEY NO DRAWER

A full featured medical records trolley including all accessories. A great addition to any healthcare facility.

### Features:

- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- CT file container

### Specifications:

Dimensions: 650 x 480 x 1000mm



Weight: 35kg

All accessories included

**ORDER CODE:  
MEDRTO**

[Click For More Information](#)



## **MEDICAL RECORDS TROLLEY ONE DRAWER**

A full featured single drawer medical records trolley including all accessories. A great addition to any healthcare facility.

### **Features:**

- One drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- CT file container

### **Specifications:**

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 100mm
- Weight: 35kg

All accessories included

**ORDER CODE:  
MEDRT1**

[Click For More Information](#)

# COMPUTER TROLLEYS



## WORKSTATION TROLLEY WITH PRINTER SHELF

A great addition to any healthcare facility or hospital. Our versatile workstation trolley is packed full of features and comes with all accessories included.

### Features:

- Display arm
- Disinfectant Holder
- Laptop box with lock & key (separate keys for each side)
- Laptop sits inside laptop box underneath
- Connects to external screen mounted on top of trolley (Note: Screen not supplied with trolley)
- Gas Lift Handle for height adjustment
- Slide keyboard holder
- Pull-out mouse tray
- Mouse holder
- Printer tray
- Twin wheel braked castors

### Specifications:

- Length: 54cm
- Width: 70cm
- Height range: 119-142cm
- Work surface height: 81-116cm
- Work surface size: 53x48cm
- Printer Tray Size: 51x68cm
- Castors: 100mm twin wheel braked

### Order Code:

**MWST**

[Click For More Information](#)

# INDEX

## A

Anaesthesia Trolley 6

## D

Difficult Airway Trolley 9

Dressing Change Trolley 6

## E

Emergency Trolley 2

## I

Infusion Trolley 7

Instrument Trolley No Shelf 8

Instrument Trolley One Shelf 8

## M

Medical Records Trolley No Drawer 9

Medical Records Trolley One Drawer 10

Medicine Trolley 3

## T

Treatment Trolley One Drawer 4

Treatment Trolley Three Drawer 5

Treatment Trolley Two Drawer 4

## W

Workstation Trolley With Printer Shelf 11



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)