



Pacific Medical Catalogue

TABLE OF CONTENTS

Medical Carts	2
Medical Carts	2
Computer Trolleys	11

MEDICAL CARTS

MEDICAL CARTS



EMERGENCY TROLLEY

A rugged portable emergency trolley constructed with steel & ABS. Along with a heavy load rating and full set of accessories included it's a best seller.

Note: Product pictured with central & seals lock. It can only come with one or the other.

Features:

- Five drawer
- Extension table
- CPR board
- Defibrillator Shelf
- Stainless steel IV pole
- Waste bin
- 304 stainless steel wire basket
- Oxygen tank holder
- External power holder
- ABS drawer dividers
- Central or seals lock
- Noise resistant castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Smooth welding

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-3 110mm, 4 - 135mm, 5 - 200mm
- Drawer load capacity: 50kg
- Overall load capacity: 200kg
- Weight: 55kg

All accessories included

ORDER CODE:

- EMT - Emergency Trolley With Central Lock**
- EMTSL - Emergency Trolley With Seals Lock**

[Click For More Information](#)

MEDICINE TROLLEY

For the safe transport of medications and supplies our four drawer medicine trolley is the perfect fit. Including all accessories and boasting a number of features this product represents great value.

Features:

- Five drawers
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 1070mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-4 125mm, 5 - 260mm
- Top loading: 75kg
- Weight: 42kg

All accessories included

ORDER CODE:
MEDT

[Click For More Information](#)





TREATMENT TROLLEY ONE DRAWER

This single drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- One drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 110mm
- Top loading: 75kg
- Weight: 35kg

All accessories included

ORDER CODE:

TTR01

[Click For More Information](#)

TREATMENT TROLLEY TWO DRAWER

This two drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- Two drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer



Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR02**

[Click For More Information](#)

**TREATMENT TROLLEY THREE DRAWER**

This three drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- Three drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR03**

[Click For More Information](#)

ANAESTHESIA TROLLEY

Anesthesia trolley with an extension table, waste bin, stainless steel wire basket and soft closing drawers.

Features:

- Five drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Storage rack

Specifications:

Dimensions: 650 x 480 x 1000mm
 Drawer Depth: 3 x 110mm, 1 x 140mm, 1 x 195mm
 Weight: 55kg

All accessories included

ORDER CODE: ANT

[Click For More Information](#)



DRESSING CHANGE TROLLEY

A popular two drawer ABS dressing trolley loaded with features such as extension table, bins, two buckets & a drawer central lock.

Features:

- Two drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Two mounted buckets

Specifications:

Dimensions: 650 x 480 x 1000mm

Drawer Depth: 110mm

Weight: 35kg

All accessories included

ORDER CODE:

DCTR

[Click For More Information](#)

INFUSION TROLLEY

A fully featured infusion trolley including bins, basket & adjustable height IV bag pole.

Features:

- One drawer
- Adjustable height rack
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 1000mm

Drawer Depth: 110mm

Weight: 22kg

All accessories included

ORDER CODE:

INFTR

[Click For More Information](#)





INSTRUMENT TROLLEY NO SHELF

A strong ABS instrument trolley for transport of goods and equipment throughout healthcare facilities.

Features:

Luxurious push handle
Strong polymer shelves
Easy to clean
Castors with brakes

Specifications:

Dimensions: 650 x 480 x 1000mm
Weight: 22kg

ORDER CODE:
INST00

[Click For More Information](#)

INSTRUMENT TROLLEY ONE SHELF

A strong ABS instrument trolley with one shelf for transport of goods throughout healthcare facilities.

Features:

Luxurious push handle
Strong polymer shelves
Easy to clean
Castors with brakes

Specifications:

Dimensions: 650 x 480 x 1000mm
Weight: 22kg

ORDER CODE:
INST01

[Click For More Information](#)





DIFFICULT AIRWAY TROLLEY

Six drawer difficult airway trolley with extension table, luxurious push handle and all accessories included.

Features:

- Six drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Lockable end cabinet

Specifications:

- Dimensions: 650 x 480 x 1070mm
- Drawer Depth: 3 x 110mm, 2 x 135mm, 1 x 195mm
- Weight: 60kg

All accessories included

ORDER CODE: DATR

[Click For More Information](#)

MEDICAL RECORDS TROLLEY NO DRAWER

A full featured medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- CT file container

Specifications:

Dimensions: 650 x 480 x 1000mm



Weight: 35kg

All accessories included

**ORDER CODE:
MEDRTO**

[Click For More Information](#)



MEDICAL RECORDS TROLLEY ONE DRAWER

A full featured single drawer medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- One drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- CT file container

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 110mm
- Weight: 35kg

All accessories included

**ORDER CODE:
MEDRT1**

[Click For More Information](#)

COMPUTER TROLLEYS



WORKSTATION TROLLEY WITH PRINTER SHELF

A great addition to any healthcare facility or hospital. Our versatile workstation trolley is packed full of features and comes with all accessories included.

Features:

- Display arm
- Disinfectant Holder
- Laptop box with lock & key (separate keys for each side)
- Laptop sits inside laptop box underneath
- Connects to external screen mounted on top of trolley (Note: Screen not supplied with trolley)
- Gas Lift Handle for height adjustment
- Slide keyboard holder
- Pull-out mouse tray
- Mouse holder
- Printer tray
- Twin wheel braked castors

Specifications:

- Length: 54cm
- Width: 70cm
- Height range: 119-142cm
- Work surface height: 81-116cm
- Work surface size: 53x48cm
- Printer Tray Size: 51x68cm
- Castors: 100mm twin wheel braked

Order Code:

MWST

[Click For More Information](#)

INDEX

A

Anaesthesia Trolley 6

D

Difficult Airway Trolley 9

Dressing Change Trolley 6

E

Emergency Trolley 2

I

Infusion Trolley 7

Instrument Trolley No Shelf 8

Instrument Trolley One Shelf 8

M

Medical Records Trolley No Drawer 9

Medical Records Trolley One Drawer 10

Medicine Trolley 3

T

Treatment Trolley One Drawer 4

Treatment Trolley Three Drawer 5

Treatment Trolley Two Drawer 4

W

Workstation Trolley With Printer Shelf 11



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au