



Pacific Medical Catalogue

TABLE OF CONTENTS

Medical Carts 2

 Medical Carts 2

 Computer Trolleys 11

MEDICAL CARTS

MEDICAL CARTS



EMERGENCY TROLLEY

A rugged portable emergency trolley constructed with steel & ABS. Along with a heavy load rating and full set of accessories included it's a best seller.

Note: Product pictured with central & seals lock. It can only come with one or the other.

Features:

- Five drawer
- Extension table
- CPR board
- Defibrillator Shelf
- Stainless steel IV pole
- Waste bin
- 304 stainless steel wire basket
- Oxygen tank holder
- External power holder
- ABS drawer dividers
- Central or seals lock
- Noise resistant castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Smooth welding

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-3 110mm, 4 - 135mm, 5 - 200mm
- Drawer load capacity: 30kg
- Overall load capacity: 200kg
- Weight: 55kg

All accessories included

ORDER CODE:

EMT - Emergency Trolley With Central Lock
EMTSL - Emergency Trolley With Seals Lock

[Click For More Information](#)

MEDICINE TROLLEY

For the safe transport of medications and supplies our four drawer medicine trolley is the perfect fit. Including all accessories and boasting a number of features this product represents great value.

Features:

- Five drawers
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 1070mm
- Drawer Dimensions: 550 x 400mm
- Drawer Depths: 1-4 125mm, 5 - 260mm
- Top loading: 75kg
- Weight: 42kg

All accessories included

ORDER CODE: MEDT

[Click For More Information](#)



TREATMENT TROLLEY ONE DRAWER

This single drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- One drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 110mm
- Top loading: 75kg
- Weight: 35kg

All accessories included

ORDER CODE:
TTR01

[Click For More Information](#)

TREATMENT TROLLEY TWO DRAWER

This two drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

- Two drawer
- Bottom shelf
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer



Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR02**

[Click For More Information](#)

**TREATMENT TROLLEY THREE DRAWER**

This three drawer treatment trolley comes with all accessories as standard to ensure you have everything you need out of the box.

Features:

Three drawer
Bottom shelf
Extension table
Waste bin
304 stainless steel wire basket
ABS drawer dividers
Central lock
Castors with brakes
Non slip mat
Luxurious push handle
Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Top loading: 75kg
Weight: 35kg

All accessories included

**ORDER CODE:
TTR03**

[Click For More Information](#)

ANAESTHESIA TROLLEY

Anesthesia trolley with an extension table, waste bin, stainless steel wire basket and soft closing drawers.

Features:

- Five drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Storage rack

Specifications:

Dimensions: 650 x 480 x 1000mm
 Drawer Depth: 3 x 110mm, 1 x 140mm, 1 x 195mm
 Weight: 55kg

All accessories included

ORDER CODE:

ANT

[Click For More Information](#)



DRESSING CHANGE TROLLEY

A popular two drawer ABS dressing trolley loaded with features such as extension table, bins, two buckets & a drawer central lock.

Features:

- Two drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Two mounted buckets

Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Weight: 35kg

All accessories included

ORDER CODE:
DCTR

[Click For More Information](#)

INFUSION TROLLEY

A fully featured infusion trolley including bins, basket & adjustable height IV bag pole.

Features:

One drawer
Adjustable height rack
Extension table
Waste bin
304 stainless steel wire basket
ABS drawer dividers
Castors with brakes
Non slip mat
Luxurious push handle
Soft closing drawer

Specifications:

Dimensions: 650 x 480 x 1000mm
Drawer Depth: 110mm
Weight: 22kg

All accessories included

ORDER CODE:
INFTR

[Click For More Information](#)





INSTRUMENT TROLLEY NO SHELF

A strong ABS instrument trolley for transport of goods and equipment throughout healthcare facilities.

Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

Specifications:

Dimensions: 650 x 480 x 1000mm
Weight: 22kg

ORDER CODE:
INST00

[Click For More Information](#)

INSTRUMENT TROLLEY ONE SHELF

A strong ABS instrument trolley with one shelf for transport of goods throughout healthcare facilities.

Features:

- Luxurious push handle
- Strong polymer shelves
- Easy to clean
- Castors with brakes

Specifications:

Dimensions: 650 x 480 x 1000mm
Weight: 22kg

ORDER CODE:
INST01

[Click For More Information](#)





DIFFICULT AIRWAY TROLLEY

Six drawer difficult airway trolley with extension table, luxurious push handle and all accessories included.

Features:

- Six drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Central lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- Lockable end cabinet

Specifications:

Dimensions: 650 x 480 x 1070mm
 Drawer Depth: 3 x 110mm, 2 x 135mm, 1 x 195mm
 Weight: 60kg

All accessories included

ORDER CODE:
DATR

[Click For More Information](#)

MEDICAL RECORDS TROLLEY NO DRAWER

A full featured medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- CT file container

Specifications:

Dimensions: 650 x 480 x 1000mm



Weight: 35kg

All accessories included

**ORDER CODE:
MEDRT0**

[Click For More Information](#)



MEDICAL RECORDS TROLLEY ONE DRAWER

A full featured single drawer medical records trolley including all accessories. A great addition to any healthcare facility.

Features:

- One drawer
- Extension table
- Waste bin
- 304 stainless steel wire basket
- ABS drawer dividers
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawer
- CT file container

Specifications:

- Dimensions: 650 x 480 x 1000mm
- Drawer Depth: 110mm
- Weight: 35kg

All accessories included

**ORDER CODE:
MEDRT1**

[Click For More Information](#)

COMPUTER TROLLEYS



WORKSTATION TROLLEY WITH PRINTER SHELF

A great addition to any healthcare facility or hospital. Our versatile workstation trolley is packed full of features and comes with all accessories included.

Features:

- Display arm
- Disinfectant Holder
- Laptop box with lock & key (separate keys for each side)
- Laptop sits inside laptop box underneath
- Connects to external screen mounted on top of trolley (Note: Screen not supplied with trolley)
- Gas Lift Handle for height adjustment
- Slide keyboard holder
- Pull-out mouse tray
- Mouse holder
- Printer tray
- Twin wheel braked castors

Specifications:

- Length: 54cm
- Width: 70cm
- Height range: 119-142cm
- Work surface height: 81-116cm
- Work surface size: 53x48cm
- Printer Tray Size: 51x68cm
- Castors: 100mm twin wheel braked

Order Code:

MWST

[Click For More Information](#)

INDEX

A

Anaesthesia Trolley 6

D

Difficult Airway Trolley 9

Dressing Change Trolley 6

E

Emergency Trolley 2

I

Infusion Trolley 7

Instrument Trolley No Shelf 8

Instrument Trolley One Shelf 8

M

Medical Records Trolley No Drawer 9

Medical Records Trolley One Drawer 10

Medicine Trolley 3

T

Treatment Trolley One Drawer 4

Treatment Trolley Three Drawer 5

Treatment Trolley Two Drawer 4

W

Workstation Trolley With Printer Shelf 11



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au