



Pacific Medical Catalogue

TABLE OF CONTENTS

New Products 2

NEW PRODUCTS



SADDLE STOOL WITH BACKREST

High quality saddle stool with backrest. It has adjustable seat and backrest angle along with backrest height and angle. Suitable for many applications.

Features:

- Adjustable seat tilt
- Height and angle adjustable backrest
- Moulded foam seat

Specifications:

- Adjustable height: 47cm - 63cm
- Seat diameter: 42cm
- Backrest size: 36.5 x 26.5cm
- Backrest height range: 23 - 38cm
- Seat angle range: -3 to +6 degrees
- Backrest angle adjustment: 25 degrees
- Foam thickness: 5cm
- Weight Capacity: 160kg

ORDER CODES:

SSNBBR - Navy Blue

SSBLBR - Black

SSGRBR - Grey

[Click For More Information](#)

ROUND STOOL WITH BACKREST

High quality round stool with backrest. It has adjustable backrest angle along with height. Suitable for many applications.

Features:

Adjustable backrest tilt
Height adjustable backrest
Comfortable foam seat

Specifications:

Adjustable height: 47cm - 63cm
Seat diameter: 40cm
Backrest size: 36.5 x 26.5cm
Backrest height range: 26.5 - 39cm
Backrest angle adjustment: 25 degrees
Foam thickness: 5cm
Weight Capacity: 160kg

ORDER CODES:

RSNBBR - Navy Blue

RSBLBR - Black

RSGRBR - Grey

[Click For More Information](#)



ROUND STOOL WITH ARMREST

High quality round stool with armrest. It has a swing around armrest with height adjustment. Suitable for many applications.

Features:

Swing around armrest
Height adjustable armrest
Comfortable foam seat

Specifications:

Adjustable height: 48cm - 64cm
Seat diameter: 40cm
Armrest size: 53 x 8 x 5cm
Armrest height range: 23 - 31cm
Foam thickness: 7cm
Weight Capacity: 160kg

ORDER CODES:

RSNBAR - Navy Blue

RSBLAR - Black



RSGRAR - Grey[Click For More Information](#)**STEALTH FLEXION DISTRACTION CHIROPRACTIC TABLE**

Our new manipulation chiropractic table is packed full with features. It gives you a wide range of treatment options from drops, breakaway, flexion and distraction with functionally adjustable headpieces. All drops are hand activated with tension controls and the removeable flexion control T-Bar will allow you to move freely around the table.

Features:

- Electric Hi Lo operation by foot switch
- Electric flexion function operated by foot switch
- Abdominal breakaway section diversifying treatment options
- Height and width adjustable padded headrest
- Removable flexion control handle
- Extendable padded footrest
- Luxurious and comfortable table
- Includes patient straps
- Sleek dark frame two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 172cm
- Width: 56-80cm
- Headrest height adjustment: 10cm
- Four tension controlled drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical to +20 Degrees
- Adjustable angle of thoracic from -15 to +18 Degrees
- Adjustable angle of pelvic (tilting up and down) -22 to +22 Degrees
- Adjustable angle of pelvic (left to right) 20 Degrees
- Height range from 56cm to 80cm
- 5cm wheels
- Padding thickness 7cm
- Foam Density: 100kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:**CTFDST - Stealth Black**

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.



[Click For More Information](#)

STEALTH MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

STEALTH PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET3BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





PROCEDURE CHAIR

Our procedure chair is a versatile multi modality product utilising fully electric functionality and a strong column lift for the rigours of modern day use.

Features:

- Fully electric operation by hand switch
- Electric seat tilt
- Strong column lift
- Small footprint frame with rounded frame legs for safety
- Gas strut adjustable face crest
- Padded armrests
- Debris tray
- Contoured backrest
- Can achieve Trendelenburg position
- 10cm central lock castors
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 198cm
- Width: 64cm
- Headrest length: 31cm
- Backrest length: 61cm
- Seat section length: 46cm
- Leg section length: 48cm
- Adjustable motorised back section from 0 to +85 Degrees
- Adjustable motorised foot section from 0 to -85 Degrees
- Adjustable motorised seat tilt from -5 to +35 Degrees
- Adjustable head rest through gas operated struts from -15 to +35 Degrees
- Height range from 59cm to 89cm in approx 30 seconds
- Max working heel height 145cm (with footrest extension and extended)
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PRCHNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

STAINLESS STEEL OXYGEN TROLLEY SINGLE DRAWER

A high grade single drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Single drawer
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE:
SST01OXY

[Click For More Information](#)



STAINLESS STEEL OXYGEN TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Two drawers
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

125mm deep drawers
95mm lockable castors
Support for oxygen bottle up to E size
Dimensions: 50 x 75 x 90cm (LxWxH)

**ORDER CODE:
SST02OXY**

[Click For More Information](#)

MOBILE PRIVACY SCREEN

A robust and lightweight antimicrobial mobile privacy screen. Mobile and suitable for simple deployment in a number of medical settings.

Features:

Lightweight aluminium frame
Antimicrobial panels
Corrosion resistant
Dirt resistant
Scratch resistant
Easy to clean & disinfect
7.5cm silent castors
Integrated handle in screen panel
Durable construction
Simple assembly

Specification:

Total Extended Length: 210cm
Height: 185cm
Number of panels: 8
Colour: White

**ORDER CODES:
PSCR**

[Click For More Information](#)



STAINLESS STEEL DEEP WASH TUB



A hardy multi purpose product manufactured from 304 stainless. Suitable for examination, grooming, bathing, dental examinations and much more.

Features:

- Durable 304 stainless steel
- Suitable for grooming and examination
- Removable and sliding top grill
- Removable bottom grill
- Removable legs for easy transport and storage

Specification:

- Length: 167cm
- Width: 58cm
- Height: 90cm
- Deep End Depth: 47.5cm
- Shallow End Depth: 6cm

ORDER CODES:

SSWTD

[Click For More Information](#)

THREE FUNCTION HOSPITAL BED

A strong, reliable and high quality three function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for aged care and hospital facilities.

Three Function Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Removeable head & foot board
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

Three Function Hospital Bed Specification:

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm



Hoist clearance: 16cm
Mattress platform length: 195cm
Mattress platform width: 90cm
Mattress size: 195x90x10cm
Backrest angle: 0-73 degrees
Knee break angle: 0-30 degrees
Leg angle: 0-16 degrees
SWL: 225KG
Warranty: 3 Years
Certification: CE, FDA, ISO13485-2016

ORDER CODES:
HB3

MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.

[Click For More Information](#)

FIVE FUNCTION HOSPITAL BED



A strong, reliable and high quality five function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for higher needs patients and longer term patient care.

Five Function Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

Five Function Hospital Bed Specification:

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm
- Hoist clearance: 16cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Mattress size: 195x90x10cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-16 degrees

Trendelenburg & reverse trendelenburg: 0-13 degrees
SWL: 225KG
Warranty: 3 Years
Certification: CE, FDA, ISO13485-2016

ORDER CODES:
HB5

MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.

[Click For More Information](#)

TILTING ICU HOSPITAL BED

A strong, reliable and high quality five function tilting ICU hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for high needs patients requiring long term treatment with the tilting function allowing treatment of the nervous system with weakness, pain, muscle spasm, and for auxiliary rehabilitation purposes such as hemiplegia, paraplegia and traumatic brain injury and other severe patients to resume standing training



Tilting ICU Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Tilt inclinometer
- Emergency stop button
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Durable patient belts
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

Tilting ICU Hospital Bed Specification:

- Overall length: 215cm
- Overall width: 110cm
- Height range: 52-76cm
- Hoist clearance: 14cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-20 degrees

Trendelenburg & reverse trendelenburg: 0-16 degrees & 0-13 degrees
Tilting angle: 0-85 degrees
SWL: 225KG
Warranty: 3 Years
Certification: CE, FDA, ISO13485-2016

ORDER CODES:
HBT

Please note these are an indent item and are manufactured to order.

MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.

[Click For More Information](#)



FINGERTIP PULSE OXIMETER

A small pulse oximeter intended for measuring the pulse oxygen saturation, pulse rate and perfusion index through a users finger. Suitable for monitoring of adult & paediatric patients in medical clinics, hospitals and at home.

Fingertip Pulse Oximeter Features:

- Special splash proof and drop resistant design
- OLED display of SpO2, PR, PI, pulse bar & waveform
- Spot check and continuous pulse oximeter measurement modes
- Pulse rate analysis for spot check measurement
- Up to 12 groups SpO2 data storage
- Audible & visual alarm with low battery indication
- Automatic power on/off
- Four direction display
- Artefact removal & anti motion
- Accurate & sensitive

Fingertip Pulse Oximeter Specification:

Spo2

Transducer: Dual-wavelength LED sensor
Measuring range: 35-100%
Measuring accuracy: $\leq 2\%$ range from 70-100%

Pulse Rate

Measuring range: 30-240bpm
Measuring accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$ (whichever is the greater)

Perfusion Index

Display range: 0-20%

Default Alarm Limit

SpO2 low limit: 90%

Pulse rate high limit: 120bpm
Pulse rate low limit: 50bpm

Measuring Mode

30 seconds spot check mode
Continuous mode

Dimensions, Weight & Power

Dimensions: 56x34x30mm (L x W x H)
Weight: 52g
Power: 2 x AAA alkaline battery

ORDER CODE:**FPO**[Click For More Information](#)

NITROSPRAY CRYOTHERAPY FLASK

Premier Nitrospray cryotherapy flask constructed of a double-walled vacuum stainless steel container to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage & disengage spray tips and other attachments as required.

Features:

- Double walled vacuum cryotherapy flask
- Easy connect luer lock connector
- Relief valve precisely regulates pressure
- Curved trigger provides unobstructed view of treatment site
- Easy to hold grooved head
- Durable stainless steel construction
- Spray tip protector prevents accidental damage and contact with patient

Specification:

- Holding Time: 6-8 Hours
- Intermittent Use: 3-6 Hours Treatment Time
- Package Inclusions: Flask, Set of five spray tips & Spray tip protector
- Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)
- Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

ORDER CODES:**LNSLITEPLUS - 300ml****LNSPLUS - 475ml**[Click For More Information](#)

CRYOTHERAPY PACKAGE ONE

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle, Extractor & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle, Extractor & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal

Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK1

[Click For More Information](#)



CRYOTHERAPY PACKAGE TWO

A great startup solution to any practice cryotherapy needs. Includes a Nitrospray 300ml, Dewar 20lt, Ladle & Gloves

Includes a Nitrospray 300ml flask, Dewar 20lt, Ladle & Gloves

Flask Specification:

Holding Time: 6-8 Hours

Intermittent Use: 3-6 Hours Treatment Time

Package Inclusions: Flask, Set of five spray tips & Spray tip protector

Capacity: 300 (Nitrospray Lite) & 475ml (Nitrospray Plus)

Construction Materials: 303 Stainless Steel, 304 Stainless Steel, Acetal



Dewar Specification:

Capacity: 20L

Weight (empty): 9.5kg

Caliber (neck opening of dewar): 50mm

Outer diameter of dewar: 394mm

Height of dewar: 672mm

Static evaporation rate: 0.12L/d

Static holding time: 168 days

Ladle outer diameter: 38mm

Ladle height: 120/276mm

No. ladles per dewar: 3~6

Order Code:

LNPACK2

[Click For More Information](#)

THREE SECTION ALL ELECTRIC TREATMENT COUCH

A high quality all electric treatment couch facilitating effective patient treatment all whilst looking after practitioner OH&S.

Features:

Electric Hi Lo operation by foot switch

Electric backrest & footrest for effortless operation

Sleek new design with rounded frame legs for safety

Head section has face hole and plug

Adjustable feet for uneven surfaces

7.5cm retractable castors for safety and mobility

Reliable TiMotion motors



Durable foam padding

Specification:

Total Length: 195cm

Width: 71cm

Head section length: 63cm

Mid section length: 63cm

Leg section length: 63cm

Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees

Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees

Height range from 47cm to 91cm in approx 30 seconds

Padding thickness 5cm

Foam Density: 36kg/m3

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

Warranty: 7 Years Construction, 3 Years Electrics

Certification: CE, FDA, ISO13485-2016

Weight Capacity: 250KG

ORDER CODES:

ET33AENB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

PML1 LED EXAMINATION LIGHT

The PML1 LED examination light is a great all-round diagnostic light. It provides a soft white light for better colour rendition and easier diagnosis and treatment. With a small light head and convenient handle with integrated on/off switch it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations and Minor Surgical Procedures.

Features:

Powerful long-life LED with up to 50,000 hours of operation

Low energy consumption

IR and heat free cold white light

Comfortable guiding handle with integrated On/Off switch

Cooling through a heatsink so no fans or moving parts



Superb colour rendering
Mobile base and wall mount provided with every light

Specifications:

Light Source: 3W LED
Light Intensity: 89,000 Lux @ 30cm / 40,000 Lux @ 50cm
Light Field Diameter: 80-150mm
Colour Temperature: 4,500-5,000K
Lamp Head Diameter: 75mm
Lamp Head Length: 110mm
Base Diameter: 530mm
Base Height: 790mm
Gooseneck Length: 760mm
Total Height: 1,830mm

Order Codes:

PML1WH - WHITE

PML1BL - BLACK

[Click For More Information](#)

**PML2 LED EXAMINATION LIGHT**

The PML2 LED examination light has 6 powerful LED's offering powerful illumination and great value. It provides white light for better colour rendition and easier diagnosis and treatment. With a small light head with integrated handle it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area. A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations & Minor Surgical Procedures.

Features:

Powerful long-life LED with up to 50,000 hours of operation
Low energy consumption
IR and heat free cold white light
Small light head with integrated handle
Cooling through a heatsink so no fans or moving parts
Superb colour rendering
Mobile base and wall mount provided with every light

Specifications:

Light Source: 6W LED
Light Intensity: 78,000 Lux @ 30cm / 40,000 Lux @ 50cm / 20,000 Lux @ 100cm
Light Field Diameter: 70-200mm
Colour Temperature: 4,500-5,500K

Lamp Head Diameter: 80mm
Base Diameter: 530mm
Base Height: 790mm
Gooseneck: Length 760mm
Total Height: 1,815mm

Order Codes:
PML2WH - WHITE
PML2BL - BLACK

[Click For More Information](#)

TWO SECTION CARDIOLOGY COUCH

A strong two section cardiology couch with comfortable sewn edge upholstery, central lock castors, electric backrest & two special cardiac cut outs allowing increased patient access. With a sleek design, great features and foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Dual cut outs for easy patient access
- Comfortable sewn edge upholstery
- Sleek design with rounded frame legs for safety
- 10cm central lock castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 67cm
- Leg section length: 123cm
- Cutout Size: 31cm x 25cm
- Cutout Mechanism: Manual opening and held in place mechanically, very strong.
- Adjustable head section through electric motor & foot switch from flat to +75 Degrees
- Height range from 58cm to 102cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:
CARNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly



to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



SINGLE SECTION TREATMENT / CHANGE COUCH WITH SIDE RAILS

This strong single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Padded side rails
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 180cm
- Width: 66cm
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CHNBSB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

VETERINARY LIFT TABLE

Features:

Electric Hi Lo operation by hand switch & footbar
Durable 304 stainless steel patient surface
Sleek new design with rounded frame legs for safety
10cm lockable castors for safety and mobility
Reliable TiMotion motors
Rechargeable battery

Specification:

Total Length: 110cm
Width: 56cm
Height range from 41cm to 96cm in approx 30 seconds
Warranty: 7 Years Construction, 3 Years Electrics
Certification: CE, FDA, ISO13485-2016
Weight Capacity: 250KG

ORDER CODES:

ETVET

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



INDEX

C

Cryotherapy Package One 15

Cryotherapy Package Two 16

F

Fingertip Pulse Oximeter 13

Five Function Hospital Bed 11

M

Mobile Privacy Screen 9

N

Nitrospray Cryotherapy Flask 14

P

PML1 LED Examination Light 17

PML2 LED Examination Light 18

Procedure Chair 7

R

Round Stool With Armrest 3

Round Stool With Backrest 3

S

Saddle Stool With Backrest 2

Single Section Treatment / Change Couch With Side Rails 20

Stainless Steel Deep Wash Tub 10

Stainless Steel Oxygen Trolley Single Drawer 8

Stainless Steel Oxygen Trolley Two Drawer 8

Stealth Flexion Distraction Chiropractic Table 4

Stealth Medical Treatment Couch 5

Stealth Physio Treatment Couch 6

T

Three Function Hospital Bed 10

Three Section All Electric Treatment Couch 16

Tilting ICU Hospital Bed 12

Two Section Cardiology Couch 19

V

Veterinary Lift Table 21



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au