



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**New Products** ..... 2

# NEW PRODUCTS



## SADDLE STOOL WITH BACKREST

High quality saddle stool with backrest. It has adjustable seat and backrest angle along with backrest height and angle. Suitable for many applications.

### Features:

- Adjustable seat tilt
- Height and angle adjustable backrest
- Moulded foam seat

### Specifications:

- Adjustable height: 47cm - 63cm
- Seat diameter: 42cm
- Backrest size: 36.5 x 26.5cm
- Backrest height range: 23 - 38cm
- Seat angle range: -3 to +6 degrees
- Backrest angle adjustment: 25 degrees
- Foam thickness: 5cm
- Weight Capacity: 160kg

### ORDER CODES:

**SSNBBR - Navy Blue**

**SSBLBR - Black**

**SSGRBR - Grey**

[Click For More Information](#)

## ROUND STOOL WITH BACKREST

High quality round stool with backrest. It has adjustable backrest angle along with height. Suitable for many applications.

### Features:

Adjustable backrest tilt  
Height adjustable backrest  
Comfortable foam seat

### Specifications:

Adjustable height: 47cm - 63cm  
Seat diameter: 40cm  
Backrest size: 36.5 x 26.5cm  
Backrest height range: 26.5 - 39cm  
Backrest angle adjustment: 25 degrees  
Foam thickness: 5cm  
Weight Capacity: 160kg

### ORDER CODES:

**RSNBBR - Navy Blue**

**RSBLBR - Black**

**RSGRBR - Grey**

[Click For More Information](#)



## ROUND STOOL WITH ARMREST

High quality round stool with armrest. It has a swing around armrest with height adjustment. Suitable for many applications.

### Features:

Swing around armrest  
Height adjustable armrest  
Comfortable foam seat

### Specifications:

Adjustable height: 48cm - 64cm  
Seat diameter: 40cm  
Armrest size: 53 x 8 x 5cm  
Armrest height range: 23 - 31cm  
Foam thickness: 7cm  
Weight Capacity: 160kg

### ORDER CODES:

**RSNBAR - Navy Blue**

**RSBLAR - Black**



**RSGRAR - Grey**[Click For More Information](#)**STEALTH FLEXION DISTRACTION CHIROPRACTIC TABLE**

Our new manipulation chiropractic table is packed full with features. It gives you a wide range of treatment options from drops, breakaway, flexion and distraction with functionally adjustable headpieces. All drops are hand activated with tension controls and the removeable flexion control T-Bar will allow you to move freely around the table.

**Features:**

- Electric Hi Lo operation by foot switch
- Electric flexion function operated by foot switch
- Abdominal breakaway section diversifying treatment options
- Height and width adjustable padded headrest
- Removable flexion control handle
- Extendable padded footrest
- Luxurious and comfortable table
- Includes patient straps
- Sleek dark frame two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

**Specification:**

- Total Length: 172cm
- Width: 56-80cm
- Headrest height adjustment: 10cm
- Four tension controlled drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical to +20 Degrees
- Adjustable angle of thoracic from -15 to +18 Degrees
- Adjustable angle of pelvic (tilting up and down) -22 to +22 Degrees
- Adjustable angle of pelvic (left to right) 20 Degrees
- Height range from 56cm to 80cm
- 5cm wheels
- Padding thickness 7cm
- Foam Density: 100kg/m<sup>3</sup>
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

**ORDER CODES:****CTFDST - Stealth Black**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**



[Click For More Information](#)

## STEALTH MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

### Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

### Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m<sup>3</sup>
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

**ET33BLST - Black**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)

## STEALTH PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

### Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

### Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m<sup>3</sup>
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

**ET3BLST - Black**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)





## PROCEDURE CHAIR

Our procedure chair is a versatile multi modality product utilising fully electric functionality and a strong column lift for the rigours of modern day use.

### Features:

- Fully electric operation by hand switch
- Electric seat tilt
- Strong column lift
- Small footprint frame with rounded frame legs for safety
- Gas strut adjustable face crest
- Padded armrests
- Debris tray
- Contoured backrest
- Can achieve Trendelenburg position
- 10cm central lock castors
- Reliable TiMotion motor
- Comfortable foam padding

### Specification:

- Total Length: 198cm
- Width: 64cm
- Headrest length: 31cm
- Backrest length: 61cm
- Seat section length: 46cm
- Leg section length: 48cm
- Adjustable motorised back section from 0 to +85 Degrees
- Adjustable motorised foot section from 0 to -85 Degrees
- Adjustable motorised seat tilt from -5 to +35 Degrees
- Adjustable head rest through gas operated struts from -15 to +35 Degrees
- Height range from 59cm to 89cm in approx 30 seconds
- Max working heel height 145cm (with footrest extension and extended)
- Padding thickness 7.5cm
- Foam Density: 36kg/m<sup>3</sup>
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

**PRCHNB - Navy Blue**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)



## STAINLESS STEEL OXYGEN TROLLEY SINGLE DRAWER

A high grade single drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

### Features:

- High grade 304 stainless steel throughout
- Single drawer
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

### Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

**ORDER CODE:**  
**SST01OXY**

[Click For More Information](#)



## STAINLESS STEEL OXYGEN TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

### Features:

- High grade 304 stainless steel throughout
- Two drawers
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

**Specifications:**

125mm deep drawers  
95mm lockable castors  
Support for oxygen bottle up to E size  
Dimensions: 50 x 75 x 90cm (LxWxH)

**ORDER CODE:  
SST02OXY**

[Click For More Information](#)

**MOBILE PRIVACY SCREEN**

A robust and lightweight antimicrobial mobile privacy screen. Mobile and suitable for simple deployment in a number of medical settings.

**Features:**

Lightweight aluminium frame  
Antimicrobial panels  
Corrosion resistant  
Dirt resistant  
Scratch resistant  
Easy to clean & disinfect  
7.5cm silent castors  
Integrated handle in screen panel  
Durable construction  
Simple assembly

**Specification:**

Total Extended Length: 210cm  
Height: 185cm  
Number of panels: 8  
Colour: White

**ORDER CODES:  
PSCR**

[Click For More Information](#)



## STAINLESS STEEL DEEP WASH TUB



A hardy multi purpose product manufactured from 304 stainless. Suitable for examination, grooming, bathing, dental examinations and much more.

### Features:

- Durable 304 stainless steel
- Suitable for grooming and examination
- Removable and sliding top grill
- Removable bottom grill
- Removable legs for easy transport and storage

### Specification:

- Length: 167cm
- Width: 58cm
- Height: 90cm
- Deep End Depth: 47.5cm
- Shallow End Depth: 6cm

### ORDER CODES:

**SSWTD**

[Click For More Information](#)

## THREE FUNCTION HOSPITAL BED

A strong, reliable and high quality three function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for aged care and hospital facilities.

### Three Function Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Removeable head & foot board
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### Three Function Hospital Bed Specification:

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm



Hoist clearance: 16cm  
Mattress platform length: 195cm  
Mattress platform width: 90cm  
Mattress size: 195x90x10cm  
Backrest angle: 0-73 degrees  
Knee break angle: 0-30 degrees  
Leg angle: 0-16 degrees  
SWL: 225KG  
Warranty: 3 Years  
Certification: CE, FDA, ISO13485-2016

**ORDER CODES:**  
**HB3**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)

## **FIVE FUNCTION HOSPITAL BED**



A strong, reliable and high quality five function hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for higher needs patients and longer term patient care.

### **Five Function Hospital Bed Features:**

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### **Five Function Hospital Bed Specification:**

- Overall length: 220cm
- Overall width: 102cm
- Height range: 44-74cm
- Hoist clearance: 16cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Mattress size: 195x90x10cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-16 degrees

Trendelenburg & reverse trendelenburg: 0-13 degrees  
 SWL: 225KG  
 Warranty: 3 Years  
 Certification: CE, FDA, ISO13485-2016

**ORDER CODES:**  
**HB5**

**MAINTENANCE:** Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.

[Click For More Information](#)

## TILTING ICU HOSPITAL BED

A strong, reliable and high quality five function tilting ICU hospital bed. Well featured including removeable boards, removeable ABS side rails, drainage bag holders, iv pole & 125mm central lock castors. An ideal product for high needs patients requiring long term treatment with the tilting function allowing treatment of the nervous system with weakness, pain, muscle spasm, and for auxiliary rehabilitation purposes such as hemiplegia, paraplegia and traumatic brain injury and other severe patients to resume standing training



### Tilting ICU Hospital Bed Features:

- Electric operation by handset with integral function lockout
- Electrically operated backrest, height & knee break angle
- Backrest manual CPR release on both sides
- Tilt inclinometer
- Emergency stop button
- Removeable head & foot board
- Patient handset & nurse control panel
- Battery backup
- Four ABS safety side rails
- Four drainage bag holders on both sides
- Durable patient belts
- Heavy duty IV pole
- Auto contour simultaneous adjustment of backrest & knee break
- 125mm heavy duty central locking castors
- Foam mattress

### Tilting ICU Hospital Bed Specification:

- Overall length: 215cm
- Overall width: 110cm
- Height range: 52-76cm
- Hoist clearance: 14cm
- Mattress platform length: 195cm
- Mattress platform width: 90cm
- Backrest angle: 0-75 degrees
- Knee break angle: 0-35 degrees
- Leg angle: 0-20 degrees

Trendelenburg & reverse trendelenburg: 0-16 degrees & 0-13 degrees  
Tilting angle: 0-85 degrees  
SWL: 225KG  
Warranty: 3 Years  
Certification: CE, FDA, ISO13485-2016

**ORDER CODES:**  
**HBT**

**Please note these are an indent item and are manufactured to order.**

**MAINTENANCE: Please make sure your bed is serviced yearly to ensure your bed functions to its optimum. Lubrication of pivot points and moving points on your bed must be done on a regular basis.**

[Click For More Information](#)



## FINGERTIP PULSE OXIMETER

A small pulse oximeter intended for measuring the pulse oxygen saturation, pulse rate and perfusion index through a users finger. Suitable for monitoring of adult & paediatric patients in medical clinics, hospitals and at home.

### Fingertip Pulse Oximeter Features:

- Special splash proof and drop resistant design
- OLED display of SpO2, PR, PI, pulse bar & waveform
- Spot check and continuous pulse oximeter measurement modes
- Pulse rate analysis for spot check measurement
- Up to 12 groups SpO2 data storage
- Audible & visual alarm with low battery indication
- Automatic power on/off
- Four direction display
- Artefact removal & anti motion
- Accurate & sensitive

### Fingertip Pulse Oximeter Specification:

#### Spo2

Transducer: Dual-wavelength LED sensor  
Measuring range: 35-100%  
Measuring accuracy:  $\leq 2\%$  range from 70-100%

#### Pulse Rate

Measuring range: 30-240bpm  
Measuring accuracy:  $\pm 2\text{bpm}$  or  $\pm 2\%$  (whichever is the greater)

#### Perfusion Index

Display range: 0-20%

#### Default Alarm Limit

SpO2 low limit: 90%

Pulse rate high limit: 120bpm  
Pulse rate low limit: 50bpm

**Measuring Mode**

30 seconds spot check mode  
Continuous mode

**Dimensions, Weight & Power**

Dimensions: 56x34x30mm (L x W x H)  
Weight: 52g  
Power: 2 x AAA alkaline battery

**ORDER CODE:**

**FPO**

[Click For More Information](#)

## THREE SECTION ALL ELECTRIC TREATMENT COUCH

A high quality all electric treatment couch facilitating effective patient treatment all whilst looking after practitioner OH&S.

**Features:**

Electric Hi Lo operation by foot switch  
Electric backrest & footrest for effortless operation  
Sleek new design with rounded frame legs for safety  
Head section has face hole and plug  
Adjustable feet for uneven surfaces  
7.5cm retractable castors for safety and mobility  
Reliable TiMotion motors  
Durable foam padding

**Specification:**

Total Length: 195cm  
Width: 71cm  
Head section length: 63cm  
Mid section length: 63cm  
Leg section length: 63cm  
Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees  
Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees  
Height range from 47cm to 91cm in approx 30 seconds  
Padding thickness 5cm  
Foam Density: 36kg/m<sup>3</sup>  
Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant  
Warranty: 7 Years Construction, 3 Years Electrics  
Certification: CE, FDA, ISO13485-2016  
Weight Capacity: 250KG

**ORDER CODES:**

**ET33AENB - Navy Blue**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)

**TWO SECTION CARDIOLOGY COUCH**

A strong two section cardiology couch with comfortable sewn edge upholstery, central lock castors, electric backrest & two special cardiac cut outs allowing increased patient access. With a sleek design, great features and foot switch as standard it is simple to operate straight from the box without any assembly required.

**Features:**

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Dual cut outs for easy patient access
- Comfortable sewn edge upholstery
- Sleek design with rounded frame legs for safety
- 10cm central lock castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

**Specification:**

- Total Length: 195cm
- Width: 66cm
- Head section length: 67cm
- Leg section length: 123cm
- Cutout Size: 31cm x 25cm
- Cutout Mechanism: Manual opening and held in place mechanically, very strong.
- Adjustable head section through electric motor & foot switch from flat to +75 Degrees
- Height range from 58cm to 102cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m<sup>3</sup>
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

**ORDER CODES:**

**CARNB - Navy Blue**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**



[Click For More Information](#)

## **SINGLE SECTION TREATMENT / CHANGE COUCH WITH SIDE RAILS**

This strong single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.



### **Features:**

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Padded side rails
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

### **Specification:**

- Total Length: 180cm
- Width: 66cm
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m<sup>3</sup>
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### **ORDER CODES:**

**CHNBSB - Navy Blue**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)

## VETERINARY LIFT TABLE



### Features:

- Electric Hi Lo operation by hand switch & footbar
- Durable 304 stainless steel patient surface
- Sleek new design with rounded frame legs for safety
- 10cm lockable castors for safety and mobility
- Reliable TiMotion motors
- Rechargeable battery

### Specification:

- Total Length: 110cm
- Width: 56cm
- Height range from 41cm to 96cm in approx 30 seconds
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

**ETVET**

**MAINTENANCE:** Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

# INDEX

## F

Fingertip Pulse Oximeter 13  
Five Function Hospital Bed 11

## M

Mobile Privacy Screen 9

## P

Procedure Chair 7

## R

Round Stool With Armrest 3  
Round Stool With Backrest 3

## S

Saddle Stool With Backrest 2  
Single Section Treatment / Change Couch With Side Rails 16  
Stainless Steel Deep Wash Tub 10  
Stainless Steel Oxygen Trolley Single Drawer 8  
Stainless Steel Oxygen Trolley Two Drawer 8  
Stealth Flexion Distraction Chiropractic Table 4  
Stealth Medical Treatment Couch 5  
Stealth Physio Treatment Couch 6

## T

Three Function Hospital Bed 10  
Three Section All Electric Treatment Couch 14  
Tilting ICU Hospital Bed 12  
Two Section Cardiology Couch 15

## V

Veterinary Lift Table 17



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)