



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**New Products** ..... 2

# NEW PRODUCTS

## PILLOW CASE WITH FLAP BLUE BX200



Easy to use disposable pillow case with flap that provide reliable protection and help maintain a clean, sanitary patient environment.

### Specifications:

75x50cm in size  
SMS 40GSM  
With flap  
Blue in colour  
Box of 200

### ORDER CODES: PCFLAP

[Click For More Information](#)

## PVC PILLOW CASE WITH ZIP LIGHT BLUE BX100

Water-resistant, easy-to-use disposable PVC pillow cases with a zip that provide reliable protection and help maintain a clean, sanitary patient environment.

### Specifications:

75x50cm in size  
Lightweight PVC  
Zip for easy closure  
Light blue in colour  
Box of 100

### ORDER CODES: PCPVCZIP

[Click For More Information](#)





## **PATIENT GOWN WITH SHORT SLEEVES BX100**

High-quality disposable patient gowns offering lightweight comfort and dependable coverage. Perfect for maintaining cleanliness and convenient patient care in healthcare environments.

### **Specifications:**

M (110x125cm) & XL (120x140cm) sizes available  
Non-woven SMS  
Neck & back closed with ties  
Dark blue in colour  
45 GSM  
Box of 100

### **ORDER CODES:**

**PGSSLM - Medium**

**PGSSLXL - Extra Large**

[Click For More Information](#)

## **MODESTY SHEET 100X70CM BX200**

A dark blue, durable & lightweight modesty sheet designed to provide comfort and privacy during treatments and procedures.

### **Specifications:**

100x70cm in size  
SMS 40GSM  
Dark blue in colour  
Box of 200

### **ORDER CODES:**

**MSH100**

[Click For More Information](#)





## MODESTY SHEET 200X120CM BX100

A dark blue, durable & lightweight modesty sheet designed to provide comfort and privacy during treatments and procedures.

### Specifications:

200x120cm in size  
SMS 25GSM  
Dark blue in colour  
Box of 100

### ORDER CODES: MSH200

[Click For More Information](#)

## PACIFIC TOWEL ROLL CTN8

Pacific towel roll an easy and economical absorbent perforated paper roll. Many applications as a hand towel, liner, underlay, surface & device protection.

### Features:

Easy & economical roll format  
Perforated  
Absorbent  
Multi Use

### Specification:

2 Ply  
100 sheets per roll  
Individually shrink wrapped rolls  
Carton of 8 rolls  
Length: 41.5m  
Width: 49cm

### ORDER CODES: PMPTL

[Click For More Information](#)





## PACIFIC BED ROLL CTN6

Pacific towel roll an easy and economical absorbent perforated paper roll. Many applications as a hand towel, liner, underlay, surface & device protection.

### Features:

- Soft yet thick & durable
- Continuous roll
- Strong & hygenic cover

### Specification:

- Individually shrink wrapped rolls
- Carton of 6 rolls
- Length: 80m
- Width: 54cm

### ORDER CODES: PMPBRL

[Click For More Information](#)

## UNDERPAD BX300

A high quality waterproof underpad suitable for use in the most demanding healthcare setting. 5 Ply for greater absorbency and waterproof backing allow better containment of fluids.

### Features:

- Waterproof
- Multi layer construction
- Light blue in colour
- Super absorbent
- Comfortable
- Environmentally friendly

### Specification:

- Fluff pulp mix
- Absorbent polymer
- Waterproof PE Film
- Non skid backing



Box of 300 pieces  
Packaged in 6 bags of 50 pieces  
Sealed on all four sides  
Length: 60cm  
Width: 40cm

**ORDER CODES:  
PMUPL**

[Click For More Information](#)



## **HALF SIZE UNDERPAD BX600**

A high quality waterproof half size underpad suitable for use in the most demanding healthcare setting. 5 Ply for greater absorbency and waterproof backing allow better containment of fluids.

### **Features:**

Waterproof  
Multi layer construction  
Light blue in colour  
Super absorbent  
Comfortable  
Environmentally friendly

### **Specification:**

Fluff pulp mix  
Absorbent polymer  
Waterproof PE Film  
Non skid backing  
Box of 600 pieces  
Packaged in 12 bags of 50 pieces  
Sealed on all four sides  
Length: 40cm  
Width: 30cm

**ORDER CODES:  
PMUPHF**

[Click For More Information](#)

## **AURA TOUCHSCREEN AUTOCLAVE 17 / 22L - CLASS B & S**

The new gold standard in benchtop sterilisation. The AURA Autoclave is loaded with features and European componentry, it will reliably serve any clinic time and time again.

### **Features:**

- Reliable European Componentry
- Sleek Finish
- Easy To Operate User Interface
- Touch Screen & Buttons Allowing Fast Operation
- Inbuilt Printer & USB Functions
- Microprocessor Control System For Precise Sensor Management
- Automated Cleaning Program & Maintenance Reminders
- Automatic RO & Manual Fill Options
- S & B Class Cycles
- Delayed Start Function
- Automatic Component Diagnosis Tool
- Water Quality Sensor
- Five Selectable Languages

### **Specification:**

- External Dimensions: 625 x 440 x 475mm
- Chamber Dimensions: 249 x 355mm (17L) / 249 x 450mm (22L)
- Weight: 52KG (17L) / 55KG (22L)
- Selectable 134°C Cycles: Universal Class B, Fast Class S, Prion, Custom
- Selectable 121°C Cycles: Universal Class B, Fast Class S, Custom
- Test Cycles: B&D Helix, Leak & Vacuum
- Jacket Heater With Six Different Settings
- Drying Settings Available: Normal, Extra Dry 5 min, Extra Dry 10 min
- Included Accessories: Tray Rack, 3 X Trays, 2 X Drain Hoses & Connectors, Tray Handle, Door Seal, Chamber & Water Tank Filter Mesh, Paper Roll, USB Drive, AC Power Cord
- Required Voltage: 120 - 230V AC
- Rated Frequency: 50 / 60 Hz
- Rated Power: 1800 VA
- Sterilising Temperature / Pressure: 121°C @ 100-130 Kpa, 134°C @ 200-240 Kpa
- Water Tank Capacity: 3.5L
- Operating Temperature: +5°C - +40°C
- Relative Humidity Range: 30% - 75%
- Warranty: 2 Years
- Certifications: CE, ISO13485-2016, ASME / Worksafe, EN13060

### **ORDER CODES:**

**AURA17 - 17L**

**AURA22 - 22L**

**NOTE: Professional installation & validation required. Please make sure your autoclave is serviced in accordance with the users manual and machine prompts to ensure it functions to its optimum. All servicing must be completed by our network of authorised agents to maintain machine warranty.**

[Click For More Information](#)



## SELF SEALING AUTOCLAVE POUCH 57X135MM BX200



A simple high quality self seal autoclave pouch 57x135mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

### ORDER CODE: SSP57x135

[Click For More Information](#)

## SELF SEALING AUTOCLAVE POUCH 70X260MM BX200

A simple high quality self seal autoclave pouch 70x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.



### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

### ORDER CODE: SSP70x260

[Click For More Information](#)

## SELF SEALING AUTOCLAVE POUCH 90X160MM BX200



A simple high quality self seal autoclave pouch 90x160mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

### ORDER CODE: SSP90x160

[Click For More Information](#)

## SELF SEALING AUTOCLAVE POUCH 90X260MM BX200

A simple high quality self seal autoclave pouch 90x260mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

### ORDER CODE: SSP90x260

[Click For More Information](#)



## SELF SEALING AUTOCLAVE POUCH 135X285MM BX200



A simple high quality self seal autoclave pouch 135x285mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:  
SSP135x285**

[Click For More Information](#)

## SELF SEALING AUTOCLAVE POUCH 190X360MM BX200

A simple high quality self seal autoclave pouch 190x360mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:  
SSP190x360**

[Click For More Information](#)





## SELF SEALING AUTOCLAVE POUCH 300X380MM BX200

A simple high quality self seal autoclave pouch 300x380mm in size. It features internal and external sterilisation indicators and a strong triple seal.

### Features:

- Internal & external sterilisation indicators
- Indicators for Eto & Steam
- Strong triple seal around pouch
- Strong & simple adhesive end
- Easy closure by hand without the need for a heat sealer
- Blue in colour film
- High quality materials

**ORDER CODE:**  
**SSP300x380**

[Click For More Information](#)

## PACIFIC LIFT CHAIR / RECLINER

The Pacific Lift Chair / Recliner is a well featured dual motor chair with a generous 170KG SWL, removable armrests and 180 degree recline to accommodate most requirements.

### Features:

- Four Castors With Two Lockable
- Two motor with lift function
- Reclining functions up to 180 degrees
- Light up remote control with USB charging port
- Quiet operation
- Removable armrests
- Removable wing headrests
- High quality charcoal grey vinyl
- Reliable TiMotion motors

### Specification:

- Seat Width: 56cm
- Seat Depth: 54cm
- Seat Height: 53cm



Overall Width (At Armrest): 74cm  
 Overall Height: 109cm  
 Total Length After Reclining: 170cm  
 SWL: 170kg

**ORDER CODE:**  
**REC - Charcoal Grey**

**MAINTENANCE: Please make sure your chair is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)



## **BARIATRIC THREE SECTION ALL ELECTRIC TREATMENT COUCH**

### **Features:**

Electric Hi Lo operation by foot switch  
 Electric backrest & footrest for effortless operation  
 Sleek new frame design  
 Head section has face hole and plug  
 12.5cm central lock castors for safety and mobility  
 Reliable TiMotion motors  
 Durable thick foam padding

### **Specification:**

Total Length: 195cm  
 Width: 80cm  
 Head section length: 63cm  
 Mid section length: 63cm  
 Leg section length: 63cm  
 Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees  
 Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees  
 Height range from 52cm to 93cm in approx 30 seconds  
 Hoist Clearance: 15cm  
 Foam Density: 36kg/m<sup>3</sup>  
 Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant  
 Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch  
 Certification: CE, FDA, ISO13485-2016  
 Weight Capacity: 350KG

**ORDER CODES:**  
**ET33BARIANB - Navy Blue**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of**

**pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)

## CHIROPRACTIC TABLE PRO SERIES

PRO series chiropractic drop table is a versatile cost effective product providing treatment flexibility for your patients. With a foot activation bar for hi-lo and improved Cervical, Thoracic, Lumbar & Pelvic drops with hand activation & tensioning.

### Features:

- Electric Hi Lo operation by foot bar
- Height and angle adjustable padded headrest
- Height adjustable & extendable padded footrest
- Sleek design with rounded frame legs and base frame cover plate
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Thick foam padding

### Specification:

- Total Length: 185-200cm
- Width: 55-65cm
- Headrest length: 35cm
- Tapered section length: 26cm
- Mid section length: 23cm
- Long section length: 70cm
- Footrest Length: 18cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable long board to +20 Degrees
- Height adjustable headrest from +20 to -20 degrees
- Height range from 46cm to 80cm in approx 30 seconds
- 7.5cm wheels
- Padding thickness 7.5cm
- Foam Density: 100kg/m<sup>3</sup>
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

### ORDER CODES:

**CTNBPRO - Navy Blue**

**CTBLPRO - Black**

**MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.**

[Click For More Information](#)





## PODIATRY CABINET SINGLE

High quality podiatry cabinet single size featuring a strong metal construction, high quality powder coated finish, glove holder, soft closing drawers & sliding drill shelf.

### Features:

- Three drawers
- Sliding drill shelf
- Soft closing magnetic drawer & drill shelf
- Smooth lockable castors
- Glove holder
- White powder coated finish

### Specifications:

- Drawer depths: 125mm
- Drill cavity: 47 x 45 x 28cm (LxWxH)
- Dimensions: 49 x 51 x 84cm (LxWxH)

### ORDER CODE:

**PCSING**

[Click For More Information](#)

## PODIATRY CABINET DOUBLE

High quality podiatry cabinet double size featuring a strong metal construction, high quality powder coated finish, glove holder, soft closing drawers & sliding drill shelf.

### Features:

- Two full width drawers
- Two half width drawers
- Sliding drill shelf
- Soft closing magnetic drawer & drill shelf
- Smooth lockable castors
- Glove holder
- White powder coated finish

### Specifications:

- Drawer depth (Half Width): 135mm
- Drawer depth (Full Width): 210mm



Drill cavity: 47 x 34 x 24xm (LxWxH)  
Overall dimensions: 49 x 71 x 84cm (LxWxH)

**ORDER CODE:**  
**PCDUBL**

[Click For More Information](#)

# INDEX

## A

AURA Touchscreen Autoclave 17 / 22L - Class B & S 7

## B

Bariatric Three Section All Electric Treatment Couch 12

## C

Chiropractic Table PRO Series 13

## H

Half Size Underpad BX600 6

## M

Modesty Sheet 100x70cm BX200 3

Modesty Sheet 200x120cm BX100 4

## P

Pacific Bed Roll CTN6 5

Pacific Lift Chair / Recliner 11

Pacific Towel Roll CTN8 4

Patient Gown With Short Sleeves BX100 3

Pillow Case With Flap Blue BX200 2

Podiatry Cabinet Double 14

Podiatry Cabinet Single 14

PVC Pillow Case With Zip Light Blue BX100 2

## S

Self Sealing Autoclave Pouch 135x285mm BX200 10

Self Sealing Autoclave Pouch 190x360mm BX200 10

Self Sealing Autoclave Pouch 300x380mm BX200 11

Self Sealing Autoclave Pouch 57x135mm BX200 8

Self Sealing Autoclave Pouch 70x260mm BX200 8

Self Sealing Autoclave Pouch 90x160mm BX200 9

Self Sealing Autoclave Pouch 90x260mm BX200 9

## U

Underpad BX300 5



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)