



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Pillow Cases** ..... 2

# PILLOW CASES

## PILLOW CASE WITH FLAP BLUE BX200



Easy to use disposable pillow case with flap that provide reliable protection and help maintain a clean, sanitary patient environment.

### Specifications:

75x50cm in size  
SMS 40GSM  
With flap  
Blue in colour  
Box of 200

### ORDER CODES: PCFLAP

[Click For More Information](#)

## PVC PILLOW CASE WITH ZIP LIGHT BLUE BX100

Water-resistant, easy-to-use disposable PVC pillow cases with a zip that provide reliable protection and help maintain a clean, sanitary patient environment.

### Specifications:

75x50cm in size  
Lightweight PVC  
Zip for easy closure  
Light blue in colour  
Box of 100

### ORDER CODES: PCPVCZIP

[Click For More Information](#)





# INDEX

## **P**

Pillow Case With Flap Blue BX200 2

PVC Pillow Case With Zip Light Blue BX100 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)