



Pacific Medical Catalogue

TABLE OF CONTENTS

Stainless Steel Products 2

 Specialty Stainless Steel Trolleys 2

 No Rail Trolleys 2

 Single Trolleys 4

 Double Trolleys 10

 Resus & Anaesthesia Trolleys 14

 IV Stands 17

 Mayo Instrument Trolleys 19

 Oxygen Trolleys 20

 Linen Skip Trolleys & Bags 23

STAINLESS STEEL PRODUCTS

SPECIALTY STAINLESS STEEL TROLLEYS

NO RAIL TROLLEYS



STAINLESS STEEL TROLLEY NO RAILS 50X50X90

A high grade stainless steel trolley 50cm wide with no draws or rails. Has lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:

SSTNR50

[Click For More Information](#)

STAINLESS STEEL TROLLEY NO RAILS 60X50X90

A high grade stainless steel trolley 60cm wide with no draws or rails.
Has lockable castors for safety.

Features:

High grade 304 stainless steel throughout
Bottom shelf
Double safety folds
Smooth welding
Smooth running lockable castors
Outstanding finish

Specifications:

95mm lockable castors
Dimensions: 60 x 50 x 90cm (LxWxH)

ORDER CODE:

SSTNR60

[Click For More Information](#)



STAINLESS STEEL TROLLEY NO RAILS 80X50X90

A high grade stainless steel trolley 80cm wide with no draws or rails.
Has lockable castors for safety.

Features:

High grade 304 stainless steel throughout
Bottom shelf
Double safety folds
Smooth welding
Smooth running lockable castors
Outstanding finish

Specifications:

95mm lockable castors
Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:

SSTNR80

[Click For More Information](#)

SINGLE TROLLEYS



STAINLESS STEEL TROLLEY NO DRAWER

A high grade stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE: SST00

[Click For More Information](#)

STAINLESS STEEL TROLLEY ONE DRAWER

A high grade single drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- One drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish



Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:**SST01**[Click For More Information](#)**STAINLESS STEEL TROLLEY TWO DRAWER**

A high grade two drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

High grade 304 stainless steel throughout
Two drawer
Bottom shelf
Double safety folds
Smooth welding
Soft closing magnetic drawer
German ball bearing drawer runners
Smooth running lockable castors
Top and bottom rails on three sides
Outstanding finish

Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:**SST02**[Click For More Information](#)

STAINLESS STEEL TROLLEY THREE DRAWER

A high grade three drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Three drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SST03

[Click For More Information](#)



STAINLESS STEEL TROLLEY FOUR DRAWER

A high grade four drawer stainless steel trolley with smooth rolling castors, shelf, rails and an outstanding finish.

Features:

- High grade 304 stainless steel throughout
- Four drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors

Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SST04

[Click For More Information](#)

STAINLESS STEEL TROLLEY ONE DRAWER WITH LOCK

A high grade single stainless steel trolley with one lockable drawer. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- One lockable drawer
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- Top locking drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:
SST01LK

[Click For More Information](#)





STAINLESS STEEL TROLLEY TWO DRAWER WITH TOP LOCKING DRAWER

High grade single stainless steel two drawer trolley. Top drawer with lock. Has top and bottom side rails and locks on the front castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- Two drawer with top drawer lockable
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- Top locking drawer
- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE: SST02LK

[Click For More Information](#)

STAINLESS STEEL TROLLEY THREE DRAWER WITH TOP LOCKING DRAWER

A high grade single stainless steel three drawers trolley. Top drawer with lock. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

- High grade 304 stainless steel throughout
- Three drawer with top drawer lockable
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides



Outstanding finish

Specifications:

125mm deep drawer

Top locking drawer

95mm lockable castors

Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:

SST03LK

[Click For More Information](#)

**STAINLESS STEEL TROLLEY FOUR DRAWER WITH TOP LOCKING DRAWER**

A high grade single stainless steel four drawer trolley. Top drawer with lock. Has top and bottom side rails and locks on the castor wheels for safety.

Features:

High grade 304 stainless steel throughout
Four drawer with top drawer lockable
Bottom shelf
Double safety folds
Smooth welding
Soft closing magnetic drawer
German ball bearing drawer runners
Smooth running lockable castors
Top and bottom rails on three sides
Outstanding finish

Specifications:

125mm deep drawer

Top locking drawer

95mm lockable castors

Dimensions: 50 x 50 x 90cm (LxWxH)

ORDER CODE:

SST04LK

[Click For More Information](#)

DOUBLE TROLLEYS



STAINLESS STEEL DOUBLE TROLLEY NO DRAWER

A high grade double stainless steel trolley with no drawer. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE: SSTD00

[Click For More Information](#)

STAINLESS STEEL DOUBLE TROLLEY ONE DRAWER (FULL WIDTH)

A high grade double stainless steel trolley with one full width drawer. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- One drawer (full width)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides



Outstanding finish

Specifications:

125mm deep drawer

95mm lockable castors

Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:**SSTD01**

[Click For More Information](#)

**STAINLESS STEEL DOUBLE TROLLEY TWO DRAWER (FULL WIDTH)**

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

High grade 304 stainless steel throughout

Two drawer (full width)

Bottom shelf

Double safety folds

Smooth welding

Soft closing magnetic drawer

German ball bearing drawer runners

Smooth running lockable castors

Top and bottom rails on three sides

Outstanding finish

Specifications:

125mm deep drawer

95mm lockable castors

Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:**SSTD02**

[Click For More Information](#)

STAINLESS STEEL DOUBLE TROLLEY TWO DRAWER (SIDE BY SIDE)

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Two drawer (side by side)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD22

[Click For More Information](#)



STAINLESS STEEL DOUBLE TROLLEY FOUR DRAWER (TWO OVER TWO)

A high grade double four drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Four drawer (two over two)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

125mm deep drawer
95mm lockable castors
Dimensions: 80 x 50 x 90cm (LxWxH)

ORDER CODE:
SSTD04

[Click For More Information](#)

RESUS & ANAESTHESIA TROLLEYS



STAINLESS STEEL ANAESTHETIC TROLLEY FOUR DRAWER

A strong high quality anaesthetic trolley with two standard and two deep drawer with rails on the top and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Four drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top rails on three sides
- Outstanding finish

Specifications:

- Drawer Sizes: 2 x 125mm & 2 x 250mm
- 95mm lockable castors
- Dimensions: 60 x 50 x 107cm (LxWxH)

ORDER CODE: SSAT04

[Click For More Information](#)

STAINLESS STEEL ANAESTHETIC TROLLEY SIX DRAWER

A strong high quality six drawer anaesthetic trolley with two side shelves with rails. It has a large top surface with rails on three sides and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Six drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top rails on three sides
- Two side shelves with rails
- Outstanding finish

Specifications:

Drawer Sizes: 1 - 75mm, 2&3 - 150mm, 4 - 75mm, 5&6 - 100mm
95mm lockable castors
Dimensions: 90 x 50 x 97cm (LxWxH)

ORDER CODE:
SSAT06

[Click For More Information](#)



STAINLESS STEEL RESUSCITATION TROLLEY FIVE DRAWER

A strong high quality five drawer resuscitation trolley with rails on both shelves and lockable castors for safety.

Features:

- High grade 304 stainless steel throughout
- Five drawer
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top & bottom rails on three sides
- Bottom shelf
- Outstanding finish

Specifications:

Drawer Sizes: 4 x half, 1 x full width
95mm lockable castors
Dimensions: 80 x 50 x 97cm (LxWxH)

ORDER CODE:
SSRT05

[Click For More Information](#)

STAINLESS STEEL RESUSCITATION TROLLEY SEVEN DRAWER

A strong high quality seven drawer resuscitation trolley with rails on top and lockable castors for safety.

Features:

High grade 304 stainless steel throughout
Seven drawer
Double safety folds
Smooth welding
Soft closing magnetic drawer
German ball bearing drawer runners
Smooth running lockable castors
Top & bottom rails on three sides
Bottom shelf
Outstanding finish

Specifications:

Drawer Sizes: 6 x 125mm half width, 1 x 250mm full width
95mm lockable castors
Dimensions: 80 x 50 x 97cm (LxWxH)

ORDER CODE:
SSRT07

[Click For More Information](#)



IV STANDS



STAINLESS STEEL IV STAND TWO HOOK

A strong stainless steel IV stand two hook with wide adjustable height range, five leg weighted base and lockable wheels for easy mobility.

Specifications:

Two hook
Stainless steel construction
Five leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm
Base diameter 54.4cm

ORDER CODE:
IVP2H

[Click For More Information](#)

STAINLESS STEEL IV STAND FOUR HOOK

A strong four hook stainless steel IV stand with a five leg weighted base and easy mobility.

Specifications:

Four hook
Stainless steel construction
Five leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm
Base diameter 54.4cm

ORDER CODE:
IVP4H

[Click For More Information](#)



MAYO INSTRUMENT TROLLEYS



MAYO INSTRUMENT TROLLEY

Strong stainless steel mayo trolley with removeable tray, height adjustment and lockable castors for safety.

Specifications:

304 stainless steel construction
Removable tray: 60 x 40 x 2cm (LxWxD)
Telescopic height adjustment: 85 - 140cm
65mm lockable castors

ORDER CODE: MIT

[Click For More Information](#)

OXYGEN TROLLEYS



STAINLESS STEEL OXYGEN TROLLEY SINGLE DRAWER

A high grade single drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Single drawer
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE: SST01OXY

[Click For More Information](#)

STAINLESS STEEL OXYGEN TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Two drawers
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawers
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE:
SST02OXY

[Click For More Information](#)



OXYGEN TROLLEY SIZE D

A strong and portable oxygen bottle trolley to suit D size cylinders.

Specifications:

- Suits D size cylinder
- 304 stainless steel throughout
- securing chain
- rugged wheels
- Dimensions: 38 x 27 x 100cm

ORDER CODE:
OTD

[Click For More Information](#)

LINEN SKIP TROLLEYS & BAGS



SINGLE LINEN SKIP

Single stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 40 x 56 x 85cm

ORDER CODE:
SLT

[Click For More Information](#)

DOUBLE LINEN SKIP

Double stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 56 x 78 x 85cm

ORDER CODE:
DLT

[Click For More Information](#)





TRIPLE LINEN SKIP

Triple stainless steel linen skip with no lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 56 x 118 x 85cm

ORDER CODE:

TLT

[Click For More Information](#)

SINGLE LINEN SKIP WITH LID

Single stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 40 x 56 x 87cm

ORDER CODE:

SLTL

[Click For More Information](#)





DOUBLE LINEN SKIP WITH LID

Double stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 56 x 78 x 87cm

ORDER CODE:

DTL

[Click For More Information](#)

TRIPLE LINEN SKIP WITH LID

Triple stainless steel linen skip with foot operated acrylic lid and lockable castors for safety.

Specifications:

304 stainless steel throughout
80mm lockable castors
Bag Frame Outer Dimension: 40 x 35cm
Bag Frame Inner Dimension: 37 x 32cm
Overall Dimensions: 56 x 118 x 87cm

ORDER CODE:

TLTL

[Click For More Information](#)





LINEN SKIP BAG STANDARD

A strong linen skip bag constructed from purpose made 300-denier polyester twill fabric known as Durotech or Sharkskin. Fits most linen skips and comes in a range of colours.

Features:

- Constructed from purpose made 300-denier polyester twill fabric known as Durotech or Sharkskin
- Heavy Duty Drawcord
- Industrial strength cord lock
- Bottom handles inside and out
- Strong double sewn bottom seams
- Four high strength polyester tabs for hanging the bag on a frame
- Fits most linen skips
- Available in a range of colours

ORDER CODES:

SLBR - Red

SLBB - Blue

SLBW - White

[Click For More Information](#)

LINEN SKIP BAG SOLUBLE SEAM BX250

A polyethylene linen skip bag with a soluble seam to eliminate the handling of soiled or infected laundry.

Features:

- Constructed from polyethylene
- Dissolvable seam
- Detachable water soluble tie
- Eliminates handling of soiled or infected laundry
- Standard size 710x990mm
- Available in a range of colours

ORDER CODES:

SSBR - Red

SSBY - Yellow

SSBP - Purple

[Click For More Information](#)



INDEX

D

Double Linen Skip 23

Double Linen Skip With Lid 25

L

Linen Skip Bag Soluble Seam BX250 26

Linen Skip Bag Standard 26

M

Mayo Instrument Trolley 19

O

Oxygen Trolley Size D 21

S

Single Linen Skip 23

Single Linen Skip With Lid 24

Stainless Steel Anaesthetic Trolley Four Drawer 14

Stainless Steel Anaesthetic Trolley Six Drawer 15

Stainless Steel Double Trolley Four Drawer (Two Over Two) 12

Stainless Steel Double Trolley No Drawer 10

Stainless Steel Double Trolley One Drawer (Full Width) 10

Stainless Steel Double Trolley Two Drawer (Full Width) 11

Stainless Steel Double Trolley Two Drawer (Side By Side) 12

Stainless Steel IV Stand Four Hook 17

Stainless Steel IV Stand Two Hook 17

Stainless Steel Oxygen Trolley Single Drawer 20

Stainless Steel Oxygen Trolley Two Drawer 21

Stainless Steel Resuscitation Trolley Five Drawer 15

Stainless Steel Resuscitation Trolley Seven Drawer 16

Stainless Steel Trolley Four Drawer 6

Stainless Steel Trolley Four Drawer With Top Locking Drawer 9

Stainless Steel Trolley No Drawer 4

Stainless Steel Trolley No Rails 50x50x90 2

Stainless Steel Trolley No Rails 60x50x90 3

Stainless Steel Trolley No Rails 80x50x90 3

Stainless Steel Trolley One Drawer 4

Stainless Steel Trolley One Drawer With Lock 7

Stainless Steel Trolley Three Drawer 6

Stainless Steel Trolley Three Drawer With Top Locking Drawer 8

Stainless Steel Trolley Two Drawer 5

Stainless Steel Trolley Two Drawer With Top Locking Drawer 8

T

Triple Linen Skip 24

Triple Linen Skip With Lid 25



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au