



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Double Trolleys ..... 2**

# DOUBLE TROLLEYS



## STAINLESS STEEL DOUBLE TROLLEY NO DRAWER

A high grade double stainless steel trolley with no drawer. Has top and bottom side rails and lockable castors for safety.

### Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

### Specifications:

- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

### ORDER CODE:

**SSTD00**

[Click For More Information](#)

## STAINLESS STEEL DOUBLE TROLLEY ONE DRAWER (FULL WIDTH)

A high grade double stainless steel trolley with one full width drawer. Has top and bottom side rails and lockable castors for safety.

### Features:

- High grade 304 stainless steel throughout
- One drawer (full width)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides



Outstanding finish

**Specifications:**

125mm deep drawer

95mm lockable castors

Dimensions: 80 x 50 x 90cm (LxWxH)

**ORDER CODE:****SSTD01**

[Click For More Information](#)

**STAINLESS STEEL DOUBLE TROLLEY TWO DRAWER (FULL WIDTH)**

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

**Features:**

High grade 304 stainless steel throughout

Two drawer (full width)

Bottom shelf

Double safety folds

Smooth welding

Soft closing magnetic drawer

German ball bearing drawer runners

Smooth running lockable castors

Top and bottom rails on three sides

Outstanding finish

**Specifications:**

125mm deep drawer

95mm lockable castors

Dimensions: 80 x 50 x 90cm (LxWxH)

**ORDER CODE:****SSTD02**

[Click For More Information](#)

## STAINLESS STEEL DOUBLE TROLLEY TWO DRAWER (SIDE BY SIDE)

A high grade double two drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

### Features:

- High grade 304 stainless steel throughout
- Two drawer (side by side)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

### Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Dimensions: 80 x 50 x 90cm (LxWxH)

**ORDER CODE:**  
**SSTD22**

[Click For More Information](#)



## STAINLESS STEEL DOUBLE TROLLEY FOUR DRAWER (TWO OVER TWO)

A high grade double four drawer stainless steel trolley. Has top and bottom side rails and lockable castors for safety.

### Features:

- High grade 304 stainless steel throughout
- Four drawer (two over two)
- Bottom shelf
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

### Specifications:

125mm deep drawer  
95mm lockable castors  
Dimensions: 80 x 50 x 90cm (LxWxH)

**ORDER CODE:**  
**SSTD04**

[Click For More Information](#)

# INDEX

## S

Stainless Steel Double Trolley Four Drawer (Two Over Two) 4

Stainless Steel Double Trolley No Drawer 2

Stainless Steel Double Trolley One Drawer (Full Width) 2

Stainless Steel Double Trolley Two Drawer (Full Width) 3

Stainless Steel Double Trolley Two Drawer (Side By Side) 4



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)