



Pacific Medical Catalogue

TABLE OF CONTENTS

IV Stands 2

IV STANDS



STAINLESS STEEL IV STAND TWO HOOK

A strong stainless steel IV stand two hook with wide adjustable height range, five leg weighted base and lockable wheels for easy mobility.

Specifications:

Two hook
Stainless steel construction
Five leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm
Base diameter 54.4cm

ORDER CODE:
IVP2H

[Click For More Information](#)

STAINLESS STEEL IV STAND FOUR HOOK

A strong four hook stainless steel IV stand with a five leg weighted base and easy mobility.

Specifications:

Four hook
Stainless steel construction
Five leg weighted base
Smooth rolling and lockable castors
Height range 135-205cm
Base diameter 54.4cm

ORDER CODE:
IVP4H

[Click For More Information](#)



INDEX

S

Stainless Steel IV Stand Four Hook 2

Stainless Steel IV Stand Two Hook 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au