



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Mayo Instrument Trolleys ..... 2**

# MAYO INSTRUMENT TROLLEYS



## MAYO INSTRUMENT TROLLEY

Strong stainless steel mayo trolley with removeable tray, height adjustment and lockable castors for safety.

### Specifications:

304 stainless steel construction  
Removable tray: 60 x 40 x 2cm (LxWxD)  
Telescopic height adjustment: 85 - 140cm  
65mm lockable castors

### ORDER CODE:

**MIT**

[Click For More Information](#)

# INDEX

## M

Mayo Instrument Trolley 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)