



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Specialty Stainless Steel Trolleys** ..... 2

    No Rail Trolleys ..... 2

# SPECIALTY STAINLESS STEEL TROLLEYS

## NO RAIL TROLLEYS



### STAINLESS STEEL TROLLEY NO RAILS 50X50X90

A high grade stainless steel trolley 50cm wide with no draws or rails. Has lockable castors for safety.

#### Features:

- High grade 304 stainless steel throughout
- Bottom shelf
- Double safety folds
- Smooth welding
- Smooth running lockable castors
- Outstanding finish

#### Specifications:

- 95mm lockable castors
- Dimensions: 50 x 50 x 90cm (LxWxH)

#### ORDER CODE: **SSTNR50**

[Click For More Information](#)

## STAINLESS STEEL TROLLEY NO RAILS 60X50X90

A high grade stainless steel trolley 60cm wide with no draws or rails.  
Has lockable castors for safety.

### Features:

High grade 304 stainless steel throughout  
Bottom shelf  
Double safety folds  
Smooth welding  
Smooth running lockable castors  
Outstanding finish

### Specifications:

95mm lockable castors  
Dimensions: 60 x 50 x 90cm (LxWxH)

### ORDER CODE:

**SSTNR60**

[Click For More Information](#)



## STAINLESS STEEL TROLLEY NO RAILS 80X50X90

A high grade stainless steel trolley 80cm wide with no draws or rails.  
Has lockable castors for safety.

### Features:

High grade 304 stainless steel throughout  
Bottom shelf  
Double safety folds  
Smooth welding  
Smooth running lockable castors  
Outstanding finish

### Specifications:

95mm lockable castors  
Dimensions: 80 x 50 x 90cm (LxWxH)

### ORDER CODE:

**SSTNR80**

[Click For More Information](#)



# INDEX

## S

Stainless Steel Trolley No Rails 50x50x90 2

Stainless Steel Trolley No Rails 60x50x90 3

Stainless Steel Trolley No Rails 80x50x90 3



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)