



Pacific Medical Catalogue

TABLE OF CONTENTS

Oxygen Trolleys 2

OXYGEN TROLLEYS



STAINLESS STEEL OXYGEN TROLLEY SINGLE DRAWER

A high grade single drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Single drawer
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE: SST01OXY

[Click For More Information](#)

STAINLESS STEEL OXYGEN TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Two drawers
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawers
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE:
SST02OXY

[Click For More Information](#)



OXYGEN TROLLEY SIZE D

A strong and portable oxygen bottle trolley to suit D size cylinders.

Specifications:

- Suits D size cylinder
- 304 stainless steel throughout
- securing chain
- rugged wheels
- Dimensions: 38 x 27 x 100cm

ORDER CODE:
OTD

[Click For More Information](#)

INDEX

O

Oxygen Trolley Size D 3

S

Stainless Steel Oxygen Trolley Single Drawer 2

Stainless Steel Oxygen Trolley Two Drawer 3



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au