



Pacific Medical Catalogue

TABLE OF CONTENTS

Thermometry 2

THERMOMETRY

INFRARED FOREHEAD THERMOMETER



This infrared forehead thermometer allows for fast, easy and accurate measurement. With a measurement time of one second at a distance of up to 5cm it is suited for use in many clinical settings such as GP, Hospital, Workplace, First Aid and more.

Specifications:

Measurement range: Human body temperature mode: 32.0°C-43.0°C
Display Resolution: 0.1°C
Measurement place: Forehead
Accuracy: 32.0°C - 34.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$ / 35.0°C - 42.0°C, maximum permissible error $\pm 0.2^{\circ}\text{C}$ / 42.1°C - 42.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$
Display screen: LCD screen with 4-digit and special icons
Operating environment: 10°C-40°C
Max. Relative humidity: <85%
Storage environment: -20°C - +55°C
Battery: 2 x 1.5V AAA Batteries

Order Code: CYT1

[Click For More Information](#)

INDEX

I

Infrared Forehead Thermometer 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au