



Pacific Medical Catalogue

TABLE OF CONTENTS

Electric Couches 2

ELECTRIC COUCHES



STEALTH FLEXION DISTRACTION CHIROPRACTIC TABLE

Our new manipulation chiropractic table is packed full with features. It gives you a wide range of treatment options from drops, breakaway, flexion and distraction with functionally adjustable headpieces. All drops are hand activated with tension controls and the removeable flexion control T-Bar will allow you to move freely around the table.

Features:

- Electric Hi Lo operation by foot switch
- Electric flexion function operated by foot switch
- Abdominal breakaway section diversifying treatment options
- Height and width adjustable padded headrest
- Removable flexion control handle
- Extendable padded footrest
- Luxurious and comfortable table
- Includes patient straps
- Sleek dark frame two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 172cm
- Width: 56-80cm
- Headrest height adjustment: 10cm
- Four tension controlled drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical to +20 Degrees
- Adjustable angle of thoracic from -15 to +18 Degrees
- Adjustable angle of pelvic (tilting up and down) -22 to +22 Degrees
- Adjustable angle of pelvic (left to right) 20 Degrees
- Height range from 56cm to 80cm
- 5cm wheels
- Padding thickness 7cm
- Foam Density: 100kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CTFDST - Stealth Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a

regular basis.

[Click For More Information](#)

STEALTH MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





STEALTH PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES: ET3BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

PROCEDURE CHAIR

Our procedure chair is a versatile multi modality product utilising fully electric functionality and a strong column lift for the rigours of modern day use.

Features:

- Fully electric operation by hand switch
- Electric seat tilt
- Strong column lift
- Small footprint frame with rounded frame legs for safety
- Gas strut adjustable face crest
- Padded armrests
- Debris tray
- Contoured backrest
- Can achieve Trendelenburg position
- 10cm central lock castors
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 198cm
- Width: 64cm
- Headrest length: 31cm
- Backrest length: 61cm
- Seat section length: 46cm
- Leg section length: 48cm
- Adjustable motorised back section from 0 to +85 Degrees
- Adjustable motorised foot section from 0 to -85 Degrees
- Adjustable motorised seat tilt from -5 to +35 Degrees
- Adjustable head rest through gas operated struts from -15 to +35 Degrees
- Height range from 59cm to 89cm in approx 30 seconds
- Max working heel height 145cm (with footrest extension and extended)
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PRCHNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





TWO SECTION TREATMENT COUCH

One of our most popular tables. This strong two section table will withstand the rigours of patient treatment whether it is used in a demanding Physiotherapy clinic or a small GP practice you can rest assured it will perform. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 67cm
- Leg section length: 123cm
- Adjustable head section through dual gas operated struts from -20 to +75 Degrees
- Height range from 47cm to 92cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

- ET2NB - Navy Blue**
- ET2BL - Black**
- ET2GR - Grey**

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33NB - Navy Blue

ET33BL - Black

ET33GR - Grey

ET33CG - Charcoal Grey

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





THREE SECTION ALL ELECTRIC TREATMENT COUCH

A high quality all electric treatment couch facilitating effective patient treatment all whilst looking after practitioner OH&S.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest & footrest for effortless operation
- Sleek new design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motors
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees
- Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees
- Height range from 47cm to 91cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33AENB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION BARIATRIC TREATMENT COUCH

Every practice should have at least one heavy duty couch to lift bariatric patients. This three equal section couch with electric backrest and dual gas strut supported footrest boasts a large lifting capacity of 280KG to perform these duties with ease.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head section through electric motor & foot switch from to -40 to +65 Degrees
- Adjustable foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 280KG

ORDER CODES:

ET32NB - Navy Blue

ET32BL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





THREE SECTION PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by foot switch or footbar
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET3NB - Navy Blue

ET3BL - Black

ET3NBSB - Navy Blue (With Surround Bar)

ET3BLSB - Black (With Surround Bar)

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

PODIATRY / MULTI PURPOSE CHAIR

A versatile podiatry and multipurpose treatment chair. This strong chair includes many features such as nose hole and plug, flip up armrests & a debris tray to suit the many disciplines of the healthcare and beauty industry.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable flip up armrests
- Padded adjustable foot rest
- Head section has face hole and headrest plug
- Debris tray included
- Adjustable feet for uneven surfaces
- 5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 220cm (Flat Table & Including Footrest)
- Width: 91cm (Including Armrests)
- Backrest length: 91cm
- Mid section length: 56cm
- Leg section length: 38cm
- Adjustable backrest through dual gas operated struts to +80 Degrees
- Motorised seat tilt to +25 Degrees
- Height range from 56.5cm to 98cm in approx 30 seconds
- Max working heel height: 165cm (with footrest extension and extended)
- Padding thickness: 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PCNB - Navy Blue

PCBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





GYNAECOLOGY TREATMENT COUCH

Everything you will need for women's health and general examinations all in one table. This flexible table includes foot pads, stirrups, debris tray, a fill section & nose hole with plug allowing for a multitude of different examinations & treatments. It features a smooth hi-lo adjustment & electric seat tilt with foot control for optimum manoeuvrability and better patient treatment.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Footpads & stirrups both included
- Debris tray
- Side rails
- Fill section allowing general use
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 185cm
- Width: 60cm
- Backrest length: 91cm
- Mid section length: 46.5cm
- Leg section length: 42.5cm
- Adjustable backrest through dual gas operated struts to +65 degrees
- Adjustable motorised seat angle from 0 to +20 degrees
- Height range from 64cm to 106cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

GCNB - Navy Blue

GCBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

GYNAECOLOGY PREMIUM TREATMENT COUCH

Our premium gynaecology couch provides for electric positioning of height & tilt through a foot control ensuring practitioner comfort. This couch has been designed with an additional headrest section transforming to a full length couch for most general examinations & minor treatments.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Electric backrest operation by foot switch
- Dual gas strut operated head rest
- Debris Tray
- Foot Pads
- Stirrups
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 181cm
- Width: 70cm
- Headrest length: 61cm
- Backrest length: 76cm
- Seat section length: 40cm
- Adjustable headrest through dual gas operated struts from -30 to +47 Degrees
- Adjustable motorised backrest from 0 to +71 Degrees
- Adjustable motorised seat tilt up to +18 Degrees
- Height range from 57cm to 118cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

GCPRNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



CHIROPRACTIC TABLE



This chiropractic drop table is a versatile cost effective table providing treatment flexibility for your practice. All Cervical, Thoracic, Lumbar & Pelvic drops by hand activation & tensioning.

Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Height adjustable & extendable padded footrest
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 180-193cm
- Width: 55-59cm
- Headrest length: 40cm
- Tapered sections length: 20.5cm & 22.5cm
- Long section length: 69.5cm
- Footrest Length: 13cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable long board to +15 Degrees
- Height adjustable headrest from +25 to -20 degrees
- Height range from 45cm to 74cm in approx 30 seconds
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

CTNB - Navy Blue

CTBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM TABLE

This Chiropractic table is our new premium table. Boasting a premium column lift eliminating any lifting mechanism. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. A premium quality & luxury product.



Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Luxurious and comfortable table
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 178cm
- Width: 54-65cm
- Headrest height adjustment: 10cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical from +25 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Adjustable angle of seat rest to +10 Degrees
- Adjustable angle of backrest to +25 Degrees
- Height range from 60cm to 80cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

CTPRNB - Navy Blue

CTPRBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM FIXED HEIGHT TABLE



This Chiropractic table is our new premium fixed height table. Very similar to our premium column but is fixed in height. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. Overall great value for any clinic.

Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 173cm
- Width: 62cm
- Headrest length: 40cm
- Adjustable angle of cervical from +15 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Fixed height of 55cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

CTFIXNB - Navy Blue

CTFIXBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH NO POSTURAL DRAINAGE

A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.



Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET51NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH WITH POSTURAL DRAINAGE



A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo and postural drainage operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.

Features:

- Electric Hi Lo operation by foot switch
- Electric postural drainage operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Motorised postural drainage adjustment to +25 Degrees
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET52NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

TWO SECTION CARDIOLOGY COUCH

A strong two section cardiology couch with comfortable sewn edge upholstery, central lock castors, electric backrest & two special cardiac cut outs allowing increased patient access. With a sleek design, great features and foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Dual cut outs for easy patient access
- Comfortable sewn edge upholstery
- Sleek design with rounded frame legs for safety
- 10cm central lock castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 67cm
- Leg section length: 123cm
- Cutout Size: 31cm x 25cm
- Cutout Mechanism: Manual opening and held in place mechanically, very strong.
- Adjustable head section through electric motor & foot switch from flat to +75 Degrees
- Height range from 58cm to 102cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CARNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



NEUROLOGICAL BOBATH TREATMENT COUCH



Designed for rehabilitation this strong two section neurological bobath has been designed with the therapist in mind. The large table top provides plenty of space for patient and therapist and three nose holes and plugs provides flexibility of treatment.

Features:

- Electric Hi Lo operation by foot bar
- Sleek design with rounded frame legs for safety
- Head section has three face holes and plugs
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 190cm
- Width: 120cm
- Backrest length: 65cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +65 Degrees
- Height range from 49cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

NT120NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CONTOURED MASSAGE COUCH

A comfortable and adjustable massage couch with contoured surface to increase patient access. The flexible padded face rest and prone arm rest ensure patient comfort and relaxation during treatment.



Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable head cushion with prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 32cm
- Long section length: 165cm
- Adjustable headrest through gas operated strut from +30 to -75 Degrees
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CMNB - Navy Blue

CMBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



TILT TABLE TWO SECTION

Our tilt table offers a simple out of the box solution for physiotherapists and in home care patients alike. This two section tilt table including patient straps, activity tray and individually adjustable foot plates will suit the most demanding rehabilitation regime and perform for the long term.

Features:

- Electric hi-lo & tilt operation by hand switch
- Tilt angle inclinometer included
- Activity tray & patient straps included
- Narrowed leg section with strap rail in centre
- Head section has face hole and plug
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 194cm
- Width: 64cm
- Backrest length: 66.5cm
- Long section length: 124cm
- Support leg width: 87cm
- Adjustable backrest through gas operated strut to +75 Degrees
- Adjustable individual foot plates from +20 to -10 Degrees
- Motorised tilt to +85 Degrees
- Height range from 59cm to 104cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

TTNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FOUR SECTION TRACTION TREATMENT COUCH

This traction couch offers a great solution for static and intermittent lumbar and cervical traction treatments. Including a flexion stool and angle adjustable machine mount this couch is an out of the box solution for many practitioners.



Features:

- Electric Hi Lo operation by foot bar
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Upholstered traction stool
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 40cm
- Long seat length: 39cm
- Short seat length: 29cm
- Backrest length: 81cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from +40 to -45 Degrees
- Max traction of 16cm between second and third sections
- Height range from 52cm to 94cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

TRACNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



BEAUTY DAY SPA COUCH

This electric hi lo couch is perfect for use in beauty, day spa and massage. It can be used in numerous comfortable positions with the head and leg section adjusted by gas struts. It has removeable padded arm rests on each side that are also height adjustable. The wings next to the head piece are removeable for patient comfort in massage allowing them to rest their arms on the prone arm rest.

Features:

- Electric Hi Lo operation by foot switch
- Luxurious soft and comfortable head cushion
- Contoured comfortable table top
- Sleek design with rounded frame legs for safety
- Removeable padded wing armrests
- Padded prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 186cm
- Width: 62cm
- Backrest length: 58cm
- Mid section length: 45cm
- Legrest length: 41cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +85 Degrees
- Adjustable footrest through gas operated strut from +5 to -75 Degrees
- Height range from 57cm to 98cm in approx 30 seconds
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

BTWH - White

BTNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION MEDICAL TREATMENT COUCH WITH SIDE RAILS

Our most popular product this three equal section table with side rails & smooth hi-lo operation represents great value.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

SRET33NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



SINGLE SECTION TREATMENT / CHANGE COUCH WITH SIDE RAILS



This strong single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Padded side rails
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 180cm
- Width: 66cm
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CHNBSB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

VETERINARY LIFT TABLE

Features:

Electric Hi Lo operation by hand switch & footbar
Durable 304 stainless steel patient surface
Sleek new design with rounded frame legs for safety
10cm lockable castors for safety and mobility
Reliable TiMotion motors
Rechargeable battery

Specification:

Total Length: 110cm
Width: 56cm
Height range from 41cm to 96cm in approx 30 seconds
Warranty: 7 Years Construction, 3 Years Electrics
Certification: CE, FDA, ISO13485-2016
Weight Capacity: 250KG

ORDER CODES:

ETVET

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



INDEX

B

Beauty Day Spa Couch 24

C

Chiropractic Premium Fixed Height Table 16

Chiropractic Premium Table 15

Chiropractic Table 14

Contoured Massage Couch 21

F

Five Section Treatment Couch No Postural Drainage 17

Five Section Treatment Couch With Postural Drainage 18

Four Section Traction Treatment Couch 23

G

Gynaecology Premium Treatment Couch 13

Gynaecology Treatment Couch 12

N

Neurological Bobath Treatment Couch 20

P

Podiatry / Multi Purpose Chair 11

Procedure Chair 5

S

Single Section Treatment / Change Couch With Side Rails 26

Stealth Flexion Distraction Chiropractic Table 2

Stealth Medical Treatment Couch 3

Stealth Physio Treatment Couch 4

T

Three Section All Electric Treatment Couch 8

Three Section Bariatric Treatment Couch 9

Three Section Medical Treatment Couch 7

Three Section Medical Treatment Couch With Side Rails 25

Three Section Physio Treatment Couch 10

Tilt Table Two Section 22

Two Section Cardiology Couch 19

Two Section Treatment Couch 6

V

Veterinary Lift Table 27



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au