



Pacific Medical Catalogue

TABLE OF CONTENTS

Treatment Couches 2

 Electric Couches 2

 Fixed Height Exam Couches 30

 Portable Massage Couches 32

 Bedsheets 33

 Treatment Couch Accessories & Spare Parts 35

TREATMENT COUCHES

ELECTRIC COUCHES



BARIATRIC THREE SECTION ALL ELECTRIC TREATMENT COUCH

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest & footrest for effortless operation
- Sleek new frame design
- Head section has face hole and plug
- 12.5cm central lock castors for safety and mobility
- Reliable TiMotion motors
- Durable thick foam padding

Specification:

- Total Length: 195cm
- Width: 80cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees
- Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees
- Height range from 52cm to 93cm in approx 30 seconds
- Hoist Clearance: 15cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 350KG

ORDER CODES:

ET33BARIANB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

STEALTH FLEXION DISTRACTION CHIROPRACTIC TABLE (REPLACED BY FLEXION DISTRACTION PRO SERIES)

Our new manipulation chiropractic table is packed full with features. It gives you a wide range of treatment options from drops, breakaway, flexion and distraction with functionally adjustable headpieces. All drops are hand activated with tension controls and the removeable flexion control T-Bar will allow you to move freely around the table.



Features:

- Electric Hi Lo operation by foot switch
- Electric flexion function operated by foot switch
- Abdominal breakaway section diversifying treatment options
- Height and width adjustable padded headrest
- Removable flexion control handle
- Extendable padded footrest
- Luxurious and comfortable table
- Includes patient straps
- Sleek dark frame two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 172cm
- Width: 56-80cm
- Headrest height adjustment: 10cm
- Four tension controlled drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical to +20 Degrees
- Adjustable angle of thoracic from -15 to +18 Degrees
- Adjustable angle of pelvic (tilting up and down) -22 to +22 Degrees
- Adjustable angle of pelvic (left to right) 20 Degrees
- Height range from 56cm to 80cm
- 5cm wheels
- Padding thickness 7cm
- Foam Density: 100kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CTFDST - Stealth Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

STEALTH FLEXION DISTRACTION CHIROPRACTIC TABLE PRO SERIES



Our new manipulation chiropractic table is packed full with features. It gives you a wide range of treatment options from drops, breakaway, flexion, distraction & electronic traction with functionally adjustable headpieces. All drops are hand activated with tension controls along with the foot bar hi-lo and removeable flexion control T-Bar will allow you to move freely around the table.

Features:

- Electric Hi Lo operation by foot bar
- Electric traction function operated by hand button
- Abdominal breakaway section diversifying treatment options
- Height and width adjustable padded headrest
- Removable flexion control handle
- Extendable padded footrest
- Luxurious and comfortable table
- Includes patient straps
- Sleek dark frame with retractable 75mm wheel system
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 172cm
- Width: 56-80cm
- Headrest height adjustment: 10cm
- Headrest tilt: +-20 Degrees
- Headrest lateral movement: +-15 Degrees
- Four tension controlled drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical to +20 Degrees
- Adjustable angle of thoracic from -13 to +15 Degrees
- Adjustable angle of pelvic (tilting up and down) -40 to +25 Degrees
- Adjustable angle of pelvic (left to right) 20 Degrees
- Height range from 56cm to 90cm
- 7.5cm wheels
- Padding thickness 7.5cm
- Foam Density: 100kg/m3
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE-MDR, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CTFDSTPRO - Stealth Black PRO Series

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of

pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

STEALTH MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm centrally retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





STEALTH PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by 360 degree foot bar
- Sleek black frame with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm centrally retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET3BLST - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

PROCEDURE CHAIR (REPLACED BY PROCEDURE CHAIR PRO SERIES)

Our procedure chair is a versatile multi modality product utilising fully electric functionality and a strong column lift for the rigours of modern day use.

PLEASE NOTE - REPLACED BY PROCEDURE CHAIR PRO SERIES



Features:

- Fully electric operation by hand switch
- Electric seat tilt
- Strong column lift
- Small footprint frame with rounded frame legs for safety
- Gas strut adjustable face crest
- Padded armrests
- Debris tray
- Contoured backrest
- Can achieve Trendelenburg position
- 10cm central lock castors
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 198cm
- Width: 64cm
- Headrest length: 31cm
- Backrest length: 61cm
- Seat section length: 46cm
- Leg section length: 48cm
- Adjustable motorised back section from 0 to +85 Degrees
- Adjustable motorised foot section from 0 to -85 Degrees
- Adjustable motorised seat tilt from -5 to +35 Degrees
- Adjustable head rest through gas operated struts from -15 to +35 Degrees
- Height range from 59cm to 89cm in approx 30 seconds
- Max working heel height 145cm (with footrest extension and extended)
- Padding thickness 7.5cm
- Foam Density: 36kg/m3
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PRCHNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



PROCEDURE CHAIR PRO SERIES

Our procedure chair pro series is a versatile multi modality product utilising fully electric functionality and a strong column lift for the rigours of modern day use.

Features:

- Fully electric operation by hand switch
- Electric seat tilt
- Strong column lift
- Small footprint frame with rounded frame legs for safety
- Gas strut adjustable face crest
- Height adjustable padded armrests
- Integrated slide out debris tray
- Contoured backrest
- Integrated paper towel holder
- Backup battery
- Can achieve Trendelenburg position
- 12.5cm central lock castors
- Reliable TiMotion motor
- Comfortable foam padding
- Optional accessories: IV Pole, U-Shaped Metal Head Support, Side Rails, Stirrups, Operating Battery

Specification:

- Total Length: 198cm
- Width: 64cm Patient Surface / 90cm At Armrests
- Headrest length: 33cm
- Backrest length: 61cm
- Seat section length: 45cm
- Leg section length: 48cm
- Adjustable motorised back section from -5 to +85 Degrees
- Adjustable motorised foot section from 0 to -85 Degrees
- Adjustable motorised seat tilt from 0 to +35 Degrees
- Adjustable head rest through gas operated struts from -15 to +35 Degrees
- Armrest height adjustment: 11cm
- Height range from 59cm to 89cm in approx 30 seconds
- Max working heel height 145cm (with footrest extension and extended)
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Optional accessories: IV Pole, U-Shaped Metal Head Support, Side Rails, Stirrups, Operating Battery
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PRCHNBPRO - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

TWO SECTION TREATMENT COUCH

One of our most popular tables. This strong two section table will withstand the rigours of patient treatment whether it is used in a demanding Physiotherapy clinic or a small GP practice you can rest assured it will perform. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 67cm
- Leg section length: 123cm
- Adjustable head section through dual gas operated struts from -20 to +75 Degrees
- Height range from 47cm to 92cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET2NB - Navy Blue

ET2BL - Black

ET2GR - Grey

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





THREE SECTION MEDICAL TREATMENT COUCH

Our most popular product this three equal section couch with dual gas strut supported sections & smooth hi-lo operation represents great value. Australia's number 1 choice for medical centres for many years.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

- ET33NB - Navy Blue**
- ET33BL - Black**
- ET33GR - Grey**
- ET33CG - Charcoal Grey**

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION ALL ELECTRIC TREATMENT COUCH

A high quality all electric treatment couch facilitating effective patient treatment all whilst looking after practitioner OH&S.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest & footrest for effortless operation
- Sleek new design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motors
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head section through electric motor & foot switch from to -40 to +75 Degrees
- Adjustable foot section through electric motor & foot switch from to -40 to +75 Degrees
- Height range from 47cm to 91cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET33AENB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





THREE SECTION PHYSIO TREATMENT COUCH

The preferred choice of physiotherapists this short head couch provides a strong platform for treatments. With a gas strut supported head and backrest coupled with a smooth hi-lo operation this table represents great value.

Features:

- Electric Hi Lo operation by foot switch or footbar
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 40cm
- Mid section length: 54cm
- Leg section length: 93cm
- Adjustable back section through dual gas operated struts from to +80 Degrees
- Adjustable head section through gas operated struts from -45 to +40 Degrees
- Height range from 47.5cm to 93cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET3NB - Navy Blue

ET3BL - Black

ET3NBSB - Navy Blue (With Surround Bar)

ET3BLSB - Black (With Surround Bar)

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

PODIATRY / MULTI PURPOSE CHAIR

A versatile podiatry and multipurpose treatment chair. This strong chair includes many features such as nose hole and plug, flip up armrests & a debris tray to suit the many disciplines of the healthcare and beauty industry.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable flip up armrests
- Padded adjustable foot rest
- Head section has face hole and headrest plug
- Debris tray included
- Adjustable feet for uneven surfaces
- 5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 220cm (Flat Table & Including Footrest)
- Width: 91cm (Including Armrests)
- Backrest length: 91cm
- Mid section length: 56cm
- Leg section length: 38cm
- Adjustable backrest through dual gas operated struts to +80 Degrees
- Motorised seat tilt to +25 Degrees
- Height range from 56.5cm to 98cm in approx 30 seconds
- Max working heel height: 165cm (with footrest extension and extended)
- Padding thickness: 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

PCNB - Navy Blue

PCBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)





GYNAECOLOGY TREATMENT COUCH

Everything you will need for women's health and general examinations all in one table. This flexible table includes foot pads, stirrups, debris tray, a fill section & nose hole with plug allowing for a multitude of different examinations & treatments. It features a smooth hi-lo adjustment & electric seat tilt with foot control for optimum manoeuvrability and better patient treatment.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Sleek design with rounded frame legs for safety
- Footpads & stirrups both included
- Debris tray
- Side rails
- Fill section allowing general use
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 185cm
- Width: 60cm
- Backrest length: 91cm
- Mid section length: 46.5cm
- Leg section length: 42.5cm
- Adjustable backrest through dual gas operated struts to +65 degrees
- Adjustable motorised seat angle from 0 to +20 degrees
- Height range from 64cm to 106cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

GCNB - Navy Blue

GCBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

GYNAECOLOGY PREMIUM TREATMENT COUCH

Our premium gynaecology couch provides for electric positioning of height & tilt through a foot control ensuring practitioner comfort. This couch has been designed with an additional headrest section transforming to a full length couch for most general examinations & minor treatments.

Features:

- Electric Hi Lo operation by foot switch
- Electric seat tilt operation by foot switch
- Electric backrest operation by foot switch
- Dual gas strut operated head rest
- Debris Tray
- Foot Pads
- Stirrups
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 181cm
- Width: 70cm
- Headrest length: 61cm
- Backrest length: 76cm
- Seat section length: 40cm
- Adjustable headrest through dual gas operated struts from -30 to +47 Degrees
- Adjustable motorised backrest from 0 to +71 Degrees
- Adjustable motorised seat tilt up to +18 Degrees
- Height range from 57cm to 118cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

GCPRNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



CHIROPRACTIC TABLE PRO SERIES



PRO series chiropractic drop table is a versatile cost effective product providing treatment flexibility for your patients. With a foot activation bar for hi-lo and improved Cervical, Thoracic, Lumbar & Pelvic drops with hand activation & tensioning.

Features:

- Electric Hi Lo operation by foot bar
- Height and angle adjustable padded headrest
- Height adjustable & extendable padded footrest
- Sleek design with rounded frame legs and base frame cover plate
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 185-200cm
- Width: 55-65cm
- Headrest length: 35cm
- Tapered section length: 26cm
- Mid section length: 23cm
- Long section length: 70cm
- Footrest Length: 18cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable long board to +20 Degrees
- Height adjustable headrest from +20 to -20 degrees
- Height range from 46cm to 80cm in approx 30 seconds
- 7.5cm wheels
- Padding thickness 7.5cm
- Foam Density: 100kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CTNBPRO - Navy Blue

CTBLPRO - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM TABLE

This Chiropractic table is our new premium table. Boasting a premium column lift eliminating any lifting mechanism. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. A premium quality & luxury product.



Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Luxurious and comfortable table
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 178cm
- Width: 54-65cm
- Headrest height adjustment: 10cm
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Adjustable angle of cervical from +25 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Adjustable angle of seat rest to +10 Degrees
- Adjustable angle of backrest to +25 Degrees
- Height range from 60cm to 80cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

CTPRNB - Navy Blue

CTPRBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CHIROPRACTIC PREMIUM FIXED HEIGHT TABLE



This Chiropractic table is our new premium fixed height table. Very similar to our premium column but is fixed in height. All Cervical, Thoracic, Lumbar and Pelvic drops are operated by foot activation with the tensioning done by hand. This allows easy activation for the practitioner. Overall great value for any clinic.

Features:

- Electric Hi Lo operation by foot switch
- Height and width adjustable padded headrest
- Extendable padded footrest
- Sleek design with two rounded frame legs and two fixed wheels
- Adjustable feet for uneven surfaces
- Reliable TiMotion motor
- Thick foam padding

Specification:

- Total Length: 173cm
- Width: 62cm
- Headrest length: 40cm
- Adjustable angle of cervical from +15 to -20 Degrees
- Adjustable angle of thoracic to +25 Degrees
- Adjustable angle of lumbar to +30 Degrees
- Adjustable angle of pelvic to +30 Degrees
- Four drop sections: Cervical, Thoracic, Lumbar, Pelvic
- Fixed height of 55cm
- 5cm wheels
- Padding thickness 8cm
- Foam Density: 110g/m²
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

CTFIXNB - Navy Blue

CTFIXBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH NO POSTURAL DRAINAGE

A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.



Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET51NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FIVE SECTION TREATMENT COUCH WITH POSTURAL DRAINAGE



A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo and postural drainage operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.

Features:

- Electric Hi Lo operation by foot switch
- Electric postural drainage operation by foot switch
- Sleek design with rounded frame legs for safety
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm to 96cm in approx 30 seconds
- Motorised postural drainage adjustment to +25 Degrees
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ET52NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

TWO SECTION CARDIOLOGY COUCH

A strong two section cardiology couch with comfortable sewn edge upholstery, central lock castors, electric backrest & two special cardiac cut outs allowing increased patient access. With a sleek design, great features and foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot switch
- Electric backrest for effortless operation
- Dual cut outs for easy patient access
- Comfortable sewn edge upholstery
- Sleek design with rounded frame legs for safety
- 10cm central lock castors for safety and mobility
- Reliable TiMotion motor
- Durable foam padding

Specification:

- Total Length: 195cm
- Width: 66cm
- Head section length: 67cm
- Leg section length: 123cm
- Cutout Size: 31cm x 25cm
- Cutout Mechanism: Manual opening and held in place mechanically, very strong.
- Adjustable head section through electric motor & foot switch from flat to +75 Degrees
- Height range from 58cm to 102cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CARNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



NEUROLOGICAL BOBATH TREATMENT COUCH



Designed for rehabilitation this strong two section neurological bobath has been designed with the therapist in mind. The large table top provides plenty of space for patient and therapist and three nose holes and plugs provides flexibility of treatment.

Features:

- Electric Hi Lo operation by foot bar
- Sleek design with rounded frame legs for safety
- Head section has three face holes and plugs
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 190cm
- Width: 120cm
- Backrest length: 65cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +65 Degrees
- Height range from 49cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

NT120NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

CONTOURED MASSAGE COUCH

A comfortable and adjustable massage couch with contoured surface to increase patient access. The flexible padded face rest and prone arm rest ensure patient comfort and relaxation during treatment.



Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Comfortable head cushion with prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 32cm
- Long section length: 165cm
- Adjustable headrest through gas operated strut from +30 to -75 Degrees
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CMNB - Navy Blue

CMBL - Black

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



TILT TABLE TWO SECTION

Our tilt table offers a simple out of the box solution for physiotherapists and in home care patients alike. This two section tilt table including patient straps, activity tray and individually adjustable foot plates will suit the most demanding rehabilitation regime and perform for the long term.

Features:

- Electric hi-lo & tilt operation by hand switch
- Tilt angle inclinometer included
- Activity tray & patient straps included
- Narrowed leg section with strap rail in centre
- Head section has face hole and plug
- Sleek design with rounded frame legs for safety
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 194cm
- Width: 64cm
- Backrest length: 66.5cm
- Long section length: 124cm
- Support leg width: 87cm
- Adjustable backrest through gas operated strut to +75 Degrees
- Adjustable individual foot plates from +20 to -10 Degrees
- Motorised tilt to +85 Degrees
- Height range from 59cm to 104cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 225KG

ORDER CODES:

TTNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

FOUR SECTION TRACTION TREATMENT COUCH

This traction couch offers a great solution for static and intermittent lumbar and cervical traction treatments. Including a flexion stool and angle adjustable machine mount this couch is an out of the box solution for many practitioners.



Features:

- Electric Hi Lo operation by foot bar
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Upholstered traction stool
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 70cm
- Headrest length: 40cm
- Long seat length: 39cm
- Short seat length: 29cm
- Backrest length: 81cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from +40 to -45 Degrees
- Max traction of 16cm between second and third sections
- Height range from 52cm to 94cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

TRACNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



BEAUTY DAY SPA COUCH

This electric hi lo couch is perfect for use in beauty, day spa and massage. It can be used in numerous comfortable positions with the head and leg section adjusted by gas struts. It has removeable padded arm rests on each side that are also height adjustable. The wings next to the head piece are removeable for patient comfort in massage allowing them to rest their arms on the prone arm rest.

Features:

- Electric Hi Lo operation by foot switch
- Luxurious soft and comfortable head cushion
- Contoured comfortable table top
- Sleek design with rounded frame legs for safety
- Removeable padded wing armrests
- Padded prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 62cm (Not Including Armrests)
- Backrest length: 58cm
- Mid section length: 45cm
- Legrest length: 41cm
- Long section length: 121cm
- Adjustable backrest through gas operated strut to +85 Degrees
- Adjustable footrest through gas operated strut from +5 to -75 Degrees
- Height range from 57cm to 98cm in approx 30 seconds
- Padding thickness 7.5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

BTWH - White

BTNB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

THREE SECTION MEDICAL TREATMENT COUCH WITH SIDE RAILS

Our most popular product this three equal section table with side rails & smooth hi-lo operation represents great value.

Features:

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Three equal sections for universal use in any room
- Head section has face hole and plug
- Dual gas strut supported back & footrest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 63cm
- Mid section length: 63cm
- Leg section length: 63cm
- Adjustable head & foot section through dual gas operated struts from to -40 to +65 Degrees
- Height range from 46.5cm to 92.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

SRET33NB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



SINGLE SECTION TREATMENT / CHANGE COUCH WITH SIDE RAILS



This strong single section couch will withstand the rigours of in clinic patient treatment, in a special needs facility or home as a change couch. It has retractable side rails for added patient safety. With a sleek design, great features and a foot switch as standard it is simple to operate straight from the box without any assembly required.

Features:

- Electric Hi Lo operation by foot surround bar
- Sleek design with rounded frame legs for safety
- Padded side rails
- Face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Comfortable foam padding

Specification:

- Total Length: 180cm
- Width: 66cm
- Height range from 48cm to 93.5cm in approx 30 seconds
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

CHNBSB - Navy Blue

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)

VETERINARY LIFT TABLE

A robust and mobile animal lift table with stainless steel surface. Suitable for exam and treatment with battery operation to suit any size clinic.

Features:

- Electric Hi Lo operation by hand switch & footbar
- Durable 304 stainless steel patient surface
- Sleek new design with rounded frame legs for safety
- 10cm lockable castors for safety and mobility
- Reliable TiMotion motors
- Rechargeable battery

Specification:

- Total Length: 110cm
- Width: 56cm
- Height range from 41cm to 96cm in approx 30 seconds
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:

ETVET

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



FIXED HEIGHT EXAM COUCHES



FREE STANDING TREATMENT COUCH

Our free standing treatment couch is the workhorse of any clinic. Flat packed with a two minute assembly and available in a number of colours choices with and without face hole this couch is a popular choice.

Features:

- Robust
- New lift and lock mechanism for easy 2 minute assembly
- Lift up backrest for assessment
- Face hole and non-face hole models available
- Flat packed to save on freight

Specification:

- Total length: 195cm
- Width: 71cm
- Height: 82cm
- Weight capacity: 280KG
- Padding thickness: 5cm
- Foam density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007

ORDER CODES:

- FSBL - BLACK**
- FSBLF - BLACK WITH FACEHOLE**
- FSGR - GREY**
- FSGRF - GREY WITH FACEHOLE**
- FSNB - NAVY BLUE**
- FSNBF - NAVY BLUE WITH FACEHOLE**
- FSNB - WHITE**
- FSNBF - WHITE WITH FACEHOLE**

[Click For More Information](#)

FREE STANDING TREATMENT COUCH WITH DRAWERS

Our three drawer free standing treatment couch is the workhorse of any clinic. Flat packed with a simple & quick assembly this couch is a popular choice.

Features:

- Robust
- New lift and lock mechanism for easy 2 minute assembly
- Lift up backrest for assessment
- Non-face hole model only
- Three drawer with simple assembly
- Flat packed to save on freight

Specification:

- Total Length: 195cm
- Width: 71cm
- Height: 82cm
- Drawer: 3
- Weight Capacity: 280KG
- Padding thickness 5cm
- Foam Density: 36kg/m³
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Certification: CE, FDA, ISO9001-2008, ISO13485-2007

ORDER CODES:

FSBLDR - BLACK

FSGRDR - GREY

FSNBDR - NAVY BLUE

[Click For More Information](#)



PORTABLE MASSAGE COUCHES

PORTABLE MASSAGE COUCH



Aluminium framed portable massage couch with an auto lock cabling system, rounded corners and a multi density 5cm foam system. A great choice for mobile practitioners.

Features:

- Aluminium frame and legs
- Manufactured without electric welding to join the aluminium
- Face crest & prone arm rest
- Auto lock cabling system
- Rounded corners and deluxe no noise handle
- Luxurious PU vinyl
- Adjustable backrest
- Multi density 5cm foam system

Specifications:

- Height Adjustable: 65-80cm
- Length: 185cm
- Width: 66cm
- Folded Size: 93x67x17cm
- Weight: 18KG
- Weight Capacity: 220KG
- Warranty: 2 Years

ORDER CODE:

PMTNB - NAVY BLUE

PMTBL - BLACK

PMCB - CARRY CASE

[Click For More Information](#)

BEDSHEETS

FITTED DISPOSABLE BEDSHEET BX100



The best form of protection for your examination couch is a disposable fitted bed sheet. Suitable for various treatment couches and stretchers, made to size and packaged in four packs of 25 for easy distribution it is the preferred solution for many healthcare facilities.

Specifications:

2000 x 750mm
Non-woven polypropylene
Available in navy blue & white
40 GSM
Box of 100 (4 x 25)

ORDER CODES:

BSHFITWH - White

BSHFITNB - Navy Blue

[Click For More Information](#)

NON-FITTED DISPOSABLE BEDSHEET BX100

The best form of protection for your examination couch is a disposable flat bed sheet. Suitable for various treatment couches, made to size and packaged in ten packs of 10 for easy distribution it is the preferred solution for many healthcare facilities.

Specifications:

2400 x 700mm
Non-woven polypropylene
Available in navy blue & white
40 GSM
Box of 100 (4 x 25)

ORDER CODE:

BSHWH - White

BSHNB - Navy Blue

[Click For More Information](#)



TREATMENT COUCH ACCESSORIES & SPARE PARTS

PAPER TOWEL HOLDER FOR ELECTRIC & FIXED HEIGHT COUCHES



A paper towel holder that fixes to electric & fixed height couches.

A paper towel holder that fixes to electric & fixed height couches.

Order Code:

PTH

[Click For More Information](#)

GAS STRUT SHORT FOR ELECTRIC COUCH

Gas strut short for hi-lo couch. Used on many headrest and footrest sections on electric couches.

Gas strut short for hi-lo couch. Used on many headrest and footrest sections on electric couches.

Order Code:

ETGSS

[Click For More Information](#)



GAS STRUT LONG FOR ELECTRIC COUCHES



Gas strut long for electric couches. Used on many backrest sections on electric couches.

Gas strut long for electric couches. Used on many backrest sections on electric couches.

Order Code:

ETGSL

[Click For More Information](#)

HAND CONTROL FOR ELECTRIC COUCH HI-LO

A hand control for activation of the couch Hi-Lo function.

A hand control for activation of the table Hi-Lo function.

Order Code:

ETHC

[Click For More Information](#)



FOOT CONTROL RED DOT FOR ELECTRIC COUCH HI-LO



A foot control used to activate the electric couch hi-lo function.

A foot control used to activate the electric couch hi-lo function.

Order Code:

ETFCRD

[Click For More Information](#)

FOOT CONTROL YELLOW DOT FOR ELECTRIC COUCH TILT

A foot control with yellow dot used to activate the electric couch tilt function. Used on many of our tables such as the Podiatry chair.

A foot control with yellow dot used to activate the electric couch tilt function. Used on many of our tables such as the Podiatry chair.

Order Code:

ETFCYD

[Click For More Information](#)





CONTROL BOX FOR ELECTRIC COUCHES

Control box for electric couches.

Control box for electric couches.

Order Code:

ETCB

[Click For More Information](#)

SPLITTER FOR ELECTRIC COUCHES

A splitter to allow two foot or hand controls to be connected to the table.

A splitter to allow two foot or hand controls to be connected to the table. Allows tables with two motorised functions to operate e.g. Podiatry seat tilt and Five Section with PD.

Order Code:

ETSP

[Click For More Information](#)



INDEX

B

Bariatric Three Section All Electric Treatment Couch 2
Beauty Day Spa Couch 26

C

Chiropractic Premium Fixed Height Table 18
Chiropractic Premium Table 17
Chiropractic Table PRO Series 16
Contoured Massage Couch 23
Control Box For Electric Couches 38

F

Fitted Disposable Bedsheet BX100 33
Five Section Treatment Couch No Postural Drainage 19
Five Section Treatment Couch With Postural Drainage 20
Foot Control Red Dot For Electric Couch Hi-Lo 37
Foot Control Yellow Dot For Electric Couch Tilt 37
Four Section Traction Treatment Couch 25
Free Standing Treatment Couch 30
Free Standing Treatment Couch With Drawers 31

G

Gas Strut Long For Electric Couches 36
Gas Strut Short For Electric Couch 35
Gynaecology Premium Treatment Couch 15
Gynaecology Treatment Couch 14

H

Hand Control For Electric Couch Hi-Lo 36

N

Neurological Bobath Treatment Couch 22
Non-Fitted Disposable Bedsheet BX100 33

P

Paper Towel Holder For Electric & Fixed Height Couches 35
Podiatry / Multi Purpose Chair 13
Portable Massage Couch 32
Procedure Chair (REPLACED BY PROCEDURE CHAIR PRO SERIES) 7
Procedure Chair PRO Series 8

S

Single Section Treatment / Change Couch With Side Rails 28
Splitter For Electric Couches 38
Stealth Flexion Distraction Chiropractic Table (REPLACED BY FLEXION DISTRACTION PRO SERIES) 3
Stealth Flexion Distraction Chiropractic Table PRO Series 4
Stealth Medical Treatment Couch 5
Stealth Physio Treatment Couch 6

T

Three Section All Electric Treatment Couch 11
Three Section Medical Treatment Couch 10
Three Section Medical Treatment Couch With Side Rails 27

Three Section Physio Treatment Couch 12
Tilt Table Two Section 24
Two Section Cardiology Couch 21
Two Section Treatment Couch 9

V

Veterinary Lift Table 29



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au