



Pacific Medical Catalogue

TABLE OF CONTENTS

Portable Massage Couches 2

PORTABLE MASSAGE COUCHES

PORTABLE MASSAGE COUCH



Aluminium framed portable massage couch with an auto lock cabling system, rounded corners and a multi density 5cm foam system. A great choice for mobile practitioners.

Features:

- Aluminium frame and legs
- Manufactured without electric welding to join the aluminium
- Face crest & prone arm rest
- Auto lock cabling system
- Rounded corners and deluxe no noise handle
- Luxurious PU vinyl
- Adjustable backrest
- Multi density 5cm foam system

Specifications:

- Height Adjustable: 65-80cm
- Length: 185cm
- Width: 66cm
- Folded Size: 93x67x17cm
- Weight: 18KG
- Weight Capacity: 220KG
- Warranty: 2 Years

ORDER CODE:

PMTNB - NAVY BLUE

PMTBL - BLACK

PMCB - CARRY CASE

[Click For More Information](#)

INDEX

P

Portable Massage Couch 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au