



Pacific Medical Catalogue

TABLE OF CONTENTS

Veterinary 2

VETERINARY



STAINLESS STEEL OXYGEN TROLLEY SINGLE DRAWER

A high grade single drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Single drawer
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawer
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE: SST01OXY

[Click For More Information](#)

STAINLESS STEEL OXYGEN TROLLEY TWO DRAWER

A high grade two drawer stainless steel trolley with extended shelf and strap to support an oxygen bottle up to E in size.

Features:

- High grade 304 stainless steel throughout
- Two drawers
- Extended bottom shelf with strap for oxygen bottle support
- Double safety folds
- Smooth welding
- Soft closing magnetic drawer
- German ball bearing drawer runners
- Smooth running lockable castors
- Top and bottom rails on three sides
- Outstanding finish

Specifications:

- 125mm deep drawers
- 95mm lockable castors
- Support for oxygen bottle up to E size
- Dimensions: 50 x 75 x 90cm (LxWxH)

ORDER CODE: SST02OXY

[Click For More Information](#)



STAINLESS STEEL DEEP WASH TUB

A hardy multi purpose product manufactured from 304 stainless. Suitable for examination, grooming, bathing, dental examinations and much more.

Features:

- Durable 304 stainless steel
- Suitable for grooming and examination
- Removable and sliding top grill
- Removable bottom grill
- Removable legs for easy transport and storage

Specification:

- Length: 167cm
- Width: 58cm
- Height: 90cm



Deep End Depth: 47.5cm
Shallow End Depth: 6cm

ORDER CODES:
SSWTD

[Click For More Information](#)

VETERINARY LIFT TABLE

A robust and mobile animal lift table with stainless steel surface. Suitable for exam and treatment with battery operation to suit any size clinic.

Features:

- Electric Hi Lo operation by hand switch & footbar
- Durable 304 stainless steel patient surface
- Sleek new design with rounded frame legs for safety
- 10cm lockable castors for safety and mobility
- Reliable TiMotion motors
- Rechargeable battery

Specification:

- Total Length: 110cm
- Width: 56cm
- Height range from 41cm to 96cm in approx 30 seconds
- Warranty: 7 Years Construction, 3 Years Electrics, 12 Months Footswitch
- Certification: CE, FDA, ISO13485-2016
- Weight Capacity: 250KG

ORDER CODES:
ETVET

MAINTENANCE: Please make sure your table is serviced yearly to ensure your table functions to its optimum. Lubrication of pivot points and moving points on your table must be done on a regular basis.

[Click For More Information](#)



INDEX

S

Stainless Steel Deep Wash Tub 3

Stainless Steel Oxygen Trolley Single Drawer 2

Stainless Steel Oxygen Trolley Two Drawer 3

V

Veterinary Lift Table 4



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au