



# Pacific Medical Catalogue

TABLE OF CONTENTS

Walmay Medical Products ..... 2

    Walmay Liquid Nitrogen Dewars ..... 2

    Walmay Liquid Nitrogen Flask ..... 3

# WALMAY MEDICAL PRODUCTS

## WALMAY LIQUID NITROGEN DEWARs



### WALMAY LIQUID NITROGEN DEWARs

Designed specifically for liquid nitrogen storage, Walmay Liquid Nitrogen Dewar's are characterised by their small calibers and low static evaporation loss.

**Features:**

- High strength aluminium alloy
- Multi layered thermal insulation
- Low static evaporation loss
- Protective jacket to prevent collision

**Specification:**

Model No.	WM-D10L	WM-D20L
Capacity	10.0 L	20.0 L
Empty Weight	6.2 kg	11.2 kg
Neck Diameter	50 mm	50 mm
Outer Diameter	303 mm	409 mm
Height	530 mm	650 mm
Static Evaporation Loss**	0.10 L/d	0.10 L/d
Static Hold Time	101 d	202 d

**ORDER CODES:**

**WMD10L - Walmay Dewar 10L**

**WMD20L - Walmay Dewar 20L**

[Click For More Information](#)

# WALMAY LIQUID NITROGEN FLASK



## WALMAY CRYOSPRAY

### Features:

- Fixed stainless steel probe with greater tensile strength
- Safety valve maintaining consistent pressure in the unit
- Increased regulator travel for better control of liquid nitrogen flow
- Ergonomic collar design
- Improved centre of gravity for easy handling
- Fine gauze strainer on the siphon tube
- Carry bag
- Regulator position can be left or right sided

### Specification:

- Capacity: 350ml
- Includes 5 spray tips and holder
- Warranty: 5 Years
- Australian designed and manufactured

### ORDER CODES:

**WM227 - Walmay Cryospray 300ml**

**MAINTENANCE:** Please make sure your flask is serviced yearly to ensure it functions to its optimum.

[Click For More Information](#)

# INDEX

## W

Walmay Cryospray 3

Walmay Liquid Nitrogen Dewars 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)