



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Walmay Liquid Nitrogen Dewars** ..... 2

# WALMAY LIQUID NITROGEN DEWARS



## WALMAY LIQUID NITROGEN DEWARS

Designed specifically for liquid nitrogen storage, Walmay Liquid Nitrogen Dewar's are characterised by their small calibers and low static evaporation loss.

### Features:

- High strength aluminium alloy
- Multi layered thermal insulation
- Low static evaporation loss
- Protective jacket to prevent collision

### Specification:

Model No.	WM-D10L	WM-D20L
Capacity	10.0 L	20.0 L
Empty Weight	6.2 kg	11.2 kg
Neck Diameter	50 mm	50 mm
Outer Diameter	303 mm	409 mm
Height	530 mm	650 mm
Static Evaporation Loss**	0.10 L/d	0.10 L/d
Static Hold Time	101 d	202 d

### ORDER CODES:

**WMD10L - Walmay Dewar 10L**

**WMD20L - Walmay Dewar 20L**

[Click For More Information](#)

# INDEX

## W

Walmay Liquid Nitrogen Dewars 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)