



Pacific Medical Catalogue

TABLE OF CONTENTS

Walmay Liquid Nitrogen Dewars 2

WALMAY LIQUID NITROGEN DEWARS



10L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 10L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

Made of a high-strength aluminium alloy
Advanced powder coating technology
Designed for low liquid nitrogen consumption
5 year vacuum warranty
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

Capacity: 10L
Weight (empty): 6.0kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 300mm
Height of dewar: 552mm
Static evaporation rate: 0.12L/d
Static holding time: 86 days
Ladle outer diameter: 38mm
Ladle height: 120mm
No. ladles per dewar: 3~6
Warranty: 5 Years Vacuum, 2 Years Container

Straws holding capacity:
Single 120mm (0.5mL): 792pc
Single 120mm (0.25mL): 1788pc

ORDER CODE: LNC10

[Click For More Information](#)

20L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 20L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

Made of a high-strength aluminium alloy
Advanced powder coating technology
Designed for low liquid nitrogen consumption
5 year vacuum warranty
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

Capacity: 20L
Weight (empty): 9.5kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 394mm
Height of dewar: 672mm
Static evaporation rate: 0.12L/d
Static holding time: 168 days
Ladle outer diameter: 38mm
Ladle height: 120/276mm
No. ladles per dewar: 3~6

Straws holding capacity:
Single 120mm (0.5mL): 792pc
Single 120mm (0.25mL): 1788pc
Double 276mm (0.5mL): 1284
Double 276mm (0.25mL): 2832

ORDER CODE:

LNC20

[Click For More Information](#)





30L LIQUID NITROGEN DEWAR

You can rest assured your liquid nitrogen will be safely stored inside our high performance 30L dewar with a safety lock lid and low liquid nitrogen consumption.

Features:

Made of a high-strength aluminium alloy
Advanced powder coating technology
Designed for low liquid nitrogen consumption
5 year vacuum warranty
Each unit comes complete with safety lock lid, neck plug, leather (PU) bag & ladle

Specifications:

Capacity: 30L
Weight (empty): 12.9kg
Caliber (neck opening of dewar): 50mm
Outer diameter of dewar: 462mm
Height of dewar: 706mm
Static evaporation rate: 0.12L/d
Static holding time: 254 days
Ladle outer diameter: 38mm
Ladle height: 120/276mm
No. ladles per dewar: 3~6

Straws holding capacity:
Single 120mm (0.5mL): 792pc
Single 120mm (0.25mL): 1788pc
Double 276mm (0.5mL): 1284
Double 276mm (0.25mL): 2832

ORDER CODE:
LNC30

[Click For More Information](#)

INDEX

1

10L Liquid Nitrogen Dewar 2

2

20L Liquid Nitrogen Dewar 3

3

30L Liquid Nitrogen Dewar 4



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au