



Pacific Medical Catalogue

TABLE OF CONTENTS

Water Distillers 2

WATER DISTILLERS

WATER DISTILLER BENCHTOP

This compact water distiller provides a suitable water supply solution for many practice autoclaves.

Specifications:

Stainless Steel Tank

750w Power

Produces 1L Per Hour

Internal Dimensions: 180x200mm

External Dimensions: 290x290x390mm

Weight: 3.5KG

CE Mark

ORDER CODE:

WDIS

[Click For More Information](#)

WATER DISTILLER CHARCOAL FILTER


Replacement carbon filter for the WDIS water distiller.

Replacement carbon filter for the WDIS water distiller.

ORDER CODE:

WDISCF

[Click For More Information](#)





WATER DISTILLER CLEANING CRYSTALS

A tub of crystals used in the cleaning of benchtop water distillers.

A tub of crystals used in the cleaning of benchtop water distillers.

Tub capacity: 150g

Tub size: 72x50x150mm

For instructions on use see water distiller user manual.

Order Code:

WDISCC

[Click For More Information](#)

WATER DISTILLER REPLACEMENT BOTTLE

Replacement water bottle for the WDIS water distiller.

Replacement water bottle for the WDIS water distiller.

Order Code:

WDISB

[Click For More Information](#)



INDEX

W

Water Distiller Benchtop 2

Water Distiller Charcoal Filter 2

Water Distiller Cleaning Crystals 3

Water Distiller Replacement Bottle 3



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au