



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**Wheelchairs ..... 2**

# WHEELCHAIRS



## STANDARD WHEELCHAIR 18"

A strong economical wheelchair with 110kg capacity. A great addition to any clinic.

### Features:

- Single cross brace powder coated frame
- Front rubber tires
- Pneumatic 24" spoke wheel rear tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate

### Specifications:

- Total Length 103cm
- Total Width 63cm
- Total Height 91cm
- Width Folded 26.5cm
- Seat Width 45cm
- Seat Depth 40cm
- Backrest Height 44cm
- Seat Height 50cm
- Weight Capacity 110kg
- Front Wheel Diameter 8" x 6"
- Rear Wheel Diameter 24"
- Gross Weight 22.5kg
- Net Weight 20Kg
- Carton Size 82 x 27 x 43cm

### ORDER CODE: 18WCS

[Click For More Information](#)

## PATIENT MOVER WHEELCHAIR 18"

Our patient mover wheelchair is a great aid to ferry patients around the largest of facilities with ease.

### Features:

- Single cross brace, powder coating
- Front PVC tires
- Rear PU tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate

### Specifications:

- Total Length 97cm
- Total Width 63.5cm
- Total Height 92cm
- Width Folded 27cm
- Seat Width 46cm
- Seat Depth 41cm
- Backrest Height 43cm
- Seat Height 49cm
- Weight Capacity 110kg
- Front Wheel Diameter 8"
- Rear Wheel Diameter 12"
- Gross Weight 16.4kg
- Net Weight 14.4kg
- Carton Size 72 x 25 x 93

### ORDER CODE: 18WCM

[Click For More Information](#)



## BARIATRIC WHEELCHAIR 22"

Every practice needs a strong wheelchair to carry heavy patients. With a foldable design, 200kg weight capacity and comfortable seat it presents great value.

### Features:

- Single cross brace powder coated frame
- Front PU tires
- Spoke rear wheel with PU tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate
- Nylon with foam seat & backrest

### Specifications:

Total Length: 103cm  
Total Width: 73cm  
Total Height: 91cm  
Width Folded: 26.5cm  
Seat Width: 55cm  
Seat Depth: 40cm  
Backrest Height: 44cm  
Seat Height: 50cm  
Weight Capacity: 200kg  
Front Wheel Diameter: 8" x 2"  
Rear Wheel Diameter: 24"  
Gross Weight: 23kg  
Net Weight: 20.5kg  
Carton Size: 82 x 27 x 93

**ORDER CODE:**  
**22WCB**

[Click For More Information](#)

# INDEX

## **B**

Bariatric Wheelchair 22" 3

## **P**

Patient Mover Wheelchair 18" 3

## **S**

Standard Wheelchair 18" 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)