



Pacific Medical Catalogue

TABLE OF CONTENTS

Wheelchairs 2

WHEELCHAIRS



STANDARD WHEELCHAIR 18"

A strong economical wheelchair with 110kg capacity. A great addition to any clinic.

Features:

- Single cross brace powder coated frame
- Front rubber tires
- Pneumatic 24" spoke wheel rear tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate

Specifications:

- Total Length 103cm
- Total Width 63cm
- Total Height 91cm
- Width Folded 26.5cm
- Seat Width 45cm
- Seat Depth 40cm
- Backrest Height 44cm
- Seat Height 50cm
- Weight Capacity 110kg
- Front Wheel Diameter 8" x 6"
- Rear Wheel Diameter 24"
- Gross Weight 22.5kg
- Net Weight 20Kg
- Carton Size 82 x 27 x 43cm

ORDER CODE: 18WCS

[Click For More Information](#)

PATIENT MOVER WHEELCHAIR 18"

Our patient mover wheelchair is a great aid to ferry patients around the largest of facilities with ease.

Features:

Single cross brace, powder coating
 Front PVC tires
 Rear PU tires
 Flip back desk length armrest
 Swing away footrest with plastic footplate

Specifications:

Total Length 97cm
 Total Width 63.5cm
 Total Height 92cm
 Width Folded 27cm
 Seat Width 46cm
 Seat Depth 41cm
 Backrest Height 43cm
 Seat Height 49cm
 Weight Capacity 110kg
 Front Wheel Diameter 8"
 Rear Wheel Diameter 12"
 Gross Weight 16.4kg
 Net Weight 14.4kg
 Carton Size 72 x 25 x 93

ORDER CODE: 18WCM

[Click For More Information](#)



BARIATRIC WHEELCHAIR 22"

Every practice needs a strong wheelchair to carry heavy patients. With a foldable design, 200kg weight capacity and comfortable seat it presents great value.

Features:

Single cross brace powder coated frame
 Front PU tires
 Spoke rear wheel with PU tires
 Flip back desk length armrest
 Swing away footrest with plastic footplate
 Nylon with foam seat & backrest

Specifications:

Total Length: 103cm
Total Width: 73cm
Total Height: 91cm
Width Folded: 26.5cm
Seat Width: 55cm
Seat Depth: 40cm
Backrest Height: 44cm
Seat Height: 50cm
Weight Capacity: 200kg
Front Wheel Diameter: 8" x 2"
Rear Wheel Diameter: 24"
Gross Weight: 23kg
Net Weight: 20.5kg
Carton Size: 82 x 27 x 93

**ORDER CODE:
22WCB**

[Click For More Information](#)

INDEX

B

Bariatric Wheelchair 22" 3

P

Patient Mover Wheelchair 18" 3

S

Standard Wheelchair 18" 2



Pacific Medical Australia

E info@pacificmedical.com.au

P 1300 076 332

www.pacificmedical.com.au