



# Pacific Medical Catalogue

# TABLE OF CONTENTS

**Wheelchairs** ..... 2

# WHEELCHAIRS



## STANDARD WHEELCHAIR 18"

A strong economical wheelchair with 110kg capacity. A great addition to any clinic.

### Features:

- Single cross brace powder coated frame
- Front rubber tires
- Pneumatic 24" spoke wheel rear tires
- Flip back desk length armrest
- Swing away footrest with plastic footplate

### Specifications:

- Total Length 103cm
- Total Width 63cm
- Total Height 91cm
- Width Folded 26.5cm
- Seat Width 45cm
- Seat Depth 40cm
- Backrest Height 44cm
- Seat Height 50cm
- Weight Capacity 110kg
- Front Wheel Diameter 8" x 6"
- Rear Wheel Diameter 24"
- Gross Weight 22.5kg
- Net Weight 20Kg
- Carton Size 82 x 27 x 43cm

### ORDER CODE: 18WCS

[Click For More Information](#)

## PATIENT MOVER WHEELCHAIR 18"

Our patient mover wheelchair is a great aid to ferry patients around the largest of facilities with ease.

### Features:

Single cross brace, powder coating  
 Front PVC tires  
 Rear PU tires  
 Flip back desk length armrest  
 Swing away footrest with plastic footplate

### Specifications:

Total Length 97cm  
 Total Width 63.5cm  
 Total Height 92cm  
 Width Folded 27cm  
 Seat Width 46cm  
 Seat Depth 41cm  
 Backrest Height 43cm  
 Seat Height 49cm  
 Weight Capacity 110kg  
 Front Wheel Diameter 8"  
 Rear Wheel Diameter 12"  
 Gross Weight 16.4kg  
 Net Weight 14.4kg  
 Carton Size 72 x 25 x 93

### ORDER CODE: 18WCM

[Click For More Information](#)



## BARIATRIC WHEELCHAIR 22"

Every practice needs a strong wheelchair to carry heavy patients. With a foldable design, 200kg weight capacity and comfortable seat it presents great value.

### Features:

Single cross brace powder coated frame  
 Front PU tires  
 Spoke rear wheel with PU tires  
 Flip back desk length armrest  
 Swing away footrest with plastic footplate  
 Nylon with foam seat & backrest

### Specifications:

Total Length: 103cm  
Total Width: 73cm  
Total Height: 91cm  
Width Folded: 26.5cm  
Seat Width: 55cm  
Seat Depth: 40cm  
Backrest Height: 44cm  
Seat Height: 50cm  
Weight Capacity: 150kg  
Front Wheel Diameter: 8" x 2"  
Rear Wheel Diameter: 24"  
Gross Weight: 23kg  
Net Weight: 20.5kg  
Carton Size: 82 x 27 x 93

**ORDER CODE:  
22WCB**

[Click For More Information](#)

# INDEX

## **B**

Bariatric Wheelchair 22" 3

## **P**

Patient Mover Wheelchair 18" 3

## **S**

Standard Wheelchair 18" 2



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)