



# Pacific Medical Catalogue

**TABLE OF CONTENTS**

**X-Ray Viewers ..... 2**

# X-RAY VIEWERS



## SLIMLINE X-RAY VIEWERS

A high quality slimline LED xray viewer with 2.4cm thickness, mitsubishi screen, high luminance, 100000 hours working time and very simple operation. A great addition to any clinic.

**Features:**

Mitsubishi Screen  
2.4cm Thickness  
No Globes Required  
Lightweight  
High Luminance

**Specifications:**

10000K Colour Temperature  
100,000 Hours Working Time  
4000cd/m2 / Average 13000 Lux Luminance

Dimensions (L x W x H): 47x2.4x54.5cm (Single), 85x2.4x54.5cm (Double), 123x4.1x54.5cm (Triple), 161x4.1x54.5cm (Quad)

**Order Codes:**

**SBXSL - Single Bay**  
**DBXSL - Double Bay**  
**TBXSL - Triple Bay**  
**QBXSL - Quad Bay**

[Click For More Information](#)

## STANDARD X-RAY VIEWERS

A high quality lightweight LED xray viewer with a mitsubishi screen, high luminance and 100000 hours working time. Great for any clinic.

**Features:**

Mitsubishi Screen  
12cm Thickness  
LED Illumination  
Lightweight

**Specifications:**

10000K Colour Temperature  
100,000 Hours Working Time  
2000cd/m2 / Average 6300 Lux Luminance

Dimensions: 44x12x59cm (Single), 80x12x59cm (Double),  
116x12x59cm (Triple), 150x12x59cm (Quad)

**Order Codes:**

**SBXST - Single Bay**  
**DBXST - Double Bay**  
**TBXST - Triple Bay**  
**QBXST - Quad Bay**

[Click For More Information](#)



# INDEX

## S

Slimline X-Ray Viewers 2

Standard X-Ray Viewers 3



**Pacific Medical Australia**

**E** [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

**P** 1300 076 332

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)