



NEW PRODUCT

INFRARED FOREHEAD THERMOMETER

CE & CFDA APPROVED

Infrared Forehead Thermometer

- Fast non-invasive measurement
- Instant one second response time
- Simple one button operation
- Measurement range 32.0°C - 43.0°C
- 0.1°C increment display
- LCD backlight display
- Can be used with or without sound
- Auto shut-off after 30 seconds
- Accuracy:
 - 32.0°C - 34.9°C: +/- 0.3°C
 - 35.0°C - 42.0°C: +/- 0.2°C
 - 42.1°C - 42.9°C: +/- 0.3°C
- Requires 2 x AA batteries

Certification: CE, CFDA, ISO13485

ORDER CODE: YT1



Pacific Medical

Phone: 1300 076 332

Email: info@pacificmedical.com.au

www.pacificmedical.com.au

