

Walmay Medical Products

presents the new and improved

CRYOSPRAY

It is the instrument of choice for all Dermatologists and Practitioners.

- Ergonomic collar design
- 350ml Flask easy for all Operators to hold.
- Fixed Stainless Steel Probe with greater tensile strength
- 5 tips, 4 Coned for accurate aim during procedure and 1 curved tip for body orifices.
- Regulator Handle may be placed on either side of the probe to suit left or right hand Operator

- Centre of gravity is towards the top of the flask, making it easy to balance and steady the instrument during use.
- Stainless Steel and Brass components for high quality instrument.
- Safety Valve maintains consistent pressure in the unit for operation.
- Carry Bag for safe and convenient transport.



5 YEAR WARRANTY MADE IN AUSTRALIA

